CE₀₁₂₃

Crosstella™ OTW

PTA Balloon Dilatation Catheter
Cathéter de dilatation à ballonnet pour PTA
PTA-Balloon dilatacijskatheter
Catéter de dilatación de balón para ATP
Cateter de Dilatação de Balão para ATP
Catetere con palloncino di dilatazione per PTA
PTA ballondilatatiekatheter
PTA ballongdilatacijskatheter
PTA ballon-dilatationskatheter
PTA ballongdilatacijskatheter
PTA-pallolaajennuskatheter
Καθετήρας Διαστολής PTA με Μπαλόνι
Периферический баллонный катетер ЧТА
Balonowy cewnik dylatacyjny do PTA
PTA Tágító Ballonkathéter
Balónkový dilatační katétr pro PTA
PTA balónkový dilatačný katéter
PTA Balon Dilatasyon Kateteri
PTA balloondilatatsiooni kateeter
PTA dilatācijas balonkathetrs
PTA kateteris su išplečiamuoju balionėliu
Balonski kateter za dilatacijs PTA
PTA balon dilatacijskatheter
Cateter de dilatare pe balon pentru PTA
PTA дилатационен балон-катетър
Дилатаційний балон-катетер для ЧТА
Balonski dilatacijski kateter za PTA

September, 2021

PM-05581
1/115 (3212-2)

REF	SN	LOT	
<p>Catalogue number Référence du catalogue Bestellnummer Número de catálogo Referência Codice prodotto Catalogusnummer Artikelnummer Varenummer Artikelnummer Tuotekoodi Κωδικός Προϊόντος Каталожный номер Kod produktu Termékkód Kód výrobku Kód výrobku Ürün kodu Kood Koda Nr. Katalogo Nr. Koda izdelka Šifra proizvoda Codul produsului Каталожен номер Номер за каталогом Kataloški broj</p>	<p>Serial number Numéro de série Seriennummer Número de serie Número de série Numero di serie Seriennummer Seriennummer Seriennummer Seriennummer Serianumero Αριθμός παραγωγής Серийный номер Numer seri Sorozatszám Sériové číslo Sériové číslo výrobku Seri Numarasi Seerianumber Sērījas Nr. Serijos Nr. Serijska številka Serijski broj Număr de serie. Сериен номер Серийный номер Serijski broj</p>	<p>Lot number Numéro de lot Chargenbezeichnung Número de lote Número de lote Numero di lotto Lotnummer Batchkod Batchnummer Batch number Eränumero Αριθμός παρτίδας Номер партии Numer serii Tételszám Číslo výrobní serie Identifikačné číslo Lot numarasi Partii nr. Sērījas Nr. Partijos kodas Številka serije Serijska oznaka Număr lotului Партиден номер Код партії Šifra serije</p>	<p>Use by A utiliser jusque Verwendbar bis Caducidad Data limite de uso Utilizzare entro Gebruiken voor Användes före Udløbsdato Brukes før-dato Eräpäivä Ημερομηνία Λήξης Годен до Użyty przed datą Felhasználható Použit do Použit' do Son kullanım tarihi Säilivusaeg Derīguma termiņš Sunaudoti iki Uporabno do Upotrebljivo do Data de expirare Срок на годност Використати до Upotrijebiti do</p>

STERILE/O		P
<p>Sterilized using ethylene oxide Stérilisé à l'oxyde d'éthylène Sterilisation mit Ethylenoxid Esterilizado con óxido de etileno Esterilizado com oxido de etileno Sterilizato con ossido di etilene Gesteriliseerd met ethyleenoxid Steriliserad med etylenoxid Ætylenoxid steriliseret Sterilisert med etylen oksid Steriloitu etyleenoksidilla Αποστειρωμένο με αιθυλινοξειδιο Стерилизовано окисью этилена Steryliżowane tlenkiem etylenu Etilénoxiddal sterilizálva Sterilizováno ethylenoxidem Sterilizované etylénoxidom Etilen Oksitle sterilize edilmiştir Steriliseeritud etüleenoksiidiga Sterilizēts ar etilēna oksīdu Sterilizuota naudojant etileno oksidą Sterilizirano z etilen oksidom Sterilisano etilen oksidom Sterilizat cu oxid de etilena Стерилизовано с этилен оксид Простерилизовано за допомогою оксиду етилену Sterilizirano etilen-oksidom</p>	<p>Do not reuse Ne pas réutiliser Nicht zur Wiederverwendung No reutilizar Não reutilizar Monouso Niet opnieuw gebruiken Ingen återanvändning Må ikke genbruges Må ikke gjenbrukes Kertakäyttöinen Μόνο για μια χρήση Не использовать повторно Nie używać powtórnie Csak egyszeri használatra Nelze použít opakovaně Na jednorázové použitie Tek kullanımlık Ainult ühekordseks kasutuseks Tikai vienreizējai lietošanai Pakartotina nenaudoti Samo za enkratno uporabo Samo za jednokratnu upotrebu Doar pentru unica folosinta Само за еднократна употреба Використовувати повторно заборонено Ne upotrebjavati ponovno</p>	<p>Pressure Pression Druck Presión Pressão Pressione Druk Tr yck Tr yk Trykk Paine Πίεση Давление Ciśnienie Nyomás Tlak Tlak Basiņç Rõhk Spiediens Slegis Tlak Pritisak Presiune Налягане Тиск Tlak</p>

	
<p>Consult instructions for use Lire le mode d'emploi Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultar instruções de uso Consultare le istruzioni per l'uso De gebruiksaanwijzing raadplegen Läs bruksanvisningen Læs brugsanvisningen Les bruksanvisning før bruk Lue käyttöohjeet Διαβάστε τις οδηγίες χρήσης Прочитайте инструкции по применению Sprawdź w instrukcji stosowania Olvassa el a használati útmutatót! Postupujte podle pokynů Prečítajte si návod na použitie Kullanma talimatına başvurunuz Tutvu kasutusjuhendiga. Pirms lietošanas izlasiet instrukciju Skaitykite naudojimo instrukciją Glej navodila za uporabo Pogledaj instrukcije za upotrebu A se consulta instrucțiunile de folosință Прочетете указанията за употреба Ознайомтеся з інструкціями для застосування Pogledati upute za upotrebu</p>	<p>Manufacturer Fabricant Hersteller Fabricante Fabricante Fabbricante Fabrikant Tiliverkare Produceret af Produsent Valmistaja Κατασκευαστής Изготовлено Wytwórca Gyártó Vyrobeno Výrobca Üretici firma Toodetud Ražots Gamintojas Izdelovalec Proizvodjač Fabricat de Произведено от Виробник Proizvođač</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;">NP</div>	<div style="border: 1px solid black; padding: 2px; display: inline-block;">RBP</div>
<p>Nominal Pressure Pression nominale Nenndruck Presión nominal Pressão nominal Pressione nominale Nominale druk Nominellt tryck Nominelt tryk Nominelt trykk Nominaalipaine Όνομαστική Πίεση Номинальное давление Ciśnienie znamionowe Névleges nyomás Jmenovitý tlak Menovitý tlak Nominal Basınç Nominaalrõhk Nominālais spiediens Minimalus slėgis Nominalni tlak Nominalni pritisk Presiune nominală Номинално налягане Номинальний тиск Nominalni tlak</p>	<p>Rated Burst Pressure Pression nominale de rupture Zerreißdruck Presión de rotura Pressão nominal de ruptura Pressione massima di sicurezza Bepaalde scheurdruk Nominellt bristningstryck Nominel øvre trykgrænse Nominelt sprengtrykk Nimellinen rikkoutumispain Όνομαστική Πίεση Ρήξης Расчетное давление разрыва Znamionowe ciśnienie rozrywające Becsült szakítónyomás Jmenovitý tlak prasknutí Menovitý tlak prasknutia Anma Patlama Basıncı Arvestuslik purunemisrõhk Aprēķinātais sprādziena spiediens Minimalus išplėtimo slėgis Naznačen delovni tlak Granični pritisk pucanja Presiunea nominală de explozie Номинално налягане на разкъсване Номинальний тиск розриву Nazivni tlak prsnuća</p>

 <p>Contents</p>	
<p>Contents Contenu Stück Unidades Conteúdo Unità Stuks Styck Indhold Innhold Sisältö Περιεχόμενα Содержание Zawartość Tartalom Obsah Obsah İçindekiler Sisu Satus Sudétis Vsebina Sadržaj Conținut Зміст Вміст Sadržaj</p>	<p>Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Inhalt beschädigter Packung nicht verwenden No usar si el paquete está dañado Não usar se a embalagem estiver danificada Non utilizzare se la confezione risulta danneggiata Niet gebruiken wanneer de verpakking beschadigd is Använd inte om förpackningen är skadad Må ikke anvendes, hvis pakningen er beskadiget Må ikke brukes hvis forpakningen er ødelagt Älä käyttää jos pakkaus on vahingoittunut Μην το χρησιμοποιήσετε εάν η συσκευασία είναι κατεστραμμένη Не использовать, если упаковка повреждена Nie stosować w przypadku uszkodzenia opakowania Ne használja, ha a csomagolás sérült! Nepoužívejte, bylo-li balení poškozeno Nepoužívajte, ak je obal poškodený Paket hasarlıysa kullanmayınız Ärge kasutage, kui pakend on kahjustatud Nelietot, ja sterilais iepakojums ir bojāts Jei pakuotė pažeista, nenaudoti Ne uporabljajte ga, če je ovojnina poškodovana. Ne upotrebljavati ukoliko je pakovanje oštećeno A nu se utiliza în cazul în care ambalajul este deteriorat Не използвайте, ако опаковка е повредена Не використовувати, якщо пакування пошкоджено Ne upotrebljavati ako je pakiranje oštećeno</p>
	
<p>Balloon diameter Diamètre du ballon Durchmesser des Ballons Diámetro del balón Diâmetro do balão Diametro del palloncino Ballondiameter Ballongdiameter Ballondiameter Ballongdiameter Pallon halkaisija Διάμετρος μπαλονιού Диаметр баллона Średnica balonu Ballon átmérő Průměr balónku Priemer balóna Balon çapı Ballooni diameeter Balona diametrs Balono diametras Premer balona Prečnik balona Diametrul balonului Диаметър на балона Диаметр балону Promjer balona</p>	<p>Balloon length Longueur du ballon Länge des Ballons Longitud de balón Comprimento do Balão Lunghezza del palloncino Ballonlengte Ballonglängd Ballonlængde Ballonglengde Pallon pituus Μήκος Μπαλονιού Длина баллона Długość balonu Ballon hossz Délka Balonku Dĺžka balóna Balon Uzunluđu Ballooni pikkus Balona garums Balono ilgis Dolžina balona Dužina balona Lungimea balonului Дължина на балона Довжина балону Duljina balona</p>

	
<p>Authorized Representative in the European Community Rep. CE (Mandataire européen) Bevollmächtigter in der Europäischen Gemeinschaft Representante autorizado de la Comunidad europea Representante autorizado na CE Rappresentante autorizzato nella comunità europea Gemachtigde vertegenwoordiger in de EU (Europese Unie) Auktoriserad återförsäljare inom EU-området Autoriseret forhandler indenfor EU området Autoriseret representant i Det europeiske Fællesmarked EY edustaja Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα Авторизованный представитель Европейского сообщества Autoryzowany przedstawiciel w Unii Europejskiej Az Európai Közösségben a képviseletre jogosult Zástupce Evropského společenství Predstaviteľ Európskej únie Yetkili AB temsilcisi Esindaja EL-s ES pārstāvis Įgaliotas atstovas Europos bendrijoje Zastopnik v Evropski Skupnosti Predstavnik u Evropskoj Uniji Reprezentant UE Упълномощен представител в ЕС Уповноважений представник в Европейському Союзу Ovlašteni predstavnik u europskoj zajednici</p>	<p>Do not reesterilize Ne pas restériliser Nicht wieder sterilisieren No reesterilizar Não re-esterilizar Non risterilizzare Niet hersteriliseren Får ej omsteriliseras Må ikke reesteriliseres Må ikke re-steriliseres Ei uudelleen steriloitua Na μην επαναποστερωθεί Не стерилизовать повторно Nie sterylizować ponownie Nem sterilizálható újra! Nesterilizujte opakovane Nesterilizujte Yeniden sterilize etmeyin Ärge steriliseerge uuesti Nesterilizēt atkārtoti Pakartotinai nesterilizuoti Ne reesterilizirati Ne reesterilisati Nu reesterilizați Да не се стерилизира повторно Повторно не стерилизувати Ne ponovo sterilizirati</p>
	
<p>Manufacturing Site Site de fabrication Herstellungsort Lugar de fabricación Site do fabricante Sito di produzione Productielocatie Tillverkningsplats Produktionssted Produksjonssted Valmistuspaikka Τόπος Κατασκευής Произведено Miejsce produkcji Gyártóhely Misto výroby Miesto výroby Üretim Yeri Tootmiskohht Ražotne Gamybos vieta Proizvodno mesto Mesto proizvodnje Locul de Productie Място на производство Виробнича дільниця Mjesto proizvodnje</p>	<p>Keep dry Conserver au sec Vor Nässe schützen Mantener seco Manter seco Conservare in luogo asciutto Droog bewaren Håll torr Holdes tørt Holdes tørt Pidä kuivassa Διατηρήστε το στεγνό Беречь от влаги Przechowywać w miejscu suchym Nedvességétől védett helyen tartandó Udržujte v suchu Udržujte ja suchom mieste Kuru tutunuz Hoida kuivas Sargāt no mitruma Laikyti sausoje vietoje Izogibajte se dežju Čuvati na suvom mestu A se păstra în locuri ferite de umezeală Съхранявайте на сухо Зберігати в сухому місці Držati suho</p>

**OTW**

Keep away from sunlight
Craint la chaleur
Vor Sonnenlicht schützen
Mantener alejado de la luz solar
Manter ao resguardo da luz solar
Tenere lontano da fonti di calore
Niet blootstellen aan zonlicht
Aktas för solljus
Må ikke udsættes for sollys
Unngå sollys
Suojattava auringonvaloita
Φυλάξτε το μακριά από το ηλιακό φως
Избегать прямого солнечного света
Chronić przed bezpośrednim działaniem promieni słonecznych
Napfényről védett helyen tartandó
Chraňte před slunečním zářením
Udržujte mimo slnečného žiarenia
Güneş ışığından uzak tutunuz
Kaitsta otsese päikesevalgus eest
Sargāt no saules stariem
Saugoti nuo saulės šviesos
Izogibajte se sončni svetlobi
Ne izlagati sunčevog svetlosti
A se pästra la umbrā
Пазете от слънчева светлина
Зберігати подалі від сонячного проміння
Držati dalje od sunčevog svjetla

Over-the-Wire type
Coaxial
Over-the-wire-Type
Coaxial
Sobre o Fio-Guia
Tipo Over-the-Wire
Over-de-voerdraad type
Över ledare
Over the wire-type
Over the wire
Over-the-wire
Τύπου Πάνω Από Το Σύρμα
Тип На Проводнике
Typ over the wire
Végig vezetődrróton futó típus
Aplicace přes vodící drát
Typ s použitím cez drôt
Tel Üstü Tipi
Kasutamiseks juhtetraadi peal
Lietošanas ar stīgu
Ant vielos įvedamas tipas
Tip OTW (preko žice)
OTW sistem dostave, koaksijalni
Tipul coaxial
тип По Жицата
тип OTW
OTW tip, koaksijalni

Rx ONLY

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
[PRODUCT DESCRIPTION]

Crosstella OTW is a Percutaneous Transluminal Angioplasty (PTA) Balloon Dilatation Catheter for peripheral indications. Crosstella OTW is an over the wire (OTW) type and the maximum diameter of the compatible guidewire is 0.018 inches (0.46 mm).

The shaft consists of an outer tube and an inner tube (co-axial lumens structure). The outer lumen is used for inflation and deflation of the balloon with contrast medium diluted with saline solution. The inner lumen (the guidewire lumen; from the distal tip to the guidewire port) is for inserting a compatible guidewire to facilitate advancing the catheter through the stenotic lesion or stent to be dilated. Inside the balloon, two radiopaque markers, which indicate the working length of the balloon, are placed to guide the physician for positioning the balloon properly in the targeted lesion under fluoroscopy. Crosstella OTW is available in various balloon sizes by its diameter and length at recommended inflation pressures (Nominal pressure). Upon inflation, the balloon diameter varies according to the inflating pressure. The balloon compliance chart of Crosstella OTW is provided in Table 1 at the end of this Instructions for Use. In the sterilization package, the device is equipped with the protective materials, a balloon protective tube and a stylet, which are to be removed before use. The balloon protective tube is placed over the balloon folded in a low profile and the stylet is inserted into the guidewire lumen to prevent collapsing the catheter shaft. The catheter's distal tip is tapered to facilitate the catheter's advancing into the stenotic lesion or stent. The surface of catheter is partially coated with hydrophilic polymer coating which generates lubricity when wet.

<NAME OF EACH PART>

Sterile and non pyrogenic in an undamaged and unopened blister pouch. This device is sterilized by ethylene oxide.

Recommended guidewire diameter:

Maximum diameter: 0.018”(0.46 mm)

Recommended inflation pressure and maximum inflation pressure

Nominal pressure: 8 atm

Rated burst pressure: 14 atm

[INDICATIONS]

The Crosstella OTW PTA balloon dilatation catheter is intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Note: The verification test for stent post-dilatation of Crosstella OTW was conducted using Complete® SE Vascular Stent System (Medtronic, Inc), S.M.A.R.T.® CONTROL® Vascular Stent System (Cordis Corporation).

[CONTRA-INDICATIONS]

1. Patients who have developed anastomotic stenosis within one month after AV shunt construction.
2. Patients having lesions which communicate with a pseudoaneurysm.

3. Patients who are pregnant.
4. Patients who cannot tolerate antiplatelet therapy or anticoagulant therapy.

[WARNINGS AND PRECAUTIONS]

[Warnings]

1. This device is for single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, the biocompatibility and the physical integrity of the device.
2. Do not use if the product or the unit packaging has been damaged.
3. Do not inflate the balloon exceeding the diameter of the blood vessel proximal and distal to the stenotic lesion.
4. Do not inflate the balloon to a pressure exceeding the rated burst pressure.
5. Do not use this catheter in the central circulatory system.

[Precautions related to procedures]

1. Do it carefully and slowly to remove this device out of the carrier tube and to remove protection parts from this device. Careless handling may damage the balloon and/or the catheter shaft of this device to impair balloon inflation and/or deflation.
2. Inflate the balloon with the contrast medium diluted with saline in the ratio of 1:1 (hereafter, inflation fluid). No gaseous medium such as air should be used for inflation. (Balloon may be inflated irregularly.)
3. Prior to use, expel all the air in the balloon and balloon inflation lumen and replace it with the inflation fluid. (In case of incomplete air removal, balloon inflation state can not be observed under fluoroscopy.)
4. In case of serious stenotic lesions such as calcified lesions, the blood vessel can not be dilated fully. Do not apply a pressure exceeding the rated burst pressure in such a case. (The balloon may burst and the debris may remain inside the body.)
5. As for the medical devices used in conjunction with this catheter, follow the instructions for use of such devices. (When the balloon catheter is used for post-dilatation of stent during peripheral stent placement procedure, refer to the manufacturer's instructions for use.)
6. For insertion into stent, withdrawal from stent or in-stent dilatation of the balloon, proceed with caution under fluoroscopy. (Doing so without fluoroscopy may damage this catheter or injure the blood vessel.)
7. Do not attempt to pass the Crosstella OTW PTA balloon dilatation catheter through a smaller sized introducer sheath than indicated on the label. Refer to product labeling.

[Precautions during usage]

1. Use the device immediately after the sterile package is opened.
2. The entire procedure should be carried out aseptically.
3. This catheter should be used only by physicians skilled in percutaneous vascular therapy.
4. Select the catheter with appropriate balloon size (balloon diameter, balloon length) according to the following (procedures) criteria.
 - 1) The inflated balloon diameter shall not exceed the inner vessel diameter proximal and distal to the lesion.
 - 2) The length of the inflated balloon shall not exceed the length of the lesion to be treated. (Inappropriate catheter selection may lead to acute vascular occlusion.)
5. Use this catheter only when emergency surgery can be performed at any time.

6. The physician in charge of the procedure should determine the duration and number of balloon inflations based on his/her past experiences.
7. Heparinized saline should be infused for anti-coagulation while this balloon catheter is inserted in the blood vessel.
8. Never inflate the balloon or advance the guidewire with the catheter bent or kinked.
9. Manipulate the catheter carefully in the blood vessel verifying the location and movement of its tip under fluoroscopy or DSA (Digital Subtraction Angiography) monitoring.
10. Always inflate the balloon while observing its inflation state under fluoroscopy. Inflate the balloon carefully while verifying its movement with the radiopaque markers. If any abnormality is found such as a little resistance, immovable catheter during manipulation or kinked distal tip of the guidewire, the procedure should be discontinued immediately and this catheter should be removed slowly together with the guidewire.
11. Do not insert or remove the catheter rapidly. (Operating rapidly may damage the catheter or injure the vascular intima.)
12. Do not move the catheter with the balloon inflated. (Moving the catheter with the balloon inflated may result in balloon burst or catheter shaft breakage.)
13. Do not inflate or deflate the balloon rapidly in the blood vessel. (Rapid inflation or deflation may damage the blood vessel or cause the balloon to burst resulting in the debris left inside the body.)
14. Do not use an injector (automatic injection device) for inflating the balloon.
15. Precautions should be taken to prevent any damage to the catheter by a surgical knife or scissors.
16. If any abnormality such as strong resistance is experienced while manipulating the catheter, the procedure should be discontinued immediately. The cause should be verified and appropriate measures should be taken. (Continuing the operation with excessive force may result in damage to the catheter or in vascular wall injury.)
17. Always pay attention to the patient's condition and monitor the body temperature, pulse, and breathing. When any abnormality is found, discontinue the use of the catheter immediately or take appropriate measures for the patient's condition on the discretion of the physician.
18. After use dispose the catheter as medical waste according to hospital procedures.
19. Do not use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of the catheter.
20. While manipulating this catheter, do not twist or rotate the catheter.
21. While inserting this device into the blood vessel or removing this device from the blood vessel, make sure that the balloon is completely deflated. (A device with larger and longer balloon requires a longer deflation time.)
22. If resistance is felt during post procedural withdrawal of this device, it is recommended to withdraw the entire system together with the introducer / guide sheath.
23. Always use the catheter along a guidewire inserted into the guidewire lumen.

- Inflation device with manometer
- Luer-lock syringe
- Three-way stopcock

[Operation method or instructions for use]

1. Selection of the catheter
Prior to use, verify the blood vessel inner diameter proximal and distal to the stenotic lesion under fluoroscopy, and select the catheter with appropriate balloon size. If two sizes are applicable to the diameter verified, select the smaller size.
2. Preparations
 - 1) After aseptically removing the catheter from the package container, remove the protective materials from the catheter.
 - 2) Replace the air in the balloon and balloon inflation lumen with inflation fluid according to the following procedures.
 - a) Attach the inflation device (not included in this kit) filled with the inflation fluid onto the balloon inflation port.
 - b) After applying negative pressure for approximately 10 seconds to the inflation device, direct the catheter shaft tip downward and release the negative pressure slowly to allow the fluid to fill the balloon and balloon inflation lumen and to expel the air.
 - c) When any air is found in the balloon and balloon inflation lumen, repeat the procedure b) until the air is expelled completely.
 - d) Detach the inflation device from the balloon inflation port and remove air from the cylinder of the inflation device.
 - e) Reattach the inflation device onto the balloon inflation port and maintain negative pressure. Verify that the air no longer returns into the syringe.
 - 3) Flush the guidewire lumen with heparinized saline and fill the lumen with the heparinized saline.
 - 4) Immerse the entire Crosstella OTW PTA balloon dilatation catheter in saline.
3. Insertion and inflation of the balloon catheter
Prior to use, please check if this catheter is compatible to the guide catheter (not included in this kit), the introducer sheath (not included in this kit), and the guidewire (not included in this kit) by the following table.

Catalogue number	Balloon diameter (mm)	Compatible sheath (Fr/mm)	Maximum diameter for the guidewire (inch/mm)
BD-F20020MW ~ BD-F60200LW	2 mm, 2.5 mm, 3 mm, 4 mm, 5 mm(excluding 120 mm, 150 mm, 200 mm Balloon Length), 6 mm (20 mm, 40 mm Balloon Length)	4 Fr / 1.35 mm	0.018 inch (0.46 mm)
	5 mm (120 mm, 150 mm, 200 mm Balloon Length), 6 mm (excluding 20 mm, 40 mm Balloon Length)	5 Fr / 1.66 mm	

When using the catheter inserted in the introducer sheath or guide catheter

- 1) When using the catheter in the guide catheter, loosen the hemostatic valve while keeping it tight enough to prevent blood leakage.
- 2) Gradually insert the catheter, with the balloon deflated completely, into the introducer sheath or guide catheter over the guidewire and advance to the distal periphery of the lesion. Under fluoroscopy, advance the catheter slowly to make its radiopaque marker in the balloon reach the target site.
- 3) Under fluoroscopy, verify the position of the radiopaque marker in the balloon and determine the dilatation site.
- 4) Attach the inflation device securely to the balloon inflation port.

[Operational Instructions]

Materials typically required for PTA with the Crosstella OTW PTA balloon dilatation catheter include:

- Guidewire(s) of appropriate diameter and length
- Appropriate introducer/guide sheath and dilator set
- Vial of contrast medium
- Vial of sterile saline

- 5) Inflate the balloon gradually to the target diameter verifying the inflation pressure and dilate the stenotic lesion (for inflated balloon diameters, consult the balloon compliance chart provided in Table 1. at the end of this Instructions for Use).
4. Removal of the catheter
 - 1) After applying negative pressure to deflate the balloon completely, remove the catheter.
5. Re-insertion of the catheter
 - 1) If a re-insertion of the removed catheter with deflated balloon is required, ensure the balloon is completely deflated and folded before re-inserting it into the patient's blood vessel. On re-insertion use the catheter as per "Insertion and inflation of the balloon catheter" section herein.
 - 2) If any resistance is encountered while re-inserting the catheter, stop advancing it and carefully remove it out of the patient's blood vessel. Prepare to use a new Crosstella OTW in accordance with the procedures provided above in this section of this Instruction for Use.

[Complications]

Device Failures:

- Balloon rupture
- Insufficient inflation / deflation of the balloon
- Breakage of the balloon and / or the catheter shaft
- Difficulty in removing the device
- Leakage of inflation fluid

Adverse Events:

- Local or systemic infections
- Local internal bleeding or hematoma
- Intimal rupture
- Vascular dissection
- Vascular perforation
- Vascular rupture
- Aneurysm
- Arrhythmia
- Acute vascular occlusion
- Venous thromboembolism
- Vasospasm
- Formation of pseudoaneurysm
- Arteriovenous fistula
- Bleeding requiring transfusion
- Allergic response to the contrast media / renal failure
- Ache or pressing pain
- Death
- Embolism
- Endocarditis
- Fever
- Hypertension / Hypotension
- Inflammation
- Myocardial Infarction
- Sepsis
- Shock
- Stroke
- Transient Ischemic Attack

[Storage and expiration date]

1. Store in a cool, dry, dark place. Avoid exposure to water, direct sunlight.

2. The expiration date is indicated on the box. Do not use after the expiration date.

[Package]

1 set/box

[Name and address of manufacturer]

Manufacturer : KANEKA CORPORATION

Address: 3-18, 2-chome, Nakanoshima, Kita-ku, Osaka-city,
OSAKA, 530-8288 JAPAN
TEL No.: (+81)-(0) 6-6226-5256
FAX No.: (+81)-(0) 6-6226-5143

Authorized representative in the European Community:

KANEKA MEDICAL EUROPE N.V.

Address: Nijverheidsstraat 16, 2260 Westerlo-Oevel, BELGIUM
TEL No.: (+32)-(0) 14-256-297
FAX No.: (+32)-(0) 14-256-298

EU Distributor : TERUMO EUROPE N.V.

Address: INTERLEUVENLAAN 40,3001 LEUVEN, Belgium
TEL No.: (+32)-(0) 16 381 211
FAX No.: (+32)-(0) 16 400 249

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.

Table 1. Crosstella OTW Balloon Compliance Chart

Pressure		Balloon Diameter (mm)					
atm	kPa	2.00	2.50	3.00	4.00	5.00	6.00
8	811	2.03*	2.50*	2.98*	4.00*	5.04*	5.93*
10	1013	2.07	2.56	3.05	4.08	5.16	6.06
12	1216	2.11	2.61	3.11	4.17	5.28	6.16
14	1419	2.15**	2.67**	3.17**	4.26**	5.41**	6.28**
16	1621	2.20	2.73	3.23	4.35	5.55	6.40

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© TERUMO CORPORATION September, 2021