CE<sub>0123</sub>

# ***Metacross™ OTW***

PTA Balloon Dilatation Catheter  
Cathéter de dilatation à ballonnet pour PTA  
PTA-Balloon dilatacijskatheter  
Catéter de dilatación de balón para ATP  
Cateter de Dilatação de Balão para ATP  
Catetere con palloncino di dilatazione per PTA  
PTA ballondilatatiekatheter  
PTA ballongdilatacijskatheter  
PTA ballon-dilatationskatheter  
PTA ballongdilatacijskatheter  
PTA-pallolaajennuskatheter  
Καθετήρας Διαστολής PTA με Μπαλόνι  
Периферический баллонный катетер ЧТА  
Balonowy sewnik dylatacyjny do PTA  
PTA Tágító Ballonkathéter  
Balónkový dilatační katétr pro PTA  
PTA balónkový dilatačný katéter  
PTA Balon Dilatasyon Kateteri  
PTA balloondilatatsiooni kateeter  
PTA dilatācijas balonkatheter  
PTA kateteris su išplečiamuoju balionėliu  
Balonski kateter za dilatacijs PTA  
PTA balon dilatacijskatheter  
Cateter de dilatare pe balon pentru PTA  
PTA дилатационен балон-катетър  
Дилатаційний балон-катетер для ЧТА  
Balonski dilatacijski kateter za PTA

---

---

September, 2021

PM-05579  
1/115 (3210-2)

<b>REF</b>	<b>SN</b>	<b>LOT</b>	
<p>Catalogue number Référence du catalogue Bestellnummer Número de catálogo Referência Codice prodotto Catalogusnummer Artikelnummer Varenummer Artikelnummer Tuotekoodi Κωδικός Προϊόντος Каталожный номер Kod produktu Termékkód Kód výrobu Kód výrobu Ürün kodu Kood Koda Nr. Katalogo Nr. Koda izdelka Šifra proizvoda Codul produsului Каталожен номер Номер за каталогом Kataloški broj</p>	<p>Serial number Numéro de série Seriennummer Número de serie Número de série Numero di serie Seriennummer Seriennummer Seriennummer Seriennummer Sarjanumero Αριθμός παραγωγής Серийный номер Numer serii Sorozatszám Sériové číslo Sériové číslo výrobu Seri Numarası Seerianumber Sērījas Nr. Serijos Nr. Serjska številka Serjski broj Număr de serie. Сериен номер Серийный номер Serjski broj</p>	<p>Lot number Numéro de lot Chargenbezeichnung Número de lote Número de lote Numero di lotto Lotnummer Batchkod Batchnummer Batch number Eränumero Αριθμός παρτίδας Номер партии Numer serii Tételszám Číslo výrobní serie Identifikačné číslo Lot numarası Partii nr. Sērījas Nr. Partijos kodas Številka serije Serjska oznaka Număr lotului Партиден номер Код партії Šifra serije</p>	<p>Use by A utiliser jusque Verwendbar bis Caducidad Data limite de uso Utilizzare entro Gebruiken voor Användes före Udløbsdato Brukes før-dato Eräpäivä Ημερομηνία λήξης Годен до Użytyc przed datą Felhasználható Použit do Použit do Son kullanim tarihi Säilivusaeg Derīguma termiņš Sunaudoti iki Uporabno do Upotrebljivo do Data de expirare Срок на годност Використати до Upotrijebiti do</p>
<b>STERILE</b>		<b>P</b>	
<p>Sterilized using ethylene oxide Stérilisé à l'oxyde d'éthylène Sterilisation mit Ethylenoxid Esterilizado con óxido de etileno Esterilizado com oxido de etileno Sterilizzato con ossido di etilene Gesteriliseerd met ethyleenoxide Steriliserad med etylenoxid /Etylenoxid steriliseret Sterilisert med etylen oksid Steriloitu etyleenioksidilla Αποστειρωμένο με αιθυλινοξειδιο Стерилизовано окисью этилена Steryliżowane tlenkiem etylenu Etilénoxidál sterilizálva Sterilizováno ethylenoxidem Sterilizované etylénoxidom Etilen Oksitle sterilize edilmiştir Steriliseeritud etüleenoksiidiga Sterilizēts ar etilēna oksīdu Sterilizuota naudojant etileno oksidą Sterilizirano z etilen oksidom Sterilisano etilen oksidom Sterilizat cu oxid de etilena Стерилизовано с этилен оксид Простерилизовано за допомогою оксиду етилену Sterilizirano etilen-oksidadom</p>	<p>Do not reuse Ne pas réutiliser Nicht zur Wiederverwendung No reutilizar Não reutilizar Monouso Niet opnieuw gebruiken Ingen återanvändning Må ikke genbruges Må ikke gjenbrukes Kertakäyttöinen Μόνο για μια χρήση Не использовать повторно Nie używać powtórnie Csak egyszeri használatra Nelze použít opakovaně Na jednorázové použitie Tek kullanılmık Ainult ühekordseks kasutuseks Tikai vienreizējai lietošanai Pakartotina! nenaudoti Samo za enkratno uporabo Samo za jednokratnu upotrebu Doar pentru unica folosinta Само за еднократна употреба Використовувати повторно заборонено Ne upotrebljavati ponovno</p>	<p>Pressure Pression Druck Presión Pressão Pressione Druk Tr yck Tr yk Trykk Paine Πίεση Давление Ciśnienie Nyomás Tlak Tlak Basiņç Rõhk Spiediens Stegis Tlak Pritisak Presiune Налягане Тиск Tlak</p>	

	
<p>Consult instructions for use  Lire le mode d'emploi  Gebrauchsanweisung beachten  Consulte las instrucciones de uso  Consultar instruções de uso  Consultare le istruzioni per l'uso  De gebruiksaanwijzing raadplegen  Läs bruksanvisningen  Læs brugsanvisningen  Les bruksanvisning før bruk  Lue käyttöohjeet  Διαβάστε τις οδηγίες χρήσης  Прочитайте инструкции по применению  Sprawdź w instrukcji stosowania  Olvassa el a használati útmutatót!  Postupujte podle pokynů  Prečítajte si návod na použitie  Kullanma talimatına başvurunuz  Tutvu kasutusjuhendiga.  Pirms lietošanas izlasiet instrukciju  Skaitykite naudojimo instrukciją  Glej navodila za uporabo  Pogledaj instrukcije za upotrebu  A se consulta instrucțiunile de folosință  Прочетете указанията за употреба  Ознайомтеся з інструкціями для застосування  Pogledati upute za upotrebu</p>	<p>Manufacturer  Fabricant  Hersteller  Fabricante  Fabricante  Fabbricante  Fabrikant  Tilverkare  Produceret af  Produzent  Valmistaja  Κατασκευαστής  Изготовлено  Wytwórca  Gyártó  Vyrobno  Výrobca  Üretici firma  Toodetud  Ražots  Gamintojas  Izdelovalec  Proizvodjač  Fabricat de  Произведено от  Виробник  Proizvođač</p>
<div style="border: 1px solid black; padding: 2px; display: inline-block;"><b>NP</b></div>	<div style="border: 1px solid black; padding: 2px; display: inline-block;"><b>RBP</b></div>
<p>Nominal Pressure  Pression nominale  Nenndruck  Presión nominal  Pressão nominal  Pressione nominale  Nominale druk  Nominellt tryck  Nominelt tryk  Nominelt trykk  Nominaalipaine  Όνομαστική Πίεση  Номинальное давление  Ciśnienie znamionowe  Névleges nyomás  Jmenovitý tlak  Menovitý tlak  Nominal Basıncı  Nominaalrõhk  Nominālais spiediens  Minimalus slėgis  Nominalni tlak  Nominalni pritisak  Presiune nominală  Номинално налягане  Номинальний тиск  Nominalni tlak</p>	<p>Rated Burst Pressure  Pression nominale de rupture  Zerreißdruck  Presión de rotura  Pressão nominal de ruptura  Pressione massima di sicurezza  Bepaalde scheurdruk  Nominellt bristningstryck  Nominel øvre trykgrænse  Nominelt sprengtrykk  Nimellinen rikkoutumispain  Όνομαστική Πίεση Ρήξης  Расчетное давление разрыва  Znamionowe ciśnienie rozrywające  Becsült szakítónyomás  Jmenovitý tlak prasknutí  Menovitý tlak prasknutia  Anma Patlama Basıncı  Arvestuslik purunemisrõhk  Aprēķinātais sprādziena spiediens  Minimalus išplėtimo slėgis  Naznačen delovni tlak  Granični pritisak pucanja  Presiunea nominală de explozie  Номинално налягане на разкъсване  Номинальний тиск розриву  Nazivni tlak prsnuća</p>

 <p>Contents</p>	
<p>Contents          Contenu          Stükk          Unidades          Conteúdo          Unità          Stuks          Styck          Indhold          Innhold          Sisäito          Περιεχόμενα          Содержание          Zawartość          Tartalom          Obsah          Obsah          İçindekiler          Sisu          Satures          Sudėtis          Vsebina          Sadržaj          Conținut          Зміст          Вміст          Sadržaj</p>	<p>Do not use if package is damaged          Ne pas utiliser si l'emballage est endommagé          Inhalt beschädigter Packung nicht verwenden          No usar si el paquete está dañado          Não usar se a embalagem estiver danificada          Non utilizzare se la confezione risulta danneggiata          Niet gebruiken wanneer de verpakking beschadigd is          Använd inte om förpackningen är skadad          Må ikke anvendes, hvis pakningen er beskadiget          Må ikke brukes hvis forpakningen er ødelagt          Älä käyttää jos pakkaus on vahingoittunut          Μην το χρησιμοποιήσετε εάν η συσκευασία είναι κατεστραμμένη          Не использовать, если упаковка повреждена          Nie stosować w przypadku uszkodzenia opakowania          Ne használnia, ha a csomagolás sérült!          Nepoužívejte, bylo-li balení poškozeno          Nepoužívajte, ak je obal poškodený          Paket hasarlıysa kullanmayınız          Ärge kasutage, kui pakend on kahjustatud          Nelietot, ja sterilais iepakojums ir bojāts          Jei pakuotė pažeista, nenaudoti          Ne uporabljajte ga, če je ovojnina poškodovana.          Ne upotrebljavati ukoliko je pakovanje oštećeno          A nu se utiliza în cazul în care ambalajul este deteriorat          Не използвайте, ако опаковка е повредена          Не використовувати, якщо пакування пошкоджено          Ne upotrebljavati ako je pakiranje oštećeno</p>
	
<p>Balloon diameter          Diamètre du ballon          Durchmesser des Ballons          Diámetro del balón          Diâmetro do balão          Diametro del palloncino          Ballondiameter          Ballongdiameter          Ballondiameter          Ballongdiameter          Pallon halkaisija          Διάμετρος μπαλονιού          Диаметр баллона          Średnica balonu          Ballon átmérő          Průměr balónku          Priemer balóna          Balon çapı          Ballooni diameeter          Balona diametrs          Baliono diametras          Premer balona          Prečnik balona          Diametrul balonului          Диаметр на балона          Диаметр балону          Promjer balona</p>	<p>&lt;L&gt;          Balloon length          Longueur du ballon          Länge des Ballons          Longitud de balón          Comprimento do Balão          Lunghezza del palloncino          Ballonlengte          Ballonglängd          Ballonlængde          Ballonglengde          Pallon pituus          Μήκος Μπαλονιού          Длина баллона          Długość balonu          Ballon hossz          Délka Balonku          Dĺžka balóna          Balon Uzunluğu          Ballooni pikkus          Balona garums          Baliono ilgis          Dolžina balona          Dužina balona          Lungimea balonului          Дължина на балона          Довжина балону          Dujšina balona</p>

	
<p>Authorized Representative in the European Community Rep. CE (Mandataire européen) Bevollmächtigter in der Europäischen Gemeinschaft Representante autorizado de la Comunidad europea Representante autorizado na CE Rappresentante autorizzato nella comunità europea Gemachtigde vertegenwoordiger in de EU (Europese Unie) Auktoriserad återförsäljare inom EU-området Autoriseret forhandler indenfor EU området Autoriseret representant i Det europeiske Fællesmarked EY edustaja Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα Авторизованный представитель Европейского сообщества Autoryzowany przedstawiciel w Unii Europejskiej Az Európai Közösségben a képviselőre jogosult Zástupce Evropského společenství Predstaviteľ Európskej únie Yetkili AB temsilcisi Esindaja EL-s ES pārstāvis Įgaliotas atstovas Europos bendrijoje Zastopnik v Evropski Skupnosti Predstavnik u Evropskoj Uniji Reprezentant UE Упълномощен представител в ЕС Уповноважений представник в Европейському Союзу Ovlašteni predstavnik u europskoj zajednici</p>	<p>Do not resterilize Ne pas restériliser Nicht wieder sterilisieren No reesterilizar Não re-esterilizar Non risterilizzare Niet hersteriliseren Får ej omsteriliseras Må ikke resteriliseres Må ikke re-steriliseres Ei uudelleen steriloitua Να μην επαναποστειρωθεί Не стерилизовать повторно Nie sterylizować ponownie Nem sterilizálható újra! Nesterilizujte opakovane Nesterilizujte Yeniden sterilize etmeyin Ärge steriliseerige uuesti Nesterilizēt atkārtoti Pakartotina nesterilizuoti Ne resterilizirati Ne resterilisati Nu resterilizați Да не се стерилизира повторно Повторно не стерилізувати Ne ponovo sterilizirati</p>
	
<p>Manufacturing Site Site de fabrication Herstellungsort Lugar de fabricación Site do fabricante Sito di produzione Productielocatie Tilverkningsplats Produktionssted Produksjonssted Valmistuspaikka Τόπος Κατασκευής Произведено Miejsce produkcji Gyártóhely Misto výroby Miesto výroby Üretim Yeri Tootmiskohht Ražotne Gamybos vieta Proizvodno mesto Mesto proizvodnje Locul de Productie Място на производство Виробнича дільниця Mjesto proizvodnje</p>	<p>Keep dry Conserver au sec Vor Nässe schützen Mantener seco Manter seco Conservare in luogo asciutto Droog bewaren Håll torr Holdes tør Holdes tørt Pidä kuivassa Διατηρήστε το στεγνό Берець от влаги Przechowywać w miejscu suchym Nedvességétől védett helyen tartandó Udržujte v suchu Udržujte ja suchom mieste Kuru tutunuz Hoida kuivas Sargāt no mitruma Laikyti sausoje vietoje Izogibajte se dežju Čuvati na suvom mestu A se păstra în locuri ferite de umezeală Съхранявайте на сухо Зберігати в сухому місці Držati suho</p>

	<div style="border: 1px solid black; padding: 5px; display: inline-block;"><b>OTW</b></div>
<p>Keep away from sunlight  Craint la chaleur  Vor Sonnenlicht schützen  Mantener alejado de la luz solar  Manter ao resguardo da luz solar  Tenere lontano da fonti di calore  Niet blootstellen aan zonlicht  Aktas för solljus  Må ikke udsættes for sollys  Unngå sollys  Suojattava auringonvalolta  Φυλάξτε το μακριά από το ηλιακό φως  Избегать прямого солнечного света  Chronić przed bezpośrednim działaniem promieni słonecznych  Napfénytől védett helyen tartandó  Chraňte před slunečním zářením  Udržujte mimo slnečného žiarenia  Güneş ışığından uzak tutunuz  Kaitsta otsese päikesevalgus eest  Sargāt no saules stariem  Saugoti nuo saulės šviesos  Izogibajte se sončni svetlobi  Ne izlagati sunčevoj svetlosti  A se păstra la umbră  Пазете от слънчева светлина  Зберігати подалі від сонячного проміння  Drżać dalej od sunčevog svjetla</p>	<p>Over-the-Wire type  Coaxial  Over-the-wire-Type  Coaxial  Sobre o Fio-Guia  Tipo Over-the-Wire  Over-de-voerdraad type  Över ledare  Over the wire-type  Over the wire  Over-the-wire  Τύπου Πάνω Από Το Σύρμα  Тип На Проводнике  Typ over the wire  Végig vezetőd්රóton futó típus  Aplikace přes vodičí drát  Typ s použitím cez drôt  Tel Üstü Tipi  Kasutamiseks juhtetraadi peal  Lietojams ar stīgu  Ant vielos įvedamas tipas  Tip OTW (preko žice)  OTW sistem dostave, koaksijalni  Tipul coaxial  тип По Жицата  тип OTW  OTW tip, koaksijalni</p>

**Rx ONLY**

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
**[PRODUCT DESCRIPTION]**

Metacross OTW is a Percutaneous Transluminal Angioplasty (PTA) Balloon Dilatation Catheter for peripheral indications. Metacross OTW is an over the wire (OTW) type and the maximum diameter of the compatible guidewire is 0.035 inches (0.89 mm).

The distal section of this catheter consists of a balloon and dual lumen shafts (co-axial). The outer lumen is used for inflation and deflation of the balloon with contrast medium diluted with saline solution. The inner lumen (the guidewire lumen; from the distal tip to the guidewire port) is for inserting a compatible guidewire to facilitate advancing the catheter through the stenotic lesion or stent to be dilated. The proximal section is a double lumen shaft with dual luer port Manifold; one for inflation/deflation of the balloon, the other for guidewire lumen access. Inside the balloon, two radiopaque markers, which indicate the working length of the balloon, are placed to guide the physician for positioning the balloon properly in the targeted lesion under fluoroscopy. Metacross OTW is available in various balloon sizes by its diameter and length at recommended inflation pressures (Nominal pressure). Upon inflation, the balloon diameter varies according to the inflating pressure. The balloon compliance chart of Metacross OTW is provided in Table 1 at the end of this Instructions for Use. In the sterilization package, the device is equipped with the protective materials, a balloon protective tube and a stylet, which are to be removed before use. The balloon protective tube is placed over the balloon folded in a low profile and the stylet is inserted into the guidewire lumen to prevent collapsing the catheter shaft. The catheter's distal tip is tapered to facilitate the catheter's advancing into the stenotic lesion or stent.

<NAME OF EACH PART>

Sterile and non pyrogenic in an undamaged and unopened blister pouch. This device is sterilized by ethylene oxide.

Recommended guidewire diameter:

Maximum diameter: 0.035”(0.89 mm)

Recommended inflation pressure and maximum inflation pressure

Nominal pressure:

- Balloon diameter: smaller than 6 mm: 10 atm
- Balloon diameter: 6 mm and larger: 8 atm

Rated burst pressure:

- Balloon diameter: smaller than 5 mm: 20 atm
- Balloon diameter: 5 mm : 18 atm
- Balloon diameter: 6mm, 7 mm, 8 mm : 14 atm
- Balloon diameter: 9 mm and larger: 12 atm

**[INDICATIONS]**

The Metacross OTW PTA balloon dilatation catheter is intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

Note: The verification test for stent post-dilatation of Metacross OTW was conducted using Complete® SE Vascular Stent System (Medtronic, Inc), S.M.A.R.T.® CONTROL® Vascular Stent System (Cordis Corporation) and Epic™ Vascular Self-Expanding Stent System (Boston Scientific Corporation).

**[CONTRA-INDICATIONS]**

1. Patients who have developed anastomotic stenosis within one month after AV shunt construction.
2. Patients having lesions which communicate with a pseudoaneurysm.
3. Patients who are pregnant.
4. Patients who cannot tolerate antiplatelet therapy or anticoagulant therapy.

**[WARNINGS AND PRECAUTIONS]**

**[Warnings]**

1. This device is for single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, the biocompatibility and the physical integrity of the device.
2. Do not use if the product or the unit packaging has been damaged.
3. Do not inflate the balloon exceeding the diameter of the blood vessel proximal and distal to the stenotic lesion.
4. Do not inflate the balloon to a pressure exceeding the rated burst pressure.
5. Do not use this catheter in the central circulatory system.

**[Precautions related to procedures]**

1. Do it carefully and slowly to remove this device out of the carrier tube and to remove protection parts from this device, Careless handling may damage the balloon and/or the catheter shaft of this device to impair balloon inflation and/or deflation.
2. Inflate the balloon with the contrast medium diluted with saline in the ratio of 1:1 (hereafter, inflation fluid). No gaseous medium such as air should be used for inflation. (Balloon may be inflated irregularly.)
3. Prior to use, expel all the air in the balloon and balloon inflation lumen and replace it with the inflation fluid. (In case of incomplete air removal, balloon inflation state can not be observed under fluoroscopy.)
4. In case of serious stenotic lesions such as calcified lesions, the blood vessel can not be dilated fully. Do not apply a pressure exceeding the rated burst pressure in such a case. (The balloon may burst and the debris may remain inside the body.)
5. As for the medical devices used in conjunction with this catheter, follow the instructions for use of such devices. (When the balloon catheter is used for post-dilatation of stent during peripheral stent placement procedure, refer to the manufacturer's instructions for use.)
6. For insertion into stent, withdrawal from stent or in-stent dilatation of the balloon, proceed with caution under fluoroscopy. (Doing so without fluoroscopy may damage this catheter or injure the blood vessel.)
7. Do not attempt to pass the Metacross OTW PTA balloon dilatation catheter through a smaller sized introducer sheath than indicated on the label. Refer to product labeling.

**[Precautions during usage]**

1. Use the device immediately after the sterile package is opened.
2. The entire procedure should be carried out aseptically.
3. This catheter should be used only by physicians skilled in percutaneous vascular therapy.
4. Select the catheter with appropriate balloon size (balloon diameter, balloon length) according to the following (procedures) criteria.
  - 1) The inflated balloon diameter shall not exceed the inner vessel diameter proximal and distal to the lesion.
  - 2) The length of the inflated balloon shall not exceed the length of the lesion to be treated.

(Inappropriate catheter selection may lead to acute vascular occlusion.)

5. Use this catheter only when emergency surgery can be performed at any time.
6. The physician in charge of the procedure should determine the duration and number of balloon inflations based on his/her past experiences.
7. Heparinized saline should be infused for anti-coagulation while this balloon catheter is inserted in the blood vessel.
8. Never inflate the balloon or advance the guidewire with the catheter bent or kinked.
9. Manipulate the catheter carefully in the blood vessel verifying the location and movement of its tip under fluoroscopy or DSA (Digital Subtraction Angiography) monitoring.
10. Always inflate the balloon while observing its inflation state under fluoroscopy. Inflate the balloon carefully while verifying its movement with the radiopaque markers. If any abnormality is found such as a little resistance, immovable catheter during manipulation or kinked distal tip of the guidewire, the procedure should be discontinued immediately and this catheter should be removed slowly together with the guidewire.
11. Do not insert or remove the catheter rapidly. (Operating rapidly may damage the catheter or injure the vascular intima.)
12. Do not move the catheter with the balloon inflated. (Moving the catheter with the balloon inflated may result in balloon burst or catheter shaft breakage.)
13. Do not inflate or deflate the balloon rapidly in the blood vessel. (Rapid inflation or deflation may damage the blood vessel or cause the balloon to burst resulting in the debris left inside the body.)
14. Do not use an injector (automatic injection device) for inflating the balloon.
15. Precautions should be taken to prevent any damage to the catheter by a surgical knife or scissors.
16. If any abnormality such as strong resistance is experienced while manipulating the catheter, the procedure should be discontinued immediately. The cause should be verified and appropriate measures should be taken. (Continuing the operation with excessive force may result in damage to the catheter or in vascular wall injury.)
17. Always pay attention to the patient's condition and monitor the body temperature, pulse, and breathing. When any abnormality is found, discontinue the use of the catheter immediately or take appropriate measures for the patient's condition on the discretion of the physician.
18. After use dispose the catheter as medical waste according to hospital procedures.
19. Do not use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of the catheter.
20. While manipulating this catheter, do not twist or rotate the catheter.
21. While inserting this device into the blood vessel or removing this device from the blood vessel, make sure that the balloon is completely deflated. (A device with larger and longer balloon requires a longer deflation time.)
22. If resistance is felt during post procedural withdrawal of this device, it is recommended to withdraw the entire system together with the introducer / guide sheath.
23. Always use the catheter along a guide wire inserted into the guide wire lumen.

**[Operational Instructions]**

Materials typically required for PTA with the Metacross OTW PTA balloon dilatation catheter include:

- Guidewire(s) of appropriate diameter and length
- Appropriate introducer/guide sheath and dilator set
- Vial of contrast medium
- Vial of sterile saline
- Inflation device with manometer
- Luer-lock syringe
- Three-way stopcock

**[Operation method or instructions for use]**

1. Selection of the catheter  
Prior to use, verify the blood vessel inner diameter proximal and distal to the stenotic lesion under fluoroscopy, and select the catheter with appropriate balloon size. If two sizes are applicable to the diameter verified, select the smaller size.
2. Preparations
  - 1) After aseptically removing the catheter from the package container, remove the protective materials from the catheter.
  - 2) Replace the air in the balloon and balloon inflation lumen with inflation fluid according to the following procedures.
    - a) Attach the inflation device (not included in this kit) filled with the inflation fluid onto the balloon inflation port.
    - b) After applying negative pressure for approximately 10 seconds to the inflation device, direct the catheter shaft tip downward and release the negative pressure slowly to allow the fluid to fill the balloon and balloon inflation lumen and to expel the air.
    - c) When any air is found in the balloon and balloon inflation lumen, repeat the procedure b) until the air is expelled completely.
    - d) Detach the inflation device from the balloon inflation port and remove air from the cylinder of the inflation device.
    - e) Reattach the inflation device onto the balloon inflation port and maintain negative pressure. Verify that the air no longer returns into the syringe.
  - 3) Flush the guidewire lumen with heparinized saline and fill the lumen with the heparinized saline.
  - 4) Immerse the entire Metacross OTW PTA balloon dilatation catheter in saline.
3. Insertion and inflation of the balloon catheter  
Prior to use, please check if this catheter is compatible to the guide catheter (not included in this kit), the introducer sheath (not included in this kit), and the guidewire (not included in this kit) by the following table.

Catalogue number	Balloon diameter (mm)	Compatible sheath (Fr/mm)	Maximum diameter for the guidewire (inch/mm)
BD-L30020SW ~ BD-L70200CW	3 mm, 4 mm, 5 mm, 6 mm (excluding 200 mm Balloon Length), 7 mm (20 mm, 40 mm Balloon Length)	5Fr / 1.66 mm	0.035 inch (0.89 mm)
	6 mm (only 200 mm Balloon Length), 7 mm (excluding 20 mm, 40 mm Balloon Length), 8 mm	6 Fr / 2.00 mm	
9 mm, 10 mm, 12 mm	7 Fr / 2.33 mm		

When using the catheter inserted in the introducer sheath or guide catheter

- 1) When using the catheter in the guide catheter, loosen the hemostatic valve while keeping it tight enough to prevent blood leakage.

- 2) Gradually insert the catheter, with the balloon deflated completely, into the introducer sheath or guide catheter over the guidewire and advance to the distal periphery of the lesion. Under fluoroscopy, advance the catheter slowly to make its radiopaque marker in the balloon reach the target site.
  - 3) Under fluoroscopy, verify the position of the radiopaque marker in the balloon and determine the dilatation site.
  - 4) Attach the inflation device securely to the balloon inflation port.
  - 5) Inflate the balloon gradually to the target diameter verifying the inflation pressure and dilate the stenotic lesion (for inflated balloon diameters, consult the balloon compliance chart provided in Table 1. at the end of this Instructions for Use).
4. Removal of the catheter
- 1) After applying negative pressure to deflate the balloon completely, remove the catheter.
5. Re-insertion of the catheter
- 1) If a re-insertion of the removed catheter with deflated balloon is required, ensure the balloon is completely deflated and folded before re-inserting it into the patient's blood vessel. On re-insertion use the catheter as per "Insertion and inflation of the balloon catheter" section herein.
  - 2) If any resistance is encountered while re-inserting the catheter, stop advancing it and carefully remove it out of the patient's blood vessel. Prepare to use a new Metacross OTW in accordance with the procedures provided above in this section of this Instruction for Use.

#### [Complications]

##### Device Failures:

- Balloon rupture
- Insufficient inflation / deflation of the balloon
- Breakage of the balloon and / or the catheter shaft
- Difficulty in removing the device
- Leakage of inflation fluid

##### Adverse Events:

- Local or systemic infections
- Local internal bleeding or hematoma
- Intimal rupture
- Vascular dissection
- Vascular perforation
- Vascular rupture
- Aneurysm
- Arrhythmia
- Acute vascular occlusion
- Venous thromboembolism
- Vasospasm
- Formation of pseudoaneurysm
- Arteriovenous fistula
- Bleeding requiring transfusion
- Allergic response to the contrast media / renal failure
- Ache or pressing pain
- Death
- Embolism
- Endocarditis
- Fever
- Hypertension / Hypotension
- Inflammation

- Myocardial Infarction
- Sepsis
- Shock
- Stroke
- Transient Ischemic Attack

#### [Storage and expiration date]

1. Store in a cool, dry, dark place. Avoid exposure to water, direct sunlight.
2. The expiration date is indicated on the box. Do not use after the expiration date.

#### [Package]

1 set/box

#### [Name and address of manufacturer]

##### Manufacturer : KANEKA CORPORATION

Address: 3-18, 2-chome, Nakanoshima, Kita-ku, Osaka-city,  
OSAKA, 530-8288 JAPAN  
TEL No.: (+81)-(0) 6-6226-5256  
FAX No.: (+81)-(0) 6-6226-5143

##### Authorized representative in the European Community:

##### KANEKA MEDICAL EUROPE N.V.

Address: Nijverheidsstraat 16, 2260 Westerlo-Oevel, BELGIUM  
TEL No.: (+32)-(0) 14-256-297  
FAX No.: (+32)-(0) 14-256-298

##### EU Distributor : TERUMO EUROPE N.V.

Address: INTERLEUVENLAAN 40,3001 LEUVEN, Belgium  
TEL No.: (+32)-(0) 16 381 211  
FAX No.: (+32)-(0) 16 400 249

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.

**Table 1. Metacross OTW Balloon Compliance Chart**

Pressure (atm)	Balloon Diameter (mm)								
	3.00	4.00	5.00	6.00	7.00	8.00	9.00	10.00	12.00
8				5.96*	6.98*	8.00*	8.98*	9.90*	11.98*
10	3.01*	3.98*	5.06*	6.11	7.14	8.19	9.18	10.16	12.28
12	3.07	4.06	5.16	6.24	7.30	8.38	9.35**	10.42**	12.61**
14	3.14	4.13	5.27	6.39**	7.48**	8.58**	9.55	10.71	12.98
16	3.21	4.21	5.38	6.55	7.66	8.79			
18	3.28	4.29	5.50**						
20	3.36**	4.38**	5.64						
22	3.46	4.49	5.79						
24	3.56	4.60							

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© TERUMO CORPORATION September, 2021