Glidesheath™

Hydrophilic Coated Introducer Sheath
Gainé D'introducer Enduit Hydrophile
Hydrophile beschichtete Einführungserschleuse
Vaina Introdutora Revestido Hidrofílico
Kit introduttore idrofilo
Inbrenghuls met hydrofiele coating
Hydrofiliskt belaggd introducerskida











TERUMO CORPORATION
 44-1, 2-CHÔME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN MADE IN JAPAN

TERUMO MEDICAL CORPORATION SOMERSET NJ 08873, U.S.A.

® : Registered Trademark M Coat and Glidesheath are trademarks of TERUMO CORPORATION.
 ® : Marque Enregistrée M Coat et Glidesheath sont des marques de TERUMO CORPORATION.
 ® : Eingetragenes Warenzeichen M Coat und Glidesheath sind Marken der TERUMO CORPORATION.
 ® : Marca Registrada M Coat y Glidesheath son marcas de TERUMO CORPORATION.
 ® : Marchio Registrato M Coat e Glidesheath sono marchi della TERUMO CORPORATION.
 ® : Geregistreerd Handelsmerk M Coat en Glidesheath zijn handelsmerken van TERUMO CORPORATION.
 ® : Registrerat Varumärke M Coat och Glidesheath är varumärken tillhörande TERUMO CORPORATION.

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ENGLISH

Glidesheath™

Read these instructions before use.

DESCRIPTION/DESIGNATION

The RADIFOCUS INTRODUCER II consists of an introducer (a sheath and a dilator), a mini guide wire, an entry needle, and a syringe. The products with "M Coat™" have Terumo Corporation's exclusive hydrophilic coating on the surface of the sheath. When the surface of the sheath is wet, it becomes very lubricous and the coefficient of friction is largely reduced. The Glidesheath Introducer is used to facilitate placing a catheter through the skin into a vein or artery.

< COMPONENT DESCRIPTION >

Refer to product labelling for appropriate system components.

All components that enter blood vessels are radiopaque.

Sheath : Incorporates a 1-way valve and a 3-way stopcock connected by a side tube. The sheath can be used with a catheter of the same Fr. size or up to two Fr. sizes smaller without blood leakage at the 1-way valve. This highly flexible sheath is designed to resist kinking.

Dilator : The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath. Mini Guide Wire : The super-elastic core wire, surrounded by a smooth, hemorepellant surface coating, prevents kinking.

Syringe : 2.5 mL syringe
 Entry Needle : Inner metallic needle and a plastic cannula
 Guide Insert

WARNING/PRECAUTIONS

- WARNINGS**
- Do not use a metal cannula as an entry needle. Withdrawing the mini guide wire through a metal cannula or advancing a metal cannula over the mini guide wire may result in shearing of the mini guide wire or scraping of its plastic coating. This may lead to damage to the blood vessel or the sheath, as well as release of fragments from the wire into the blood stream.
 - Manipulate the mini guide wire slowly and carefully not to damage the vessel wall, while monitoring the tip position and movement under fluoroscopy.
- PRECAUTIONS**
- This kit must be used by a trained physician under fluoroscopy.
 - This kit is for single use only. Do not resterilize or reuse.
 - Do not use if the unit package or the product has been damaged or soiled.
 - Use immediately after opening the unit packaging.
 - Dispose of safely after single use to avoid risk of infection.
 - Before use, make sure the sheath(Fr.) and dilator size are appropriate for the access vessel and the system to be used.
 - The entire procedure from skin incision to sheath removal must be carried out aseptically.
 - Do not use power injector for contrast media injection from the side tube.
 - Use of alcohol, antiseptic solutions or other solvents must be avoided, as they may adversely affect the surface of the sheath.

CAUTION : Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

- Make a small skin incision at the puncture site with a surgical knife.
- Insert an entry needle to puncture the blood vessel. Then, remove the inner metallic needle (Fig.1) and advance the plastic cannula into the vessel.

CAUTION Do not insert the inner needle of the entry needle into the plastic cannula. Dispose the removed inner metallic needle safely being careful to avoid infections.

- Before inserting angled/spring mini guide wire, set the guide inserter at the plastic cannula hub for easy insertion. Insert the flexible end of the mini guide wire through the plastic cannula into the vessel (Fig. 2).

CAUTION Advance or withdraw the mini guide wire slowly. If resistance is met, do not advance or withdraw the mini guide wire until the cause of resistance is determined.

- Remove the plastic cannula over the mini guide wire.
- Connect a flushing line to the 3-way stopcock of the introducer sheath. Fill the sheath assembly completely with heparinized saline, removing all air. For the products with "M Coat™", also wet the surface of the sheath with heparinized saline solution.
- Prime the dilator using the syringe with heparinized saline.
- Insert the vessel dilator fully into the sheath. The female hub of the sheath connects with the male hub of the dilator, and locks in place by means of grip (Fig. 3).

CAUTIONS

- Insert the dilator into the valve center of sheath. Forced dilator insertion missing the valve center may cause valve damage, resulting in blood leakage.
- Lock the dilator hub into the sheath hub securely. If the dilator hub is not locked into the sheath hub, only the sheath will advance into the vessel and the tip of the sheath may damage the vessel.
- For the products with "M Coat™", use it in a wet state. If the surface has dried, wet before use.

- Insert the dilator and sheath together over the mini guide wire, and into the blood vessel (Fig. 4).
- "Unlock" the dilator hub from the sheath hub by bending the dilator hub downward (Fig. 5).
- Slowly remove the dilator and mini guide wire, leaving the outer sheath in the vessel (Fig. 6). If injection or sampling is necessary at this point, remove the mini guide wire, and use the dilator hub as an injection port before removing it.

CAUTION Slowly remove the dilator from the outer sheath. Rapid withdrawal of the dilator may result in incomplete closing of the 1-way valve, resulting in blood flow through the valve. If this occurs, replace the dilator into the sheath and remove again slowly. After removing the dilator and mini guide wire, be careful for advancing the sheath. It may cause damage to the vessel.

- Insert the catheter into the valve center of sheath.
- Insert a catheter through the sheath and into the blood vessel, and advance to the desired location (Fig. 7). When exchanging catheters, remove the used catheter and repeat step12.

CAUTIONS

- Before removing or inserting the catheter through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath (Fig. 8).
- Be careful for rapid aspiration by a syringe. It may cause air to be drawn through the valve.
- When puncturing, suturing, or incising the tissue near the sheath, be careful not to damage the sheath. Do not put a clamp on the sheath nor bind it with a thread.

The flushing line may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock (Fig. 9).

- When inserting, manipulating or withdrawing a catheter from the sheath always hold it in place. To temporarily suture the sheath (for continued access) use the suture eye.

CAUTION Do not place suture on the sheath tubing since this may restrict access / flow through the sheath.

- After the intended procedure is completed, remove the catheter and then the sheath.

CAUTION FOR STORAGE

DO NOT STORE AT EXTREME TEMPERATURES AND HUMIDITY.
 AVOID DIRECT SUNLIGHT.

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© TERUMO CORPORATION May 2007 07.30