


# *Pinnacle<sup>®</sup> Precision*

Access System

## SMOOTH VASCULAR ACCESS COMPARATIVE BENCH TOP RESULTS



# THE STANDARD FOR SMOOTH VASCULAR ACCESS

## Bench top testing of PINNACLE PRECISION Access System Sheath compared to Cook® Medical Check-Flo® Performer® Introducer Set.

### Results show PINNACLE PRECISION Access System Sheath provides superior performance



- Needle Penetration Force
- Transition from Guidewire-to-Dilator and Dilator-to-Sheath
- Guidewire Rail Strength
- Sheath Strength Test
- Sheath Penetration Force

#### NEEDLE PENETRATION FORCE

##### Test Protocol

Needle penetration was measured by the amount of force it took for the needle to penetrate through a 50-micron polyethylene film.

Needle Penetration Force



#### CONCLUSION

The PINNACLE PRECISION needle technology requires less penetration force compared to the competitive needle.

##### PINNACLE PRECISION

##### Cook Micropuncture

###### Top View



###### Side View



###### Bottom View



#### Construction

The PINNACLE PRECISION needle features:

- Tapered transition from 21G tip to 19G shaft to provide stability upon insertion
- Back bevel cuts on the tip of the needle help facilitate a straighter entry

#### Clinical Benefit

PINNACLE PRECISION provides access with less force designed to provide a sharper more precise entry through the skin.


## TRANSITIONS FROM GUIDEWIRE-TO-DILATOR AND DILATOR-TO-SHEATH

### Test Protocol

The sheath transition test measured the gap between the distal tip of the sheath and the dilator by bending the sheath at a 45 degree angle with a 7.5mm radius curve.

#### PINNACLE PRECISION



Guidewire-to-Dilator Transition



Dilator-to-Sheath Transition



#### Cook Micropuncture



## CONCLUSION



PINNACLE PRECISION sheaths have Total Integrated Fit technology which offers a smoother transition between guidewire-to-dilator and dilator-to-sheath providing smaller gaps between transitions.

## GUIDEWIRE RAIL STRENGTH

### Test Protocol

The rail strength of the starter guidewire was measured by testing the stiffness at 18 different points along the guidewire starting at the distal floppy tip down to the back end.

#### Guidewire Rigidity




Guidewire  
Downward Force

## CONCLUSION

PINNACLE PRECISION starter guidewire has a floppy distal tip which progressively gets stronger towards the proximal end resulting in greater rail strength of the guidewire.

#### PINNACLE PRECISION

Cross Sectional Magnified View



Total Integrated Fit Technology

Sheath  
Dilator

#### Cook Micropuncture

Cross Sectional Magnified View



Sheath  
Dilator


### Construction

PINNACLE PRECISION Total Integrated Fit tapered sheath technology features a seamless transition between the dilator and sheath.

### Clinical Benefit

PINNACLE PRECISION smooth transitions help to reduce the risk of gaps and sharp edges catching on tissue and damaging artery walls upon sheath entry, advancement or removal.

#### Guidewire Rigidity



6.5 cm - Floppy Palladium Coil with Tapered Nitinol Core  
7.0 cm - Uniform Nitinol Core  
4.0 cm - Tapered Nitinol Core

### Construction

The PINNACLE PRECISION 0.021" guidewire features:

- Floppy 6.5cm taper on the distal tip which is designed to enhance trackability and visibility
- Solid Nitinol back end to provide a strong rail for the sheath to track over

### Clinical Benefit



PINNACLE PRECISION has the only 0.021" wire on market which provides greater rail strength to help reduce the risk of buckling through scarred or calcified arteries.

## SHEATH STRENGTH TEST

### Test Protocol

The sheath strength test measured the flexural strength of the introducer sheath tube construction without the dilator by recording the amount of force required to bend the sheath in a three point bending test set-up.

Sheath Strength Test



## CONCLUSION



Pinnacle Precision has a higher sheath tube bending strength than the competitive sheaths.

## SHEATH PENETRATION FORCE

### Test Protocol


The sheath penetration test measured the amount of force needed for the dilator and sheath to penetrate a 50-micron polyethylene film. The sheaths were tracking over the starter guidewire that are provided in the respective kits.

Penetration Force



## CONCLUSION

Competitive sheath requires 77% higher insertion force to penetrate through the polyethylene film than Pinnacle Precision.



### Construction

Pinnacle Precision kink resistant material and design helps to maintain lumen patency throughout the procedure.

### Clinical Benefit

Pinnacle Precision Sheaths are designed for difficult insertion through calcified or scarred arteries to provide smooth and reliable vascular access.



### Construction

Pinnacle Precision has a combination of features which enable smooth and reliable vascular access that include:

- Greater rail strength with 0.021" starter guidewire
- Total Integrated Fit tapered sheath technology which is designed to provide smooth transitions from guidewire-to-dilator and dilator-to-sheath

### Clinical Benefit

Pinnacle Precision provides access with less force retaining a smooth dilator-to-sheath transition with less risk of damage to the sheath upon penetration.

# PUSHING BOUNDARIES

Terumo Interventional Systems **is committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**





## PINNACLE PRECISION Access System includes:

Sheath length, 10cm | Wire length, 43cm (Nitinol), 45cm (Stainless Steel) | Tapered needle gauge, 21G/19G | 10 per shelf box

CODE	SHEATH SIZE (Fr)	WIRE	NEEDLE
70-4160	4	Nitinol	Echogenic
70-4130	4	Stainless Steel	Echogenic
70-5160	5	Nitinol	Echogenic
70-5130	5	Stainless Steel	Echogenic
70-6160	6	Nitinol	Echogenic
70-6130	6	Stainless Steel	Echogenic
70-7160	7	Nitinol	Echogenic
70-7130	7	Stainless Steel	Echogenic
70-8160	8	Nitinol	Echogenic
70-8130	8	Stainless Steel	Echogenic
70-4165	4 Stiff*	Nitinol	Echogenic
70-4135	4 Stiff*	Stainless Steel	Echogenic
70-5165	5 Stiff*	Nitinol	Echogenic
70-5135	5 Stiff*	Stainless Steel	Echogenic

\*Includes a stiff dilator for greater pushability.

**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](https://www.terumo.com)

**RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.**

#### References:

Data on file at Terumo: 20170018

Note: Bench testing results may not be indicative of clinical performance.

©2023 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-07330

**TERUMO**  
INTERVENTIONAL  
SYSTEMS