**R2P RADIFOCUS GLIDEWIRE ADVANTAGE**  
(3カ国語マルチ)

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
## ENGLISH

### INSTRUCTION FOR USE

Read the following warnings, precautions and directions for use carefully.

### INDICATION FOR USE

The R2P RADIFOCUS GLIDEWIRE ADVANTAGE is designed to direct a catheter to the desired anatomical location in the peripheral vasculature during diagnostic or interventional procedures. This device is not intended for neurovascular or coronary interventions.

### CONTRAINDICATIONS

No-contraindications are known when using this device in peripheral vascular procedures. Carefully read the list of Contraindications in the instructions for use accompanying the other interventional devices to be used with the guide wire.

### COMPLICATIONS

Complications that may be related to the use of guide wire can include but may not be limited to the following:  
 • vessel perforation • intimal tear • vessel damage • vessel spasm • hematoma at the puncture site • infections

### WARNINGS/PRECAUTIONS

Failure to abide by the following warnings might result in damage to the vessel, shearing of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE, and release of plastic fragments from the R2P RADIFOCUS GLIDEWIRE ADVANTAGE. Such pieces or fragments from the wire may have to be removed from the vessel.

#### WARNINGS

- **Do not** manipulate or withdraw the R2P RADIFOCUS GLIDEWIRE ADVANTAGE through a metal entry needle or a metal dilator. Manipulation and/or withdrawal through a metal entry needle or a metal dilator may result in destruction and/or separation of the outer polyurethane coating requiring retrieval. A plastic entry needle is recommended when using this wire for initial placement.
- **Do not** use the R2P RADIFOCUS GLIDEWIRE ADVANTAGE with devices which contain metal parts such as atherectomy catheters, laser catheter, or metal introduction devices as they may cause the GLIDEWIRE ADVANTAGE plastic coating to shear and/or sever the wire.
- **Do not** reshape the R2P RADIFOCUS GLIDEWIRE ADVANTAGE by any means. Attempting to reshape the wire may cause damage, resulting in the release of wire fragments into the vessel.
- When exchanging or withdrawing a catheter over the R2P RADIFOCUS GLIDEWIRE ADVANTAGE, secure and maintain the guide wire in place under fluoroscopy to avoid unexpected guide wire advancement; otherwise damage to the vessel wall by the wire's tip may occur.
- A retrieving device, such as a gripper or basket forceps, can only be used after the R2P RADIFOCUS GLIDEWIRE ADVANTAGE has been removed from the patient's vessel. Using a retrieving device while the R2P RADIFOCUS GLIDEWIRE ADVANTAGE is in the vessel may cause the R2P RADIFOCUS GLIDEWIRE ADVANTAGE to break.
- Manipulate the R2P RADIFOCUS GLIDEWIRE ADVANTAGE slowly and carefully in the vessel while confirming the behavior and location of the wire's tip under fluoroscopy. Improper manipulation of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE without fluoroscopic confirmation may result in vessel perforation.
- **Do not** apply repetitive bending force to one specific point of the device as this may cause damage to the R2P RADIFOCUS GLIDEWIRE ADVANTAGE.
- If any resistance is felt or if the tip's behavior and/or location seems improper, stop manipulating the R2P RADIFOCUS GLIDEWIRE ADVANTAGE and/or the catheter and determine the cause by fluoroscopy. Continuing to manipulate or rotate the R2P RADIFOCUS GLIDEWIRE ADVANTAGE or failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the catheter, or damage to the vessel.
- **Do not** attempt to use the R2P RADIFOCUS GLIDEWIRE ADVANTAGE if it has been bent, kinked or damaged. Use of a damaged wire may result in damage to the vessel or the release of wire fragments into the vessel.
- Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE.

#### PRECAUTIONS

- The R2P RADIFOCUS GLIDEWIRE ADVANTAGE should be used by physicians, who are well trained in manipulation and observation of guide wires under fluoroscopy.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. **Do not** use if the unit package or the guide wire is broken or soiled. The R2P RADIFOCUS GLIDEWIRE ADVANTAGE should be used immediately after opening the package and be disposed of safely and properly after use, following local regulations for medical waste management.
- When using a drug or a device concurrently with the R2P RADIFOCUS GLIDEWIRE ADVANTAGE, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the R2P RADIFOCUS GLIDEWIRE ADVANTAGE. For example when using the R2P RADIFOCUS GLIDEWIRE ADVANTAGE with any device that emits energy (laser, pressure, ultrasound, etc.) confirm that the R2P RADIFOCUS GLIDEWIRE ADVANTAGE is retracted into a position where it **will not** be impacted by the energy.
- Consider the use of systemic heparinization.
- The surface of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE is **not** lubricous unless it is wet. Before taking it out of its holder and inserting it through a catheter, fill the holder and the catheter with heparinized physiological saline solution.
- **Do not** use a metal torque device with the R2P RADIFOCUS GLIDEWIRE ADVANTAGE. Use of a metal torque device may result in damage to the R2P RADIFOCUS GLIDEWIRE ADVANTAGE.
- **Do not** slip a tightened up torque device or Y-connector over the wire, as this may result in damage to the wire.
- Due to the slippery nature of the hydrophilic coating on the R2P RADIFOCUS GLIDEWIRE ADVANTAGE, the operator may encounter some difficulties in handling the wire. A RADIFOCUS TORQUE DEVICE, sold separately, is recommended for easier handling/manipulation of the wire.
- Due to variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, it is advisable to stop using such catheters.
- **Do not** manipulate the R2P RADIFOCUS GLIDEWIRE ADVANTAGE through a tightened up rotating hemostasis valve, as this may result in damage to the wire.
- After removal from the patient's vessel, and prior to reinserting it into the same patient during the same catheterization, the R2P RADIFOCUS GLIDEWIRE ADVANTAGE should be rinsed in a bowl full of heparinized physiological saline solution. Any blood residues still adhering to the wire can be removed by wiping once with a gauze moistened with heparinized physiological saline solution. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE.
- The entire operation should be carried out aseptically.
- The R2P RADIFOCUS GLIDEWIRE ADVANTAGE has been sterilized by ethylene oxide gas. For single use only. **Do not** resterilize or reuse.
- The R2P RADIFOCUS GLIDEWIRE ADVANTAGE is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.
- The R2P RADIFOCUS GLIDEWIRE ADVANTAGE is intended to be used under fluoroscopic guidance, determine if therapeutic benefit is outweighed by risks in individuals who are or may be pregnant.
- **Do not** use under MRI. This device is **not** MRI compatible.
- Prior to beginning radial artery access, conduct screening, such as an Allen test, to ensure that radial access is appropriate for the patient.
- The safety and effectiveness of the device is unknown in vascular regions other than those specifically indicated.

### DIRECTIONS FOR USE

#### Guide Wire


1. Remove the R2P RADIFOCUS GLIDEWIRE ADVANTAGE and the holder together from the package.
2. Fill the holder with heparinized physiological saline solution through the hub of the holder using a syringe.
3. Remove the R2P RADIFOCUS GLIDEWIRE ADVANTAGE from the holder and inspect the R2P RADIFOCUS GLIDEWIRE ADVANTAGE prior to use, to verify that it is lubricated. If the R2P RADIFOCUS GLIDEWIRE ADVANTAGE **can not** be easily removed from the holder, inject more heparinized physiological saline solution into the holder and try again.
4. Prime the catheter with heparinized physiological saline solution to ensure smooth movement of the R2P RADIFOCUS GLIDEWIRE ADVANTAGE within

the catheter.

5. The R2P RADIFOCUS GLIDEWIRE ADVANTAGE may slide entirely into the catheter or slip out of the catheter because of its low sliding friction.
6. Keep at least 5cm of the wire extended out of the hub of the catheter during introduction.

#### Guide wire clip

1. Remove the guide wire clip from its holder.
2. Shape the R2P RADIFOCUS GLIDEWIRE ADVANTAGE into a loop (Fig. 1).
3. Secure the coiled R2P RADIFOCUS GLIDEWIRE ADVANTAGE with the guide wire clip (Fig. 2).
4. When removing the guide wire clip from the coiled R2P RADIFOCUS GLIDEWIRE ADVANTAGE, repeat above steps in reverse (Fig. 3).

### PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperatures and high humidity during storage. Store under controlled room temperature.

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All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.