
**NaviCross™**

Support Catheter  
Cathéter de soutien  
Catéter de apoyo  
Cateter de Apoio

2023-02\_Rev.4

REF

#

STERILE

LOT

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TERUMO CORPORATION  
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN  
MADE IN JAPAN

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.  
Tous les noms de marques sont des marques de commerce ou des marques déposées de TERUMO CORPORATION et de leurs propriétaires respectifs.  
Todas las marcas son marcas comerciales o marcas registradas de TERUMO CORPORATION y sus respectivos dueños.  
Todos os nomes de empresas e marcas são marcas comerciais ou marcas registradas da TERUMO CORPORATION e seus respectivos proprietários.

TERUMO

© Terumo Corporation 2023-02

## ENGLISH

Please read all instructions prior to use.

## INDICATIONS

NaviCross 0.018" is intended to guide and support a guide wire during access of the peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline or diagnostic contrast agents.

## CONTRAINDICATIONS

• Patients with a contraindication to anti-platelet and/or anti-coagulation therapy.

• Women who are pregnant or are supposed to be pregnant.

• Do not use within the cerebral vasculature and the coronary arteries.

## COMPLICATIONS

• Possible complications associated to catheterization include, but are not limited to, the following.

• Bleeding complications + Allergies to drugs + Distal embolization + Hypotension

• Infection and pain at puncture site - Arteriovenous fistula + Hematoma

• Bradycardia and Palpitation - Intra-vessel thrombosis

• Femoral pseudoaneurysm/Pseudoneurysm formation - Arterial spasms

• Arterial injury: perforation/dissection - Arterial embolism/occlusion - Arrhythmia

• Hemorrhage and hemorrhagic shock - Fever /chill - Nausea and vomiting

## WARNINGS AND PRECAUTIONS

## WARNINGS

• Do not heat or bend the catheter tip. Abrasion of the hydrophilic coating or damage to the catheter may result.

• Carefully handle the product under fluoroscopy. If any resistance is felt while handling the product, immediately stop the manipulation and find out the cause to avoid damage to blood vessels and separation or breakage of the product.

• The catheter should be carried out aseptically.

• Before inserting/ withdrawing the product, clean the surface of the guide wire with gauze moistened with saline solution. Advancing/ withdrawing the product over a guide wire with residual blood on its surface or a guide wire which is not fully wet may result in separation or breakage of the product.

• Do not torque the product excessively while the distal part of the product crosses the stenosis or is in the stent.

• Be careful when inserting/ withdrawing the product through an opening of the stent struts to avoid damage to the product. A catch by the stent struts may result in separation or breakage of the product.

• In case of separation if any resistance is felt, stop the manipulation and remove the product together with the guide wire to avoid damage to the vessel. Carefully insert a guide wire into the vessel.

• After removing the catheter, the catheter should be rinsed with heparinized saline solution to remove blood residues from the catheter's surface.

• If blood residues do not come off, wipe gently once with a gauze soaked with heparinized saline solution.

• Do not use disinfectants, which can compromise the integrity of the hydrophilic coating. If necessary, flush the lumen of the catheter with heparinized saline solution to remove blood residues.

• This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties /characteristics of the drug or device and exercise due caution to avoid damage to the catheter.

• Carefully inspect the unit package to make sure that no damage to the sterile barrier has occurred.

• Consider the use of systemic heparinization.

• Infusion pressure for this product at maximum injection pressures for saline and iohexol (10.6 mPa s 37°C) are provided below. Do not exceed the maximum permissible injection pressure. Contrast media should be injected at 37°C.

## PRECAUTIONS

• This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties /characteristics of the drug or device and exercise due caution to avoid damage to the catheter.

• Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

• Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.

• Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.

• The catheter should be disposed of after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood.

• Use the product prior to the expiry date indicated on the package.

• The catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.

• Be careful when manipulating, advancing and/or withdrawing this device through stents or other devices with sharp edges, or through tortuous or calcified blood vessels. These may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vessel or device damage.

• Avoid wiping the device with dry gauze as this may damage the device coating.

• Avoid excessive wiping of the coated device.

• Do not use under MR.

• This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.

## PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature or high humidity during storage.

This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize.

Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.

Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.

The catheter should be disposed of after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood.

Use the product prior to the expiry date indicated on the package.

The catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.

Be careful when manipulating, advancing and/or withdrawing this device through stents or other devices with sharp edges, or through tortuous or calcified blood vessels. These may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vessel or device damage.

Avoid wiping the device with dry gauze as this may damage the device coating.

Avoid excessive wiping of the coated device.

Do not use under MR.

This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.

• This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties /characteristics of the drug or device and exercise due caution to avoid damage to the catheter.

• Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

• Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.

Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.

The catheter should be disposed of after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood.

Use the product prior to the expiry date indicated on the package.

The catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.

Be careful when manipulating, advancing and/or withdrawing this device through stents or other devices with sharp edges, or through tortuous or calcified blood vessels. These may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vessel or device damage.

Avoid wiping the device with dry gauze as this may damage the device coating.

Avoid excessive wiping of the coated device.

Do not use under MR.

This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.

• This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties /characteristics of the drug or device and exercise due caution to avoid damage to the catheter.

• Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

• Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.

Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.

The catheter should be disposed of after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood.

Use the product prior to the expiry date indicated on the package.

The catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.

Be careful when manipulating, advancing and/or withdrawing this device through stents or other devices with sharp edges, or through tortuous or calcified blood vessels. These may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vessel or device damage.

Avoid wiping the device with dry gauze as this may damage the device coating.

Avoid excessive wiping of the coated device.

Do not use under MR.

This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.

• This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties /characteristics of the drug or device and exercise due caution to avoid damage to the catheter.

• Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

• Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.

Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.

The catheter should be disposed of after opening the package

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