
TERUMO

2020-06 Rev.4

RADIFOCUS®

## OPTITORQUE™

Angiographic Catheter / Cathéter d' Angiographie  
Catéter Angiográfico / Cateter AngiográficoSoft tip  
Extrémité souple  
Punta suave  
Ponta flexívelBraided shaft  
Armature tressée  
Parte mallada  
Exio aramadoWing hub  
Embase à ailettes  
Cono de entrada  
Base em forma de asaWhite catheter inserter /  
Green strain relief  
Redresseur blanc /  
Manchon protecteur vert  
Insertor blanco /  
Descarga de tensión verde  
Introdutor do cateter branco /  
Redutor verde de pressãoDrawing not to scale  
Ce dessin n'est pas à l'échelle  
El dibujo no está a escala  
Desenho não à escala

## Rx ONLY

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.  
ATTENTION: Conformément au droit fédéral des Etats-Unis, ce dispositif ne peut être vendu que par ou sur prescription d'un médecin.  
ATENCIÓN: La Ley Federal (USA) restringe la venta de este dispositivo a un médico o bien por orden de éste.  
ATENÇÃO: A Lei Federal dos E.U.A. restringe a venda deste dispositivo por médicos ou por sua prescrição.Stacking limit by n  
Ne pas empiler plus de n  
No apilar más de n unidades  
Não apilar a mas de nFragile  
Fragile  
Frágil  
FracoKeep dry  
Craint l'humidité  
Mantener seco  
Proteger contra umidadeKeep Away From Sunlight  
Protéger De La Lumière Directe Du Soleil  
Mantener Alejado De La Luz Solar Directa  
Evite a luz direta do solMax guidewire outer diameter  
Longueur externe maximal du guide  
Diamètre guida  
D.E. Máximo do fio-guiaUsable length  
Longueur  
Longitud  
Comprimento utilizávelMaximum injection pressure  
Pression maximale d'injection  
Presión máxima de inyección  
Pressão maxima de injeção

## Side holes

Side holes  
Trous latéraux  
Orificios laterales  
Orificios laterais

## Tip curve L:

Tip curve length  
Longueur de l'extrémité courbe  
Longitud del extremo curvo  
Comprimento da ponta da curvaContents  
Contenu  
Contenuto  
Unidades  
Conteúdo

## REF

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TERUMO CORPORATION 44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN MADE IN JAPAN

TERUMO

© TERUMO CORPORATION 2020-06

## ENGLISH

## INSTRUCTION FOR USE

CAUTION: The RADIFOCUS OPTITORQUE Angiographic Catheter should be used by a physician who is well trained in manipulation and observation under fluoroscopy. This device is intended for single use only. Do not resterilize and/or reuse this device. Read all instructions prior to use. To avoid complications, observe all warnings and precautions throughout these instructions.

## PACKAGING / STORAGE

Terumo's catheters are sterilized by ethylene oxide. Sterile and non-pyrogenic in an unopened and undamaged packaging.

- CAUTIONS**
- Do not use if the unit package or the product has been damaged or soiled. Use immediately after opening the unit package. Dispose of safely after single use to avoid risk of infection.
  - Do not store at extreme temperatures and humidity. Avoid direct sunlight.

## INDICATION FOR USE

The RADIFOCUS OPTITORQUE Angiographic Catheter is indicated for use in cardiac and vascular procedures. It is designed to deliver radiopaque media, guide wires, catheters, and therapeutic agents to selected sites in the vascular system. The different shapes are designed to selectively engage arteries from access sites such as the femoral, radial, and brachial artery.

## DESCRIPTION

- The RADIFOCUS OPTITORQUE has a large lumen to allow contrast media to be injected at high flow rates and is designed to exhibit good torque control.
- The white catheter inserter on the pigtail catheter, can be used to straighten the catheter tip for easier insertion into the sheath.
- Fixed, green strain relief, prevents the catheter from kinking at its proximal end.

## WARNINGS AND PRECAUTIONS

- WARNINGS**
- Failure to observe these warnings may result in damage to the vessel, damage to or breakage/ separation of the catheter that may necessitate retrieval.
  - Do not heat or bend the catheter tip. Damage to the catheter may result.
  - Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the catheter.
  - Never advance the guide wire briskly and/or force it into the catheter when the catheter is bent or twisted. This may cause breakage/separation of the catheter, resulting in damage to the vessel.

## PRECAUTIONS

- Select the catheter of optimal tip shape and size, taking into account the site in which it is to be advanced to, as well as the patient's anatomy.
- When using a drug or a device with the RADIFOCUS OPTITORQUE, the operator should have a complete understanding of the properties/characteristics of the drug or device and exercise due caution to avoid damage to the catheter.
- Consider the use of systemic heparinization.
- An access device such as a sheath should be placed prior to use of this device.
- Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Due to the potential for clot formation in the absence of anticoagulation with this or any device, this catheter should always be flushed with heparinized saline solution prior to use. Systemic heparinization should be used.
- Infusion flow rates for the RADIFOCUS OPTITORQUE at maximum injection pressures for saline and iohexol (10.6mPa s 37°C) are provided below. Do not exceed the maximum injection pressure. Contrast media should be injected at 37°C.

French size	Flow rate (mL / sec)		Maximum Injection Pressure (psi/kPa)
	saline	iohexol	
4 Fr. (1.40 mm)	19	12	750psi / 5,171kPa
5 Fr. (1.70 mm)	32	21	1,000psi / 6,895kPa
6 Fr. (2.00 mm)	37	25	1,000psi / 6,895kPa

All measurements were taken using 100cm length RADIFOCUS OPTITORQUE catheters at 37°C.

## CONTRAINDICATIONS

In certain medical conditions, angiography may be contraindicated or, special precautionary measures may need to be taken prior to and/or during the procedure. The following is a list of some of these medical conditions.

- The acute phase of myocardial infarction •A serious heart failure •A serious arrhythmia •A serious systemic infection or fever •A serious disease other than coronary disease
- Serious serum electrolyte imbalance •An allergy to or reaction to from the contrast medium •Renal dysfunction •Blood coagulopathies •Some Respiratory Disorders •Mental Disease
- Pregnancy •This device should not be used in patients who cannot be systemically anticoagulated.

## COMPLICATIONS

Catheterization may be accompanied by but not limited to the following:

- Arterial embolism or obstruction •Artery dissection •Artery injury •Acute myocardial infarction •Unstable angina pectoris •False aneurysm •Ventricular fibrillation/arrhythmia
- Artery perforation •Arteriovenous fistula •Spasm •Vascular thrombosis •Distal embolism •Infection and pain at the puncture site. •Hematoma •Bradycardia •Hemorrhage

## DIRECTIONS FOR USE

- Carefully open the sterile pouch and gently remove the catheter from the package. Flush it by injecting heparinized saline solution through the catheter hub using a syringe.
 

**CAUTIONS**

  - Do not use if the catheter has been damaged or any other anomaly is observed.
  - Employ an aseptic technique during removal from the package and use.
- Insert a guide wire of appropriate size into the catheter through its hub and advance the wire to approximately 5cm beyond the catheter's distal tip.
 

**CAUTIONS**

  - The RADIFOCUS OPTITORQUE can accommodate a maximum guide wire size of 0.038" (0.97 mm) diameter.
- Gain access to the artery by using a percutaneous or cut down technique.
- Insert the guide wire alone into the artery. Then proceed to advance the catheter into the artery over the guide wire.
 

**CAUTIONS**

  - To avoid damage to the catheter after it has been advanced into the vessel, manipulate the guide wire carefully, particularly when negotiating a bend in the catheter and/or when passing through the catheter's tip.
- Manipulate the catheter slowly and carefully in the artery.
 

**WARNING**

  - Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. Failure to exercise proper caution may result in damage to the vessel or catheter. Separation of the catheter may occur requiring retrieval in some cases.
- When the catheter tip has reached the branch of the desired vessel, remove the guide wire through the catheter.
- While confirming the location of the catheter tip under fluoroscopy, advance the catheter to the desired site and perform angiography.
 

**WARNINGS**

  - Reinforce the catheter shaft with a guide wire while torque is given to the catheter. Excessive torque without use of a guide wire may result in kinking or twisting of the catheter, which may cause difficulty in catheter retrieval.
  - Before starting infusion, verify that the catheter has not been kinked or blocked. Failure to abide by this warning may cause the catheter to break/rupture/separate, resulting in damage to the vessel.

**CAUTIONS**

  - When injecting contrast media, do not exceed the maximum injection pressure of 750 psi (5,171 kPa) for the 4Fr.(1.40mm) catheter or of 1,000psi (6,895 kPa) for the 5Fr. (1.70mm) and 6Fr. (2.00mm) catheter.
- When using the catheter to guide a micro catheter, insert the micro catheter carefully into the lumen of the catheter.
 

**CAUTIONS**

  - Use a micro catheter that is below the maximum guide wire size compatible to this catheter.
- After the procedure is completed, draw the catheter back from the site. Insert the guide wire into the catheter until it extends slightly beyond the distal end of the catheter. Carefully remove the catheter and guide wire together.

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Fabricado porAll brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.