REF	Code No. Code n N° de código Código N°	STERILE EO	Sterilized Using ETO Stérilisé à l'oxyde d'éthylène Esterilizado mediante ETO (óxido de etileno) Esterilizado por ETO
LOT	Lot No. Lot n N° de lote Lote N°		Read Instructions Before Use Lire le mode d'emploi avant utilisation Leer las instrucciones antes del uso Leia as Instruções antes de Utilizar
	Exp. Date Date d'exp. Fecha de caducidad Data de Vencimento		Single Use Only A usage unique seulement Para un sólo uso Para uma Única Utilização

PINNACLE® R/O II HiFlo

INTRODUCER SHEATH WITH RADIOPAQUE MARKER
INTRODUCTEUR AVEC MARQUEUR RADIO-OPAQUE
VAINA INTRODUTORA CON MARCADOR RADIOPACO
BAINHA INTRODUTORA COM MARCADOR RADIOPACO

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The device is intended for single use only. Do not reuse. Do not resterilize. Do not reprocess as reprocessing may compromise the biocompatibility and functional integrity of the product.

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TERUMO MEDICAL CORPORATION
ELKTON MD, 21921 USA
MADE IN USA

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


DESCRIPTION AND INDICATIONS FOR USE

The Pinnacle R/O II HiFlo Introducer Sheath consists of an introducer (a sheath incorporating a radiopaque marker and a dilator), a mini guidewire and a guide inserter.

The Pinnacle R/O II HiFlo Introducer Sheath is used to facilitate placing a catheter through the skin into a vein or artery. The mini guidewire is an accessory device which is used for placement of the sheath into the vein or artery.

COMPONENT DESCRIPTION

Refer to product labeling for appropriate system components. All components that enter blood vessels are radiopaque.

Sheath

Incorporates a radiopaque marker located approximately 5 mm proximal to the tip, 1-way valve and a detachable 3-way stopcock connected by a side tube. The side tube incorporates a flow occluding pinch clamp. The sheath can be used with a catheter of the same Fr. size or up to two Fr. sizes smaller without blood leakage at the 1-way valve. This highly flexible sheath is designed to resist kinking.


Dilator

The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath.

Spring Coil Mini Guidewire

45 cm double-ended (J and straight) stainless steel guidewire. (Some catalog numbers may not include a mini guidewire.)

Guide Inserter

The contents of this package may differ from that shown in the instructions for use. Please consult available catalog information for proper contents.

PRECAUTIONS

- When using metal needle cannula, do not withdraw the guidewire back into the cannula, as shearing of the guidewire may result.
- This kit must be used by a trained physician.
- This kit is for single use only. Do not resterilize or reuse.
- Contents are sterile, non-toxic, and non-pyrogenic in unopened, undamaged package.
- Do not use if the package or product is stained or damaged. Use the introducer kit immediately after opening the package and dispose of the kit after use.
- Before use, make sure the sheath size (Fr.) is appropriate for the access vessel and the catheter to be used.
- The entire procedure from skin incision to sheath removal must be carried out aseptically.
- When inflating a balloon at, or close to, the sheath's tip, be sure not to inflate it inside the distal end of sheath. The radiopaque marker is located approx. 5 mm proximal to the tip and therefore does not actually mark the true distal tip of the sheath (Fig. 9).

- Do not use a power injector through the side tube and 3-way stopcock.
- If using for high flow rate applications, first remove the 3-way stopcock and connect directly to the female luer fitting of the side tube.

CAUTION




- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

1. Make a small skin incision at the puncture site with a surgical knife.
2. Insert a cannula into the vessel (Fig. 1).
3. Insert the selected flexible end of the mini guidewire through the cannula into the vessel (Fig. 2).

CAUTION





- Advance or withdraw the mini guidewire slowly. If resistance is met, do not advance or withdraw the mini guidewire until the cause of resistance is determined.
4. Remove the cannula over the mini guidewire.
 5. Connect a flushing line to the 3-way stopcock or the sidetubing female luer of the introducer sheath. Fill the sheath assembly completely with heparinized saline, removing all air.
 6. Prime the dilator using a syringe with heparinized saline.
 7. Insert the vessel dilator fully into the sheath. The female hub of the sheath connects with the male hub of the dilator, and locks in place by means of grip (Fig. 3).

CAUTIONS

- Insert the dilator into the center of the sheath valve. Forced insertion of the dilator which misses the center of the sheath valve may cause damage, and result in blood leakage.
- Lock the dilator hub into the sheath hub securely. Ensure the skirt of the dilator hub fits on the sides of the side tube well and not directly on top of the side tube well (Fig. 3a and 3b). If the dilator hub is not locked into the sheath hub, only the sheath will advance into the

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

Fig. 3b

Incorreto

vessel and the tip of the sheath may damage the vessel. Advancing the sheath alone may cause damage to the vessel.

8. Insert the dilator and sheath together over the mini guidewire, and into the blood vessel (Fig. 4 and 4a).
9. “Unlock” the dilator hub from the sheath hub by bending the dilator hub downward (Fig. 5).
10. Slowly remove the dilator and mini guidewire together, leaving the sheath in the vessel (Fig. 6). If injection or sampling is necessary at this point, remove the mini guidewire only, and use the dilator hub as an injection port before removing it.

CAUTION

- Slowly remove the dilator from the sheath. Rapid withdrawal of the dilator may result in the incomplete closing of the 1-way valve, resulting in blood flow through the valve. If this occurs, replace the dilator into the sheath and remove again slowly.
11. The radiopaque marker will provide sheath tip location under fluoroscopy (the radiopaque marker is located approx. 5 mm proximal to the sheath’s tip – Fig. 4a). Insert a catheter through the sheath and into the blood vessel, and advance to the desired location (Fig. 7).

Fig. 3a

Correct


Fig. 3b

Incorrect


Fig. 4

Fig. 4a

Approx. 5 mm





Fig. 5

CAUTIONS

- Before removing or inserting the catheter through the sheath, aspirate blood from the 3-way stopcock or the sidetubing female luer to remove any fibrin deposition which may have accumulated in or on the tip of the sheath (Fig. 8).
 - When puncturing, suturing, or incising the tissue near the sheath, be careful not to damage the sheath. Do not put a clamp on the sheath nor bind it with a thread.
 - When inflating a balloon at, or close to, the sheath's tip, be sure not to inflate it inside the distal end of the sheath. **The radiopaque marker is located approx. 5 mm proximal to the tip** (Fig. 4a) and therefore does not actually mark the true distal tip of the sheath (Fig. 9). The flushing line may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock or directly to the sidetubing female luer (Fig. 10).
12. When inserting, manipulating or withdrawing a catheter from the sheath, always hold it in place. To temporarily suture the sheath (for continued access) use the suture eye.

CAUTIONS

- Do not place suture on the sheath tubing since this may restrict access/flow through the sheath.
 - Do not use a power injector through the side tube and 3-way stopcock.
 - If using for high flow rate applications, first remove the 3-way stopcock and connect directly to the sidetubing female luer.
13. When exchanging catheters, remove the used catheter and repeat Step 11.
 14. The side tube may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock or directly to the female luer of the side tube (Fig. 10).
 15. The 3-way stopcock is detachable through the use of the integrated luer lock. When attaching, ensure the luer tapers are fully mated and locking collar is securely tightened.
 16. When the 3-way stopcock is removed, flow through the side tube may be occluded by depressing the top arm of the pinch clamp until the top arm is engaged by the lower lip.
 17. When the 3-way stopcock is removed and flow is occluded by the pinch clamp, flow through the side tube may be resumed by pressing the lower lip outward to release the top arm.
 18. After the intended procedure is completed, remove the catheter and then the sheath.

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© Terumo Medical Corporation
- Ⓣ Pinnacle and TERUMO are registered in the U.S. Patent and Trademark Office by TERUMO CORPORATION.
- Ⓣ Pinnacle et TERUMO sont enregistrés auprès de l'U.S. Patent and Trademark Office par TERUMO CORPORATION.
- Ⓣ Pinnacle y TERUMO son marcas registradas en la Oficina de marcas comerciales y patentes de EE.UU. (U.S. Patent and Trademark Office) por TERUMO CORPORATION.
- Ⓣ Pinnacle e TERUMO estão registradas no U.S. Patent and Trademark Office (Escritório para Registro de Marcas e Patentes dos Estados Unidos) pela TERUMO CORPORATION.