

Progreat Lambda

Peripheral Microcatheter

HydroPearl

Microspheres

AZUR CX

Peripheral Coil System

PROSTATIC ARTERY EMBOLIZATION

TAKE AIM AT BPH



TERUMO
INTERVENTIONAL
SYSTEMS

REACH A HIGHER LEVEL

Progreat Lambda

Peripheral Microcatheter

Progreat Alpha

Peripheral Microcatheter

HydroPearl

Microspheres

Indicated for Prostatic
Artery Embolization

AZUR CX

Peripheral Coil System

Available in 2-3 mm sizes

PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**

LEARN FROM THE LEADERS

EXPLORE ADVANCED TRAINING PROGRAMS WITH PAE EXPERTS

Experience immersive peer-to-peer education and training programs from **Terumo Learning Edge™** to help you expand your PAE procedural know-how through hands-on instructional techniques, including simulated learning and personalized course curricula.

TERUMO
INTERVENTIONAL
SYSTEMS

REACH A HIGHER LEVEL OF PERFORMANCE

With small vessel access and embolization solutions designed for **the way you work**

GO LONG IN RADIAL PAE PROCEDURES

Progreat Lambda

Peripheral Microcatheter

Optimal length for radial access in pelvic procedures with the ability to deliver embolic materials^{1,2}

- 1.9Fr peripheral microcatheter available in 165 cm
- Capable of delivering 0.018" detachable coils and microspheres up to 500 µm
- 900 psi for superior flow rates compared to similar sized microcatheters

Glidecath

Hydrophilic Coated Catheter

GLIDECATH® R.A.V.I. MG 2 shape is uniquely designed for internal iliac engagement via radial access

Available in 135 cm and 150 cm

Progreat Lambda

Peripheral Microcatheter

Progreat Alpha

Peripheral Microcatheter

HydroPearl


Microspheres

Indicated for Prostatic Artery Embolization

AZUR CX

Peripheral Coil System

Available in 2-3 mm sizes



Achieve distal branch selectivity with improved flexibility and multiple tip shape options²

SMALL VESSEL SOLUTIONS

Progreat Alpha

Peripheral Microcatheter

Glidewire GT

Super-Selective Hydrophilic Coated Guidewire

The standard for most small vessel embolization procedures

Compatible with GLIDEWIRE® GT 0.016"

Pre-shaped tip options for greater vessel selectivity²

TERUMO
INTERVENTIONAL
SYSTEMS

ORDERING INFORMATION

PROGREAT LAMBDA™ Peripheral Microcatheter

Fr Size	Distal ID (in)	Product Code	Length (cm)	Tip Shape	Max Pressure (psi)	RO Markers	Hydrophilic Coating Length (cm)
1.9	0.019	CC*M1916SN	165	Straight	900	1	105
		CC*M1916AN		35° Angle			
		CC*M1916DN		Triple Angle			

PROGREAT ALPHA™ Peripheral Microcatheter

Fr Size	Distal ID (in)	Product Code	Length (cm)	Tip Shape	Max Pressure (psi)	RO Markers	Hydrophilic Coating Length (cm)
2.0	0.019	MC*PC2011Y	110	Straight	750	1	50
		MC*PC2013Y	130	Straight			70
		MC*PC2013ZRA	130	70° Angle			70
		MC*PC2013ZRC	130	J Curve			70
		MC*PC2015Y	150	Straight			90
		MC*PC2015ZRA	150	70° Angle			90
		MC*PC2015ZRC	150	J Curve			90

GLIDECATH® Hydrophilic Coated Catheter MG2

Fr Size	Product Code	Length (cm)	Tip Shape	Max Pressure (psi)	Guidewire Compatibility	Hydrophilic Coating Length (cm)
4 Fr	CG441	135	MG2	750	0.038"	40
	CG442	150				
5 Fr XP	RF*ZUF511DJM	135	MG2	1000	0.038"	40
	RF*ZUF5115JM	150				

GLIDEWIRE® GT Super-Selective Hydrophilic Coated Guidewire 0.016"

Diameter	Product Code	Length (cm)	Tip Shape	Grind Length (cm)	Radiopaque Length (cm)
0.016"	RG*GW1618SP	180	Double Angle	25	2
	RG*ES1618SP		Straight/Shapeable		



FIND OUT MORE



Phone: 800.888.3786



terumo.com



Fax: 800.411.5870

INTENDED USE: The HydroPearl® Microspheres are intended for the embolization of arteriovenous malformations and hypervascular tumors, including uterine fibroids, and for embolization of prostatic arteries (PAE) for symptomatic benign prostatic hyperplasia (BPH).^{1,2}

References: 1. Progreat Lambda Microcatheter IFU 22700BZX00078_US_1_01. Revised 2021-12.
2. Data on file.

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2022 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-05654

