

# AZUR® Peripheral Embolization Systems

## Detachable Systems

### Detachable 18 System/Pack of 1 (hydrogel coil)

Product Code	Loop Diameter	Length*
45-480305	3 mm	5 cm
45-480310	3 mm	10 cm
45-480405	4 mm	5 cm
45-480410	4 mm	10 cm
45-480415	4 mm	15 cm
45-480420	4 mm	20 cm
45-480610	6 mm	10 cm
45-480615	6 mm	15 cm
45-480620	6 mm	20 cm
45-480810	8 mm	10 cm
45-480815	8 mm	15 cm
45-480820	8 mm	20 cm
45-481010	10 mm	10 cm
45-481015	10 mm	15 cm
45-481020	10 mm	20 cm
45-481215	12 mm	15 cm
45-481220	12 mm	20 cm
45-481515	15 mm	15 cm
45-481520	15 mm	20 cm
45-481530	15 mm	30 cm
45-482020	20 mm	20 cm
45-482030	20 mm	30 cm

### Detachable 18 Framing Coil System/Pack of 1 (non-hydrogel coil)

Product Code	Loop Diameter	Length*
45-680820	8 mm	20 cm
45-681026	10 mm	26 cm
45-681434	14 mm	34 cm
45-682050	20 mm	50 cm

### Detachable 35 System/Pack of 1 (hydrogel coil)

Product Code	Loop Diameter	Length*
45-450405	4 mm	5 cm
45-450410	4 mm	10 cm
45-450415	4 mm	15 cm
45-450610	6 mm	10 cm
45-450615	6 mm	15 cm
45-450620	6 mm	20 cm
45-450815	8 mm	15 cm
45-450820	8 mm	20 cm
45-451015	10 mm	15 cm
45-451020	10 mm	20 cm
45-451215	12 mm	15 cm
45-451220	12 mm	20 cm
45-451230	12 mm	30 cm
45-451520	15 mm	20 cm
45-451530	15 mm	30 cm
45-452020	20 mm	20 cm
45-452030	20 mm	30 cm

### Detachment Controller For Use With Detachable Systems/Pack of 5

Product Code	Product Description
45-4001	AZUR Detachment Controller

#### Terumo Product Code Nomenclature:

45-YXAABB  
 45 = Terumo Interventional Systems  
 X = 2 for Pushable; 4 for Detachable; 6 for Framing Coil  
 Y = 8 for 0.018"; 5 for 0.035"  
 AA = Loop diameter in mm  
 BB = Coil length in cm

#### Examples:

45-250510 stands for a pushable coil 0.035" with a loop diameter of 5 mm and a length of 10 cm.  
 45-480815 stands for a detachable coil 0.018" with a loop diameter of 8 mm and a length of 15 cm.

\*Length is calculated from the tip of the coil to the point of connection with the pusher wire when the coil is straight.

Please see reverse side for **Pushable Systems**

# AZUR® Peripheral Embolization Systems

## Pushable Systems

### Pushable 18 System/Pack of 3 (hydrogel coil)

Product Code	Loop Diameter	Length*
45-280202	2 mm	2 cm
45-280302	3 mm	2 cm
45-280304	3 mm	4 cm
45-280402	4 mm	2 cm
45-280404	4 mm	4 cm
45-280406	4 mm	6 cm
45-280504	5 mm	4 cm
45-280506	5 mm	6 cm
45-280510	5 mm	10 cm
45-280606	6 mm	6 cm
45-280610	6 mm	10 cm
45-280614	6 mm	14 cm
45-280810	8 mm	10 cm
45-280814	8 mm	14 cm
45-281014	10 mm	14 cm

### Pushable 35 System/Pack of 3 (hydrogel coil)

Product Code	Loop Diameter	Length*
45-250404	4 mm	4 cm
45-250406	4 mm	6 cm
45-250506	5 mm	6 cm
45-250510	5 mm	10 cm
45-250610	6 mm	10 cm
45-250810	8 mm	10 cm
45-250814	8 mm	14 cm
45-251014	10 mm	14 cm
45-251514	15 mm	14 cm

**Terumo Product Code Nomenclature:**

45-XYAABB

45 = Terumo Interventional Systems

X = 2 for Pushable; 4 for Detachable; 6 for Framing Coil

Y = 8 for 0.018"; 5 for 0.035"

AA = Loop diameter in mm

BB = Coil length in cm

**Examples:**

45-250510 stands for a pushable coil 0.035" with a loop diameter of 5 mm and a length of 10 cm.

45-480815 stands for a detachable coil 0.018" with a loop diameter of 8 mm and a length of 15 cm.

\*Length is calculated from the tip of the coil to the point of connection with the pusher wire when the coil is straight.

Please see reverse side for **Detachable Systems**

**Terumo – A Global Leader in Medical Devices**

For more information, call 800.862.4143 to speak to an Inside Sales Specialist or visit [www.terumo/is.com](http://www.terumo/is.com)  
Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2010 Terumo Medical Corporation. All rights reserved.

All brand names are trademarks or registered trademarks of Terumo. Lit. TIS111US r2 08/10