

# Terumo *Glidecath*® 5 Fr

## No other catheter tracks over a guidewire like a Terumo *Glidecath*



- Softer construction for superb trackability and navigation of even the most tortuous vessels

### Tapered tip

- Smoother transition between the guidewire and *Glidecath* enhances trackability

### Softer construction\*

- Single-braided stainless steel (16 strands) for even better trackability



### Terumo *Glide* Technology™ hydrophilic coating on distal portion

- Unique hydrophilic coating has been imitated but never equaled
- Enables smooth navigation and more efficient handling

**Glidecath and Glidewire® together — unsurpassed performance**

















\*Compared to the Terumo *Glidecath* 4 Fr and *Glidecath* XP

## Hydrophilic Coated Catheters

# Terumo *Glidecath*® 5 Fr

### Hydrophilic Coated Catheters

0.038" guidewire compatible and packed 5 per box.

Product Code	Catheter Size (Fr)	Total Length (in)	Flow Rate* (mL/sec)	Hydrophilic Coating Length (cm)	Tip Shape	
CG505	5	65	21	25	Straight Taper	
CG507	5	65	21	25	Angle Taper	
CG501	5	65	21	25	C1	
CG502	5	65	21	25	C2	
CG504	5	65	21	25	J-Tip	
CG506	5	100	18	40	Straight Taper	
CG508	5	100	18	40	Angle Taper	
CG503	5	100	18	40	C2	
CG510	5	100	18	40	Sim 1	
CG511	5	100	18	40	Sim 2	
CG512	5	100	18	40	Sim 3	
CG513	5	100	18	40	H1	
CG514	5	100	18	40	JB1	
CG515	5	100	18	40	JB2	
CG516	5	100	18	40	Mani	
CG509	5	100	18	40	MP	

\*5 Fr- All measurements taken at 1000 psi and 37°C using Renografin® 76.

For more information, call 800.862.4143 to speak to a Customer Care Representative, or visit [www.terumomedical.com/interventional](http://www.terumomedical.com/interventional)  
 Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2006 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit. # TIS012 04/06