

TERUMO *Glidewire*[®] Hydrophilic Coated Guidewires

The Glidewire product line you know and trust...
now expanded to give you more options for procedure success

The Terumo Glidewire product line offers an expanded range of tip shapes and configurations, making it possible for you to find the right tools for your procedures – even when you're faced with very challenging situations or vasculature.

1. Stiff Glidewire

- Extra support and guidewire pushability helps straighten tortuous vasculature for better device delivery

2. 1 cm Taper Glidewire

- Less floppy tip is ideal for crossing stenotic lesions

3. Long Taper Glidewire

- Long taper helps to maintain catheter tip shape and gain deeper vessel purchase to support guidewire advancement

4. J-Tip Glidewire


- J-curve provides atraumatic tip and avoids selection of vessel sub-branches

5. Shapeable Glidewire

- Tip can be custom shaped to match unusual vessel take-offs

6. Bolia Curve Glidewire

- Unique double angle facilitates selection of unusual vessel angles and take-offs



Consistent, reliable Terumo Glidewire performance that gets you there

With the Terumo Glide Technology™ hydrophilic coating for unmatched lubricity and a nitinol alloy core for kink resistance, the Terumo Glidewire offers the reliable performance you know and trust.

	Product Code	Wire Type	Diameter (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
Stiff Shaft Glidewire	GS3505	Stiff	0.035	80	3	Angle
	GS3506	Stiff	0.035	150	3	Angle
	GS3508	Stiff	0.035	180	3	Angle
	GS3509	Stiff	0.035	260	3	Angle
	GS3501	Stiff	0.035	150	3	Straight
	GS3503	Stiff	0.035	180	3	Straight
	GS3504	Stiff	0.035	260	3	Straight
	GS3804	Stiff	0.038	80	3	Angle
	GS3805	Stiff	0.038	150	3	Angle
	GS3807	Stiff	0.038	180	3	Angle
	GS3801	Stiff	0.038	150	3	Straight
	GS3803	Stiff	0.038	180	3	Straight
Long Taper Glidewire	GS3507	Stiff	0.035	150	5	Angle
	GS3510	Stiff	0.035	260	5	Angle
	GS3502	Stiff	0.035	150	5	Straight
	GS3806	Stiff	0.038	150	5	Angle
	GS3802	Stiff	0.038	150	5	Straight
	GR3507	Standard	0.035	150	8	Angle
	GR3502	Standard	0.035	150	8	Straight
	GR3807	Standard	0.038	150	8	Angle
	GR3802	Standard	0.038	150	8	Straight
1cm Taper Glidewire	GR3515	Standard	0.035	150	1	Angle
	GR3516	Standard	0.035	180	1	Angle
	GR3517	Standard	0.035	260	1	Angle
	GR3518	Standard	0.035	150	1	Straight
	GR3519	Standard	0.035	180	1	Straight
	GR3520	Standard	0.035	260	1	Straight
	GS3511	Stiff	0.035	150	1	Angle
Shapeable Tip Glidewire	GE3501	Standard	0.035	150	3	Straight/Shapeable
	GE3502	Standard	0.035	180	3	Straight/Shapeable
	GE3801	Standard	0.038	150	3	Straight/Shapeable
J-Tip Glidewire	GR3510	Standard	0.035	150	3	J-Tip
	GR3810	Standard	0.038	150	3	J-Tip
Bolia Curve Glidewire	GR3521	Standard	0.035	180	5	Bolia
	GR3522	Standard	0.035	180	8	Bolia

Packed 5 per box.

Terumo – A Global Leader in Medical Devices

For more information, call 800.862.4143 to speak to an Inside Sales Customer Care Representative, or visit www.terumomedical.com/interventional
Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2007 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit. #TIS028 12/07