

Terumo *Glidewire*® Access Systems

Success begins with smooth access



Glidewire Access Systems

• **0.018" Hydrophilic Coated Glidewire**

- Allows smooth passage through challenging access sites
- Super-elastic nitinol core virtually eliminates kinking

• **SURFLO® I.V. Catheter**

- 22-gauge entry needle provides a smaller puncture
- Design is compatible with Glidewire

• **Coaxial Dilators**

- In 4 Fr and 5 Fr
- Designed for smooth transitions

Also available— GlideAccess™ Systems

- Choice of 4 Fr and 5 Fr coaxial dilators
- Choice of stainless steel or platinum-tipped nitinol 0.018" guidewire
- Choice of 21-gauge needles
 - Non-echogenic
 - Echogenic for higher visibility

When access makes all the difference, This is TERUMO Territory.™

Terumo *Glidewire*® Access Systems

Micro Accessing Systems with SURFLO® I.V. Catheter

Packed 10 per box.

	Product Code	Coaxial Diameter (Fr)	Wire	Wire Length (cm)	Wire Diameter (in)
	RSM09	4	Glidewire	52.5	0.018
	RSM10	5	Glidewire	52.5	0.018

Terumo *GlideAccess*™ Systems

Micro Accessing Systems

Packed 10 per box.

	Product Code	Coaxial Diameter (Fr)	Wire	Wire Length (cm)	Wire Diameter (in)
GlideAccess Systems with 21-gauge Non-echogenic Needle	RSM01	4	Stainless Steel	40	0.018
	RSM02	5	Stainless Steel	40	0.018
	RSM03	4	Platinum-Tipped Nitinol	40	0.018
	RSM04	5	Platinum-Tipped Nitinol	40	0.018
	Product Code	Coaxial Diameter (Fr)	Wire	Wire Length (cm)	Wire Diameter (in)
GlideAccess Systems with 21-gauge Echogenic Needle	RSM05	4	Stainless Steel	40	0.018
	RSM06	5	Stainless Steel	40	0.018
	RSM07	4	Platinum-Tipped Nitinol	40	0.018
	RSM08	5	Platinum-Tipped Nitinol	40	0.018

For more information, call 800.862.4143 to speak to a Customer Care Representative, or visit www.terumomedical.com/interventional
 Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2006 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit. # TIS024 04/06