

Terumo Pinnacle®

Achieve effortless vessel access with a smooth transition

Small outside diameter and large inside diameter

- Minimizes the size of the access site without compromising ease of device passage

Kink-resistant design

- Helps keep the lumen open throughout the procedure

Smooth transition of guidewire-to-dilator and dilator-to-sheath

- Allows easy access




Unique Terumo Cross-Cut Valve maintains uncompromised hemostasis

- “Snap-on” dilator lock prevents backout during insertion
- Hub has suture eye for securing sheath to skin and is color-coded for accurate sizing

Compare Pinnacle's penetration with other sheaths¹

Average Sheath Penetration - 5 Fr Introducer Sheath by Manufacturer



Ease of insertion is indicated by the lowest amount of pressure

When access makes all the difference, This is TERUMO Territory.™

1. Data on file. Terumo Medical Corporation.

Standard Sheaths for Coronary Applications (4 Fr-11 Fr)

Terumo Pinnacle®

Introducer Sheaths

Packed 10 per box and include dilator.

Product Code	Length (cm)	Inner Diameter (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)	Mini Guidewire Length (cm)
RSS402	10	4	2.5	0.035*	45
RSS502	10	5	2.5	0.038*	45
RSS602	10	6	2.5	0.038*	45
RSS702	10	7	2.5	0.038*	45
RSS802	10	8	2.5	0.038*	45
RSS902	10	9	2.5	0.038*	45
RSS002	10	10	2.5	0.038*	45
RSS102	10	11	2.5	0.038*	45
RSS404	10	4	2.5	—	—
RSS504	10	5	2.5	—	—
RSS604	10	6	2.5	—	—
RSS704	10	7	2.5	—	—
RSS804	10	8	2.5	—	—
RSS904	10	9	2.5	—	—
RSS004	10	10	2.5	—	—
RSS104	10	11	2.5	—	—
RSS506	25	5	2.5	—	—
RSS606	25	6	2.5	—	—
RSS706	25	7	2.5	—	—
RSS806	25	8	2.5	—	—
RSS906	25	9	2.5	—	—
RSS006	25	10	2.5	—	—
RSS106	25	11	2.5	—	—

*Mini stainless steel spring coil guidewire included.

Terumo Pinnacle®

Introducer Sheaths

Packed 10 per box and include dilator.

Product Code	Length (cm)	Inner Diameter (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)	Mini Guidewire Length (cm)
RSS401	10	4	2.5	0.035*	45
RSS501	10	5	2.5	0.035*	45
RSS601	10	6	2.5	0.035*	45
RSS701	10	7	2.5	0.035*	45
RSS801	10	8	2.5	0.035*	45
RSS901	10	9	2.5	0.035*	45
RSS001	10	10	2.5	0.035*	45
RSS101	10	11	2.5	0.035*	45
RSS403	10	4	2.5	–	
RSS503	10	5	2.5	–	
RSS603	10	6	2.5	–	
RSS703	10	7	2.5	–	
RSS803	10	8	2.5	–	
RSS903	10	9	2.5	–	
RSS003	10	10	2.5	–	
RSS103	10	11	2.5	–	
RSS505	25	5	2.5	–	
RSS605	25	6	2.5	–	
RSS705	25	7	2.5	–	
RSS805	25	8	2.5	–	
RSS905	25	9	2.5	–	
RSS005	25	10	2.5	–	
RSS105	25	11	2.5	–	

*Mini stainless steel spring coil guidewire included.

For more information, call 800.862.4143 to speak to a Customer Care Representative, or visit www.terumomedical.com/interventional
Terumo Interventional Systems · 2101 Cottontail Lane · Somerset, NJ 08873 · Fax: 800.411.5870

©2006 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. Lit. # TIS020 04/06



TERUMO INTERVENTIONAL SYSTEMS