

This is TERUMO Territory.™



**Terumo Interventional Systems
Product Portfolio**

Sheaths

Pinnacle® Introducer Sheaths for Coronary Applications.....1
Pinnacle® Introducer Sheaths for Peripheral Applications2
Pinnacle TIF Tip™ Introducer Sheaths for Coronary Applications2
Pinnacle TIF Tip™ Introducer Sheaths for Peripheral Applications3
Pinnacle® R/O II Radiopaque Marker Introducer Sheaths3
Pinnacle® R/O II HiFlo Radiopaque Marker Introducer Sheaths4
Glidesheath™ Hydrophilic Coated Introducer Sheaths4
Glidesheath™ Access Kits.....4
Pinnacle® Destination® Guiding Sheaths5
GlideAccess™ Systems Micro Accessing Systems6
TR Band™ Radial Compression Devices6

Guidewires

Glidewire® Hydrophilic Coated Guidewires7-9
Glidewire® GT Super-Selective Hydrophilic Coated Guidewires9
Glidewire® Gold Hydrophilic Coated Guidewires10
Runthrough® NS Coronary Guidewires10

Angiographic Catheters

Glidecath® 4 Fr Hydrophilic Coated Catheters11
Glidecath® XP 5 Fr Hydrophilic Coated Catheters11
Glidecath® 5 Fr Hydrophilic Coated Catheters12

Coronary Catheters

Optitorque™ Coronary Diagnostic Catheters13
FineCross™ MG Coronary Micro-Guide Catheters13

Microcatheters

Progreat™ Ω Microcatheters with R/O Marker Band14
Progreat™ Microcatheters.....14
Progreat™ 2.4 Fr Microcatheters with R/O Marker Band.....14
Progreat™ Ω Coaxial Microcatheter Systems with R/O Marker Band..15
Progreat™ Coaxial Microcatheter Systems15

Embolics

AZUR™ Peripheral HydroCoil® Embolization Systems16-17
AZUR™ Peripheral Framing Coil Embolization Systems.....17
Bead Block™ Compressible Microspheres.....17

Accessories

Obturator and Torque Device.....18

Pinnacle®

Introducer Sheaths for Coronary Applications

Packed 10 per box and includes dilator.

Product Code	Length (cm)	French Size (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)
RSS402	10	4	2.5	0.035*
RSS502	10	5	2.5	0.038*
RSS602	10	6	2.5	0.038*
RSS702	10	7	2.5	0.038*
RSS802	10	8	2.5	0.038*
RSS902	10	9	2.5	0.038*
RSS002	10	10	2.5	0.038*
RSS102	10	11	2.5	0.038*
RSS404	10	4	2.5	–
RSS504	10	5	2.5	–
RSS604	10	6	2.5	–
RSS704	10	7	2.5	–
RSS804	10	8	2.5	–
RSS904	10	9	2.5	–
RSS004	10	10	2.5	–
RSS104	10	11	2.5	–
RSS506	25	5	2.5	–
RSS606	25	6	2.5	–
RSS706	25	7	2.5	–
RSS806	25	8	2.5	–
RSS906	25	9	2.5	–
RSS006	25	10	2.5	–
RSS106	25	11	2.5	–

*45 cm mini stainless steel spring coil guidewire included.
 Maximum guidewire compatible size: 0.035" for 4 Fr., 0.038" for all others.

Pinnacle®

Introducer Sheaths for Peripheral Applications

Packed 10 per box and includes dilator.

Product Code	Length (cm)	French Size (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)
RSS401	10	4	2.5	0.035*
RSS501	10	5	2.5	0.035*
RSS601	10	6	2.5	0.035*
RSS701	10	7	2.5	0.035*
RSS801	10	8	2.5	0.035*
RSS901	10	9	2.5	0.035*
RSS001	10	10	2.5	0.035*
RSS101	10	11	2.5	0.035*
RSS403	10	4	2.5	–
RSS503	10	5	2.5	–
RSS603	10	6	2.5	–
RSS703	10	7	2.5	–
RSS803	10	8	2.5	–
RSS903	10	9	2.5	–
RSS003	10	10	2.5	–
RSS103	10	11	2.5	–
RSS505	25	5	2.5	–
RSS605	25	6	2.5	–
RSS705	25	7	2.5	–
RSS805	25	8	2.5	–
RSS905	25	9	2.5	–
RSS005	25	10	2.5	–
RSS105	25	11	2.5	–

*45 cm mini stainless steel spring coil guidewire included.

Maximum guidewire compatible size: 0.035" for 4 Fr, 0.038" for all others.

Pinnacle TIF Tip™

Introducer Sheaths for Coronary Applications

Packed 10 per box and includes dilator and guidewire.

Product Code	Length (cm)	French Size (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)
10-2531	10	5	2.5	0.038*
10-2532	10	6	2.5	0.038*
10-2533	10	7	2.5	0.038*
10-2534	10	8	2.5	0.038*

*45 cm mini stainless steel spring coil guidewire included.

Maximum guidewire compatible size: 0.038".

Pinnacle TIF Tip™

Introducer Sheaths for Peripheral Applications

Packed 10 per box and includes dilator and guidewire.

Product Code	Length (cm)	French Size (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)
10-2031	10	5	2.5	0.035*
10-2032	10	6	2.5	0.035*
10-2033	10	7	2.5	0.035*
10-2034	10	8	2.5	0.035*

*45 cm mini stainless steel spring coil guidewire included.

Maximum guidewire compatible size: 0.038".

Pinnacle® R/O II

Radiopaque Marker Introducer Sheaths

Packed 10 per box and includes dilator.

Product Code	Length (cm)	French Size (Fr)	Dilator Protruding Length (cm)	Mini Guidewire Diameter (in)
RSB401	6	4	2.5	–
RSB501	6	5	2.5	–
RSB601	6	6	2.5	–
RSB701	6	7	2.5	–
RSB801	6	8	2.5	–
RSB901	6	9	2.5	–
RSB402	10	4	2.5	–
RSB502	10	5	2.5	–
RSB602	10	6	2.5	–
RSB702	10	7	2.5	–
RSB802	10	8	2.5	–
RSB902	10	9	2.5	–
RSB403	25	4	2.5	–
RSB503	25	5	2.5	–
RSB603	25	6	2.5	–
RSB703	25	7	2.5	–
RSB803	25	8	2.5	–
RSB903	25	9	2.5	–
RSB411	6	4	2.5	0.035*
RSB511	6	5	2.5	0.038*
RSB611	6	6	2.5	0.038*
RSB711	6	7	2.5	0.038*
RSB811	6	8	2.5	0.038*
RSB911	6	9	2.5	0.038*
RSB412	10	4	2.5	0.035*
RSB512	10	5	2.5	0.038*
RSB612	10	6	2.5	0.038*
RSB712	10	7	2.5	0.038*
RSB812	10	8	2.5	0.038*
RSB912	10	9	2.5	0.038*

*45 cm mini stainless steel spring coil guidewire included.

Maximum guidewire compatible size: 0.035" for 4 Fr, 0.038" for all others.

Pinnacle® R/O II HiFlo

Radiopaque Marker Introducer Sheaths

Packed 10 per box and includes dilator.

Product Code	Length (cm)	French Size (Fr)	Short Guidewire
RSB604H	4	6	Spring 0.035"/40 cm*
RSB704H	4	7	Spring 0.035"/40 cm*
RSB804H	4	8	Spring 0.035"/40 cm*
RSB605H	4	6	–
RSB705H	4	7	–
RSB805H	4	8	–

Maximum guidewire compatible size is 0.038" *Guidewire included.

Glidesheath™

Hydrophilic Coated Introducer Sheaths

Packed 5 per box and includes a dilator and a spring wire.*

Product Code	Length (cm)	Sheath Inner Dia. (in)	Tip Inner Dia. (in)	French Size (Fr)	Guidewire Length (cm)
RM*BF4F10PA	10	0.061	0.055	4	45*
RM*BF5F10PA	10	0.074	0.068	5	45*
RM*BF6F10PA	10	0.087	0.082	6	45*
RM*BF6M10MA	10	0.087	0.082	6	45**
RM*BF4F16PA	16	0.061	0.055	4	80*
RM*BF5F16PA	16	0.074	0.068	5	80*
RM*BF6F16PA	16	0.087	0.082	6	80*
RM*BF6M16MA	16	0.087	0.082	6	80**
RM*BF4F25PA	25	0.061	0.055	4	80*
RM*BF5F25PA	25	0.074	0.068	5	80*
RM*BF6F25PA	25	0.087	0.082	6	80*
RM*BF6M25MA	25	0.087	0.082	6	80**

*Mini stainless steel spring coil guidewire (0.021") with straight tip.

**Mini stainless steel spring coil guidewire (0.035") with J-tip.

Glidesheath™ Access Kits

Packed 5 per box. Includes Glidesheath, SURFLO® IV catheter, scalpel, dilator, 2.5 mL syringe, nitinol/plastic mini guidewire.

Product Code	Length (cm)	Sheath Inner Dia. (in)	Tip Inner Dia. (in)	French Size (Fr)	I.V. Cath Size (G)	Guidewire Length (cm)
RM*AF4G10HAT	10	0.061	0.055	4	20	45
RM*AF5G10HAT	10	0.074	0.068	5	20	45
RM*AF6G10HAT	10	0.087	0.082	6	20	45
RM*AF4D10HAU	10	0.061	0.055	4	22	45
RM*AF5D10HAU	10	0.074	0.068	5	22	45
RM*AF6D10HAU	10	0.087	0.082	6	22	45

4

Pinnacle® Destination®

Guiding Sheaths

Packed 1 per box and includes dilator. All dilators are 0.038" wire compatible.*

Product Code	Length (cm)	French Size (Fr)	Coating Length (cm)	Valve Type	Tip Shape
54-54502	45	5	5	Cross-Cut	Hockey Stick
54-54503	45	5	5	Cross-Cut	Multi-Purpose
54-54504	45	5	5	Cross-Cut	RDC
54-54505	45	5	5	Cross-Cut	LIMA
54-54501	45	5	35	Cross-Cut	Straight
54-64501	45	6	35	Cross-Cut	Straight
54-64506	45	6	35	Tuohy-Borst	Straight
54-74501	45	7	35	Cross-Cut	Straight
54-74506	45	7	35	Tuohy-Borst	Straight
RSR01	45	6	5	Cross-Cut	Straight
RSR02	45	6	5	Cross-Cut	Hockey Stick
RSR03	45	6	5	Cross-Cut	Multi-Purpose
RSR04	45	7	5	Cross-Cut	Straight
RSR05	45	7	5	Cross-Cut	Hockey Stick
RSR06	45	7	5	Cross-Cut	Multi-Purpose
RSR07	45	6	5	Tuohy-Borst	Straight
RSR08	45	6	5	Tuohy-Borst	Hockey Stick
RSR09	45	6	5	Tuohy-Borst	Multi-Purpose
RSR10	45	7	5	Tuohy-Borst	Straight
RSR11	45	7	5	Tuohy-Borst	Hockey Stick
RSR12	45	7	5	Tuohy-Borst	Multi-Purpose
RSR13	45	6	5	Cross-Cut	RDC
RSR14	45	6	5	Cross-Cut	LIMA
RSR15	45	7	5	Cross-Cut	RDC
RSR16	45	7	5	Cross-Cut	LIMA
RSR17	45	6	5	Tuohy-Borst	RDC
RSR18	45	6	5	Tuohy-Borst	LIMA
RSR19	45	7	5	Tuohy-Borst	RDC
RSR20	45	7	5	Tuohy-Borst	LIMA
RSP01	65	6	35	Cross-Cut	Straight
RSP02	65	7	35	Cross-Cut	Straight
RSP03	65	6	35	Tuohy-Borst	Straight
RSP04	65	7	35	Tuohy-Borst	Straight
RSC01	90	6	15	Tuohy-Borst	Straight
RSC02	90	7	15	Tuohy-Borst	Straight
RSC03	90	6	15	Tuohy-Borst	Multi-Purpose
RSC04	90	7	15	Tuohy-Borst	Multi-Purpose
RSC05	90	6	15	Cross-Cut	Straight
RSC06	90	7	15	Cross-Cut	Straight
RSC07	90	6	15	Cross-Cut	Multi-Purpose
RSC08	90	7	15	Cross-Cut	Multi-Purpose

*Guidewire not included.

5

GlideAccess™ Systems**Micro Accessing Systems**

Packed 10 per box.

GlideAccess Systems with 21 Gauge Non-Echogenic Needle

Product Code	Coaxial Diameter (Fr)	Wire	Wire Length (cm)	Wire Diameter (in)
RSM03	4	Platinum-tipped Nitinol	40	0.018
RSM04	5	Platinum-tipped Nitinol	40	0.018

GlideAccess Systems with 21 Gauge Echogenic Needle

Product Code	Coaxial Diameter (Fr)	Wire	Wire Length (cm)	Wire Diameter (in)
RSM05	4	Stainless Steel	40	0.018
RSM06	5	Stainless Steel	40	0.018
RSM07	4	Platinum-tipped Nitinol	40	0.018
RSM08	5	Platinum-tipped Nitinol	40	0.018

TR Band™**Radial Compression Devices**

Packed 5 per box. Includes TR Band inflator.

Product Code	TR Band Type	Length (cm)
XX*RF06	Regular	24
XX*RF06L	Long	29

Peripheral Guidewires (0.032" - 0.038")

Glidewire®**Hydrophilic Coated Guidewires**

Packed 5 per box.

Standard Glidewire

Product Code	Wire Type	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
GR3202	Standard	0.032	150	3	Angle
GR3201	Standard	0.032	150	3	Straight
GR3505	Standard	0.035	120	3	Angle
GR3506	Standard	0.035	150	3	Angle
GR3508	Standard	0.035	180	3	Angle
GR3509	Standard	0.035	260	3	Angle
GR3501	Standard	0.035	150	3	Straight
GR3503	Standard	0.035	180	3	Straight
GR3504	Standard	0.035	260	3	Straight
GR3513	Standard	0.035	80	3	Angle
GR3805	Standard	0.038	120	3	Angle
GR3806	Standard	0.038	150	3	Angle
GR3808	Standard	0.038	180	3	Angle
GR3809	Standard	0.038	260	3	Angle
GR3811	Standard	0.038	80	3	Angle
GR3801	Standard	0.038	150	3	Straight
GR3803	Standard	0.038	180	3	Straight
GR3804	Standard	0.038	260	3	Straight

Shapeable Tip Glidewire

Product Code	Wire Type	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
GE3501	Shapeable	0.035	150	3	Straight/Shapeable
GE3502	Shapeable	0.035	180	3	Straight/Shapeable
GE3801	Shapeable	0.038	150	3	Straight/Shapeable

Long Taper Glidewire

Product Code	Wire Type	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
GR3507	Long Taper	0.035	150	8	Angle
GR3502	Long Taper	0.035	150	8	Straight
GR3807	Long Taper	0.038	150	8	Angle
GR3802	Long Taper	0.038	150	8	Straight

Peripheral Guidewires (0.035" - 0.038")

Glidewire®

Hydrophilic Coated Guidewires

Packed 5 per box.

Stiff Shaft Glidewire

Product Code	Wire Type	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
GS3505	Stiff	0.035	80	3	Angle
GS3506	Stiff	0.035	150	3	Angle
GS3508	Stiff	0.035	180	3	Angle
GS3509	Stiff	0.035	260	3	Angle
GS3501	Stiff	0.035	150	3	Straight
GS3503	Stiff	0.035	180	3	Straight
GS3504	Stiff	0.035	260	3	Straight
GS3804	Stiff	0.038	80	3	Angle
GS3805	Stiff	0.038	150	3	Angle
GS3807	Stiff	0.038	180	3	Angle
GS3801	Stiff	0.038	150	3	Straight
GS3803	Stiff	0.038	180	3	Straight

Stiff Shaft Long Taper Glidewire

GS3507	Stiff	0.035	150	5	Angle
GS3510	Stiff	0.035	260	5	Angle
GS3502	Stiff	0.035	150	5	Straight
GS3806	Stiff	0.038	150	5	Angle
GS3802	Stiff	0.038	150	5	Straight

1 cm Taper Glidewire

GR3515	Standard	0.035	150	1	Angle
GR3516	Standard	0.035	180	1	Angle
GR3517	Standard	0.035	260	1	Angle
GR3518	Standard	0.035	150	1	Straight
GR3519	Standard	0.035	180	1	Straight
GR3520	Standard	0.035	260	1	Straight
GS3511	Stiff	0.035	150	1	Angle

J-Tip Glidewire

GR3510	Standard	0.035	150	3	J-Tip
GR3810	Standard	0.038	150	3	J-Tip

Bolia Curve Glidewire

GR3521	Standard	0.035	180	5	Bolia
GR3522	Standard	0.035	180	8	Bolia

Small Vessel Guidewires (0.018" - 0.025")

Glidewire®

Hydrophilic Coated Guidewires

Packed 5 per box.

Standard Glidewire

Product Code	Wire Type	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Tip Shape
GR1804	Standard	0.018	80	3	Angle
GR1805	Standard	0.018	150	3	Angle
GR1806	Standard	0.018	180	3	Angle
GR1802	Standard	0.018	150	3	Straight
GR1803	Standard	0.018	180	3	Straight
GR1807	Standard	0.018	260	3	Angle
GR1808	Standard	0.018	260	3	Straight
GR2504	Standard	0.025	150	3	Angle
GR2505	Standard	0.025	180	3	Angle
GR2506	Standard	0.025	260	3	Angle
GR2501	Standard	0.025	150	3	Straight
GR2502	Standard	0.025	180	3	Straight
GR2503	Standard	0.025	260	3	Straight

Shapeable Tip Glidewire

GE1801	Shapeable	0.018	150	3	Straight/Shapeable
--------	-----------	-------	-----	---	--------------------

Glidewire® GT

Super-Selective Hydrophilic Coated Guidewires

Packed 1 per box and includes a torque device.

Product Code	Dia. (in)	Total Length (cm)	Grind Length (cm)	Radiopaque Length (mm)	Tip Shape
RG*GA1818SA	0.018	180	25	20	45° Angle ↙
RG*GE1818SA	0.018	180	25	20	90° Angle ↖
RG*GW1818SA	0.018	180	25	20	90°/60° Dble Angle† ↗
RG*ES1818SA	0.018	180	25	20	Straight/Shapeable

Note: Flexible Tip Type is standard for the above product codes.

†Distal/proximal.

Small Vessel Guidewires (0.018")

Glidewire® Gold

Hydrophilic Coated Guidewires

Packed 1 per box and includes a torque device.

Product Code	Dia. (in)	Total Length (cm)	Flexible Tip Length (cm)	Radiopaque Length (mm)	Tip Shape
GM1811	0.018	180	3	2	45° Angle
GM1813	0.018	180	3	2	70° Angle
GM1812	0.018	180	5	2	45° Angle
GM1814	0.018	180	5	2	70° Angle

Coronary Guidewire (0.014")

Runthrough® NS

Coronary Guidewires

Packed 5 per box and includes a re-shaper.

Product Code	Dia. (in)	Length (cm)	Tip Radiopacity (cm)	Tip Shape
25-1011	0.014	180	3	Straight/Shapeable
25-1013	0.014	300	3	Straight/Shapeable

Packed 1 per box. Does not include a re-shaper.

25-9092	0.014	150	–	Extension Wire
---------	-------	-----	---	----------------

Glidecath® 4 Fr

Hydrophilic Coated Catheters

Packed 5 per box. 0.038" guidewire compatible.

Product Code	Total Length (cm)	Flow Rate* (mL/sec)	Hydrophilic Coating Length (cm)	Tip Shape
CG412	65	15	15	Straight
CG415	65	15	15	Angle
CG408	65	15	15	C1
CG409	65	15	15	C2
CG411	65	15	15	J-Tip
CG422 †	70	13	45	Yashiro
CG413	100	13	15	Straight
CG416	100	13	15	Angle
CG410	100	13	15	C2
CG418	100	13	15	MP
CG401	100	13	15	Sim 1
CG402	100	13	15	Sim 2
CG403	100	13	15	Sim 3
CG404	100	13	15	H1
CG405	100	13	15	JB1
CG406	100	13	15	JB2
CG407	100	13	15	Mani
CG414	120	11	15	Straight
CG417	120	11	15	Angle

*All measurements taken at 750 psi and 37°C using Renografin® 76.

† CG422 is packed 3 per box.

Glidecath® XP 5 Fr

Hydrophilic Coated Catheters

Packed 5 per box. 0.038" guidewire compatible.
















Product Code	Total Length (cm)	Flow Rate* (mL/sec)	Hydrophilic Coating Length (cm)	Tip Shape
RF*ZM7510GM	65	24	15	Straight
RF*ZV9510GM	65	24	15	Angle
RF*ZB5510GM	65	24	15	C2
RF*ZM75110M	100	21	15	Straight
RF*ZV95110M	100	21	15	Angle
RF*WA15110M	100	21	15	Sim 1
RF*WA25110M	100	21	15	Sim 2
RF*WB15110M	100	21	15	H1

*All measurements taken at 1000 psi and 37°C using Renografin® 76.

Glidecath® 5 Fr

Hydrophilic Coated Catheters

Packed 5 per box. 0.038" guidewire compatible.

Product Code	Total Length (cm)	Flow Rate* (mL/sec)	Hydrophilic Coating Length (cm)	Tip Shape
CG505	65	21	25	Straight Taper 
CG507	65	21	25	Angle Taper 
CG501	65	21	25	C1 
CG502	65	21	25	C2 
CG504	65	21	25	J-Tip 
CG506	100	18	40	Straight Taper 
CG508	100	18	40	Angle Taper 
CG503	100	18	40	C2 
CG510	100	18	40	Sim 1 
CG511	100	18	40	Sim 2 
CG512	100	18	40	Sim 3 
CG513	100	18	40	H1 
CG514	100	18	40	JB1 
CG515	100	18	40	JB2 
CG516	100	18	40	Mani 
CG519 †	70	18	25	Yashiro 
CG509	100	18	40	MP 

* All measurements taken at 1000 psi and 37°C using Renografin® 76.

† CG519 is packed 3 per box.

Optitorque™

Coronary Diagnostic Catheters

Packed 5 per box.

Product Code by Size		Tip Shape	Length (cm)	Side Hole
5 Fr/1.70 mm	6 Fr/2.00 mm			
40-5011	40-6011	Radial Tiger 4.0	100	1
40-5012	40-6012	Radial Tiger 4.5	100	1
40-5021	40-6021	Jacky Radial 3.5	100	2
40-5022	40-6022	Sarah Radial 4.0	100	2
40-5030	40-6030	Straight Pigtail	110	6
40-5031	40-6031	Angled Pigtail 155°	110	6
40-5040	40-6040	Judkins Left 3.5 cm	100	—
40-5041	40-6041	Judkins Left 4.0 cm	100	—
40-5042	40-6042	Judkins Left 5.0 cm	100	—
40-5043	40-6043	Judkins Right 3.5 cm	100	—
40-5044	40-6044	Judkins Right 4.0 cm	100	—
40-5045	40-6045	Judkins Right 5.0 cm	100	—
40-5050	40-6050	Internal Mammary	100	—
40-5060	40-6060	Amplatz Left Small Type	100	—
40-5061	40-6061	Amplatz Left Middle Type	100	—

Note: Maximum Pressure: 1000 psi.

FineCross™ MG

Coronary Micro-Guide Catheters

Packed 1 per box.

Product Code	Catheter Length	Distal Outer Diameter	Proximal Outer Diameter	Distal Inner Diameter	Proximal Inner Diameter	Guidewire* Size
35-1430	130 cm	1.8 Fr. (0.60 mm)	2.6 Fr. (0.87 mm)	0.018" (0.45 mm)	0.021" (0.55 mm)	0.014" (0.36 mm)
35-1450	150 cm					

*Guidewire not included.

Progreat™ Ω

- Flexible Tip: 17 cm
- Max GW: 0.021 in
- Max Pressure: 900 psi

Microcatheters with R/O Marker Band

Packed 1 per box and includes one hemostatic valve and one shaping mandrel.

Product Code	Total Length (cm)	O.D. D/P (Fr)	I.D. (in)	Coating Length (cm)	Flow Rate† (mL/sec)
MC*PB2811Y	110	2.8/3.0	0.027	50	4.7
MC*PB2813Y	130	2.8/3.0	0.027	70	4.3

Progreat™ Ω Microcatheters have a 0.7 mm platinum marker band.

†Flow rates taken at 900 psi and 37°C using Isovue® 300.

Progreat™

- Flexible Tip: 17 cm
- Max GW: 0.021 in
- Max Pressure: 750 psi

Microcatheters

Packed 1 per box and includes one hemostatic valve and one shaping mandrel.

Product Code	Total Length (cm)	O.D. D/P (Fr)	I.D. (in)	Coating Length (cm)	Flow Rate† (mL/sec)
MC*PC2711Y	110	2.7/2.9	0.025	50	3.5
MC*PC2713Y	130	2.7/2.9	0.025	70	3.2

†Flow rates taken at 750 psi and 37°C using Isovue 300.

Progreat™ 2.4 Fr

- Flexible Tip: 17 cm
- Max GW: 0.018 in
- Max Pressure: 750 psi

Microcatheters with R/O Marker Band

Packed 1 per box and includes one hemostatic valve and one shaping mandrel.

Product Code	Total Length (cm)	O.D. D/P (Fr)	I.D. (in)	Coating Length (cm)	Flow Rate† (mL/sec)
MC*PB2411Y	110	2.4/2.9	0.022	50	2.3
MC*PB2413Y	130	2.4/2.9	0.022	70	2.0

Progreat™ 2.4 Fr Microcatheters have a 0.7 mm platinum marker band.

†Flow rates taken at 750 psi and 37°C using Isovue 300.

Key to abbreviations:	O.D. = Outer Diameter
	I.D. = Inner Diameter
	D = Distal, P = Proximal

Progreat™ Ω

- Flexible Tip: 17 cm
- Max GW: 0.021 in
- Max Pressure: 900 psi

Coaxial Microcatheter Systems with R/O Marker Band

Each system includes a Progreat Microcatheter, pre-loaded Glidewire® GT 0.021" guidewire, one wire stopper, one guidewire introducer, one hemostatic valve, one 2.5 mL syringe and one shaping mandrel.

Product Code	Total Length (cm)	O.D. D/P (Fr)	I.D. (in)	Coating Length (cm)	Flow Rate† (mL/sec)
MC*PE28111YB	110	2.8/3.0	0.027	50	4.7
MC*PE28131YB	130	2.8/3.0	0.027	70	4.3

Progreat Coaxial Microcatheters have a 0.7 mm platinum marker band.

†Flow rates taken at 900 psi and 37°C using Isovue 300.

Progreat™

- Flexible Tip: 17 cm
- Max GW: 0.021 in
- Max Pressure: 750 psi

Coaxial Microcatheter Systems

Each system includes a Progreat Microcatheter, pre-loaded Glidewire® GT 0.021" guidewire, one wire stopper, one guidewire introducer, one hemostatic valve, one 2.5 mL syringe and one shaping mandrel.

Product Code	Total Length (cm)	O.D. D/P (Fr)	I.D. (in)	Coating Length (cm)	Flow Rate† (mL/sec)
MC*PE27111Y	110	2.7/2.9	0.025	50	3.5
MC*PE27131Y	130	2.7/2.9	0.025	70	3.2

†Flow rates taken at 750 psi and 37°C using Isovue 300.

Key to abbreviations:	O.D. = Outer Diameter
	I.D. = Inner Diameter
	D = Distal, P = Proximal

AZUR™

Peripheral HydroCoil® Embolization Systems**Detachable 18 System.** Hydrogel Coil. Packed 1 per box.

Product Code	Loop Diameter (mm)	Length* (cm)
45-480305	3	5
45-480405	4	5
45-480310	3	10
45-480410	4	10
45-480610	6	10
45-480810	8	10
45-481010	10	10
45-480415	4	15
45-480615	6	15
45-480815	8	15
45-481015	10	15
45-481215	12	15
45-481515	15	15
45-480420	4	20
45-480620	6	20
45-480820	8	20
45-481020	10	20
45-481220	12	20
45-481520	15	20
45-482020	20	20
45-481530	15	30
45-482030	20	30

Detachment Controller for use with Detachable Systems

Packed 5 per box.

Product Code	Product Description
45-4001	AZUR Detachment Controller

Pushable 18 System. Hydrogel Coil. Packed 3 per box.

Product Code	Loop Diameter (mm)	Length* (cm)
45-280202	2	2
45-280302	3	2
45-280402	4	2
45-280304	3	4
45-280404	4	4
45-280504	5	4
45-280406	4	6
45-280506	5	6
45-280606	6	6
45-280510	5	10
45-280610	6	10
45-280810	8	10
45-280614	6	14
45-280814	8	14
45-281014	10	14

16

*When the coil is straight, length is calculated from the tip of the coil to the point of connection with the pusher wire.

AZUR™

Pushable 35 System. Hydrogel Coil. Packed 3 per box.

Product Code	Loop Diameter (mm)	Length* (cm)
45-250404	4	4
45-250406	4	6
45-250506	5	6
45-250510	5	10
45-250610	6	10
45-250810	8	10
45-250814	8	14
45-251014	10	14
45-251514	15	14

Peripheral Framing Coil Embolization Systems**Detachable 18 System.** Non-hydrogel coil. Packed 1 per box.

Product Code	Loop Diameter (mm)	Length* (cm)
45-680820	8	20
45-681026	10	26
45-681434	14	34
45-682050	20	50

Note: Use in conjunction with AZUR Detachment Controller.

BeadBlock™
**Compressible Microspheres** Packed 1 per box. 2 mL 1 Pack.

Product Code	Size Range (µm)	Color Code
EB2S103	100 – 300	Yellow
EB2S305	300 – 500	Blue
EB2S507	500 – 700	Red
EB2S709	700 – 900	Green
EB2S912	900 – 1200	Purple

20 mL syringe contains 2 mL volume of Bead Block Microspheres suspended in physiological buffered saline.

*When the coil is straight, length is calculated from the tip of the coil to the point of connection with the pusher wire.

17

**For more information, please call 800.862.4143 to
speak to an Inside Sales Customer Care Representative,
or visit www.terumo.is.com**

Terumo Interventional Systems

2101 Cottontail Lane • Somerset, NJ 08873 • Fax: 800.411.5870

© 2009 Terumo Medical Corporation. All rights reserved. HydroCoil is a registered trademark of MicroVention Inc. Bead Block is a trademark of Biocompatibles UK Ltd. Renografin and Isovue are registered trademarks of Bracco Diagnostics Inc. All other brand names are trademarks or registered trademarks of Terumo. TIS 400 6/09



TERUMO INTERVENTIONAL SYSTEMS