

# DIMENSION REFERENCE | 18 System

	Product Code	Product Description	Loop Diameter (mm)	Length (cm)	Number of Loops*	OD (in)*	Coin Stack Length (mm)*	Volume (mm3)*	
FRAMING 18	45-680820	AZUR 18D Framing 8mm x 20cm	8	20	6	0.014		19.86	
	45-681026	AZUR 18D Framing 10mm x 26cm	10	26	6	0.014		25.82	
	45-681434	AZUR 18D Framing 14mm x 34cm	14	34	6	0.015		38.76	
	45-682050	AZUR 18D Framing 20mm x 50cm	20	50	6	0.015		57.00	
DETACHABLE 18	45-480202	AZUR 18D System 2mm x 2cm	2	2	3.18	0.034	2.75	11.72	
	45-480204	AZUR 18D System 2mm x 4cm	2	4	6.36	0.034	5.49	23.43	
	45-480302	AZUR 18D System 3mm x 2cm	3	2	2.12	0.034	1.83	11.72	
	45-480305	AZUR 18D System 3mm x 5cm	3	5	5.31	0.034	4.59	29.29	
	45-480310	AZUR 18D System 3mm x 10cm	3	10	10.61	0.034	9.16	58.58	
	45-480405	AZUR 18D System 4mm x 5cm	4	5	3.98	0.034	3.44	29.29	
	45-480410	AZUR 18D System 4mm x 10cm	4	10	7.96	0.034	6.87	58.58	
	45-480415	AZUR 18D System 4mm x 15cm	4	15	11.94	0.034	10.31	87.86	
	45-480420	AZUR 18D System 4mm x 20cm	4	20	15.92	0.034	13.75	117.15	
	45-480505	AZUR 18D System 5 mm x 5cm	5	5	3.18	0.034	2.75	29.29	
	45-480510	AZUR 18D System 5mm x 10cm	5	10	6.37	0.034	5.5	58.58	
	45-480515	AZUR 18D System 5mm x 15cm	5	15	9.55	0.034	8.25	87.86	
	45-480520	AZUR 18D System 5mm x 20cm	5	20	12.73	0.034	11.00	117.15	
	45-480610	AZUR 18D System 6mm x 10cm	6	10	5.31	0.034	4.59	58.58	
	45-480615	AZUR 18D System 6mm x 15cm	6	15	7.96	0.034	6.87	87.86	
	45-480620	AZUR 18D System 6mm x 20cm	6	20	10.61	0.034	9.16	117.15	
	45-480810	AZUR 18D System 8mm x 10cm	8	10	3.98	0.034	3.44	58.58	
	45-480815	AZUR 18D System 8mm x 15cm	8	15	5.97	0.034	5.16	87.86	
	45-480820	AZUR 18D System 8mm x 20cm	8	20	7.96	0.034	6.87	117.15	
	45-481010	AZUR 18D System 10mm x 10cm	10	10	3.18	0.034	2.75	58.58	
	45-481015	AZUR 18D System 10mm x 15cm	10	15	4.77	0.034	4.12	87.86	
	45-481020	AZUR 18D System 10mm x 20cm	10	20	6.37	0.034	5.50	117.15	
	45-481215	AZUR 18D System 12mm x 15cm	12	15	3.98	0.034	3.44	87.86	
	45-481220	AZUR 18D System 12mm x 20cm	12	20	5.31	0.034	4.59	117.15	
	45-481515	AZUR 18D System 15mm x 15cm	15	15	3.18	0.034	2.75	87.86	
	45-481520	AZUR 18D System 15mm x 20cm	15	20	4.24	0.034	3.66	117.15	
	45-481530	AZUR 18D System 15mm x 30cm	15	30	6.37	0.034	5.50	175.73	
	45-482020	AZUR 18D System 20mm x 20cm	20	20	3.18	0.034	2.75	117.15	
	45-482030	AZUR 18D System 20mm x 30cm	20	30	4.77	0.034	4.12	175.73	
	CROSS SECTIONAL OCCLUDER 18	45-780202	AZUR CX 18D System 2mm x 2cm	2	2	3	0.015		2.28
		45-780204	AZUR CX 18D System 2mm x 4cm	2	4	6	0.015		4.56
		45-780304	AZUR CX 18D System 3mm x 4cm	3	4	4	0.015		4.56
		45-780308	AZUR CX 18D System 3mm x 8cm	3	8	8	0.015		9.12
45-780413		AZUR CX 18D System 4mm x 13cm	4	13	9	0.015		14.82	
45-780516		AZUR CX 18D System 5mm x 16cm	5	16	9	0.015		18.24	
45-780620		AZUR CX 18D System 6mm x 20cm	6	20	9	0.015		22.80	
45-780724		AZUR CX 18D System 7mm x 24cm	7	24	9	0.015		23.84	
45-780828		AZUR CX 18D System 8mm x 28cm	8	28	9	0.015		31.92	
45-780928		AZUR CX 18D System 9mm x 28cm	9	28	8	0.015		31.92	
45-781032		AZUR CX 18D System 10mm x 32cm	10	32	8	0.015		36.48	
45-781238		AZUR CX 18D System 12mm x 38cm	12	38	8	0.015		43.32	
45-781434		AZUR CX 18D System 14mm x 34cm	14	34	6	0.015		38.76	
45-781639		AZUR CX 18D System 16mm x 39cm	16	39	6	0.015		44.46	
45-781836		AZUR CX 18D System 18mm x 36cm	18	36	5	0.015		41.04	
45-782040		AZUR CX 18D System 20mm x 40cm	20	40	5	0.015		45.60	
PUSHABLE 18	45-280202	AZUR 18P System 2mm x 2cm	2	2	3.18	0.034	2.75	11.72	
	45-280302	AZUR 18P System 3mm x 2cm	3	2	2.12	0.034	1.83	11.72	
	45-280304	AZUR 18P System 3mm x 4cm	3	4	4.24	0.034	3.66	23.43	
	45-280402	AZUR 18P System 4mm x 2cm	4	2	1.59	0.034	1.37	11.72	
	45-280404	AZUR 18P System 4mm x 4cm	4	4	3.18	0.034	2.75	23.43	
	45-280406	AZUR 18P System 4mm x 6cm	4	6	4.77	0.034	4.12	35.15	
	45-280504	AZUR 18P System 5mm x 4cm	5	4	2.55	0.034	2.20	23.43	
	45-280506	AZUR 18P System 5mm x 6cm	5	6	3.82	0.034	3.30	35.15	
	45-280510	AZUR 18P System 5mm x 10cm	5	10	6.37	0.034	5.50	58.58	
	45-280514	AZUR 18P System 5mm x 14cm	5	14	8.91	0.034	7.70	82.01	
	45-280606	AZUR 18P System 6mm x 6cm	6	6	3.18	0.034	2.75	35.15	
	45-280610	AZUR 18P System 6mm x 10cm	6	10	5.31	0.034	4.59	58.58	
	45-280614	AZUR 18P System 6mm x 14cm	6	14	7.43	0.034	6.42	82.01	
	45-280620	AZUR 18P System 6mm x 20cm	6	20	10.61	0.034	9.16	117.15	
	45-280810	AZUR 18P System 8mm x 10cm	8	10	3.98	0.034	3.44	58.58	
	45-280814	AZUR 18P System 8mm x 14cm	8	14	5.57	0.034	4.81	82.01	
	45-280820	AZUR 18P System 8mm x 20cm	8	20	7.96	0.034	6.87	117.15	
	45-281014	AZUR 18P System 10mm x 14cm	10	14	4.46	0.034	3.85	82.01	
	45-281020	AZUR 18P System 10mm x 20cm	10	20	6.37	0.034	5.50	117.15	

# AZUR® Peripheral Embolization System

## DIMENSION REFERENCE | 35 System

# AZUR®

Peripheral HydroCoil Embolization System

	Product Code	Product Description	Loop Diameter (mm)	Length (cm)	Number of Loops*	OD (in)*	Coin Stack Length (mm)*	Volume (mm3)*
FRAMING 35	45-650820	AZUR 35D Framing 8mm x 20cm	8	20	6	0.022		49.05
	45-651026	AZUR 35D Framing 10mm x 26cm	10	26	6	0.022		63.76
	45-651434	AZUR 35D Framing 14mm x 34cm	14	34	6	0.022		83.38
	45-652050	AZUR 35D Framing 20mm x 50cm	20	50	6	0.022		122.62
DETACHABLE 35	45-450405	AZUR 35D System 4mm x 5cm	4	5	3.98	0.048	4.85	58.37
	45-450410	AZUR 35D System 4mm x 10cm	4	10	7.96	0.048	9.71	116.75
	45-450415	AZUR 35D System 4mm x 15cm	4	15	11.94	0.048	14.56	175.12
	45-450610	AZUR 35D System 6mm x 10cm	6	10	5.31	0.048	6.47	116.75
	45-450615	AZUR 35D System 6mm x 15cm	6	15	7.96	0.048	9.71	175.12
	45-450620	AZUR 35D System 6mm x 20cm	6	20	10.62	0.048	12.94	233.49
	45-450815	AZUR 35D System 8mm x 15cm	8	15	5.97	0.048	7.28	175.12
	45-450820	AZUR 35D System 8mm x 20cm	8	20	7.96	0.048	9.71	233.49
	45-451015	AZUR 35D System 10mm x 15cm	10	15	4.78	0.048	5.82	175.12
	45-451020	AZUR 35D System 10mm x 20cm	10	20	6.37	0.048	7.77	233.49
	45-451215	AZUR 35D System 12mm x 15cm	12	15	3.98	0.048	4.85	175.12
	45-451220	AZUR 35D System 12mm x 20cm	12	20	5.31	0.048	6.47	233.49
	45-451230	AZUR 35D System 12mm x 30cm	12	30	7.96	0.048	9.71	350.24
	45-451520	AZUR 35D System 15mm x 20cm	15	20	4.25	0.048	5.18	233.49
	45-451530	AZUR 35D System 15mm x 30cm	15	30	6.37	0.048	7.77	350.24
	45-452020	AZUR 35D System 20mm x 20cm	20	20	3.18	0.048	3.88	233.49
45-452030	AZUR 35D System 20mm x 30cm	20	30	4.77	0.048	5.82	350.24	
CROSS SECTIONAL OCCLUDER 35	45-750407	AZUR CX Detachable 035 4mm x 7cm	4	7	6	0.029		29.83
	45-750511	AZUR CX Detachable 035 5mm x 11cm	5	11	7	0.029		46.88
	45-750609	AZUR CX Detachable 035 6mm x 9cm	6	9	4	0.029		38.35
	45-750617	AZUR CX Detachable 035 6mm x 17cm	6	17	8	0.029		72.44
	45-750812	AZUR CX Detachable 035 8mm x 12cm	8	12	4	0.029		51.14
	45-750824	AZUR CX Detachable 035 8mm x 24cm	8	24	8	0.029		102.27
	45-751019	AZUR CX Detachable 035 10mm x 19cm	10	19	5	0.029		80.97
	45-751324	AZUR CX Detachable 035 13mm x 24cm	13	24	5	0.029		102.27
	45-751632	AZUR CX Detachable 035 16mm x 32cm	16	32	5	0.029		136.37
45-752039	AZUR CX Detachable 035 20mm x 39cm	20	39	5	0.029		166.19	
PUSHABLE 35	45-250404	AZUR 35P System 4mm x 4cm	4	4	3.18	0.048	3.88	46.70
	45-250406	AZUR 35P System 4mm x 6cm	4	6	4.77	0.048	5.82	70.05
	45-250506	AZUR 35P System 5mm x 6cm	5	6	3.82	0.048	4.66	70.05
	45-250510	AZUR 35P System 5mm x 10cm	5	10	6.37	0.048	7.77	116.75
	45-250610	AZUR 35P System 6mm x 10cm	6	10	5.31	0.048	6.47	116.75
	45-250614	AZUR 35P System 6mm x 14cm	6	14	7.43	0.048	9.06	163.44
	45-250810	AZUR 35P System 8mm x 10cm	8	10	3.98	0.048	4.85	116.75
	45-250814	AZUR 35P System 8mm x 14cm	8	14	5.57	0.048	6.79	163.44
	45-250820	AZUR 35P System 8mm x 20cm	8	20	7.96	0.048	9.70	233.49
	45-251014	AZUR 35P System 10mm x 14cm	10	14	4.46	0.048	5.44	163.44
	45-251020	AZUR 35P System 10mm x 20cm	10	20	6.37	0.048	7.76	233.49
	45-251514	AZUR 35P System 15mm x 14cm	15	14	2.97	0.048	3.62	163.44
	45-251520	AZUR 35P System 15mm x 20cm	15	20	4.24	0.048	5.17	233.49
	45-251620	AZUR 35P System 16mm x 20cm	16	20	3.98	0.048	4.85	233.49



### AZUR® Detachable 18 Framing Coil System:

Diameters: 8mm – 20mm  
Lengths: 20cm – 50cm

### AZUR® Detachable 18 System:

Diameters: 2mm – 20mm  
Lengths: 2cm – 40cm

### AZUR® Pushable 18 System:

Diameters: 2mm – 10mm  
Lengths: 2cm – 20cm

### AZUR® Detachable 35 Framing Coil System:

Diameters: 8mm – 20mm  
Lengths: 20cm – 50cm

### AZUR® Detachable 35 System:

Diameters: 4mm – 20mm  
Lengths: 5cm – 30cm

### AZUR® Pushable 35 System:

Diameters: 4mm – 16mm  
Lengths: 4cm – 20cm

**FIND OUT MORE** ☎ Phone: 800.862.4143 🌐 [terumo.com](http://terumo.com) 📠 Fax: 800.411.5870

\* Values are estimated based on engineering calculations and should only be used as reference. The actual measurements may vary based on multiple factors such as: blood flow, landing zone, sizing, packing density, technique, etc.

©2017 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. TIS-493-07272017/DF