


# GLIDEWIRE Tip Load Weight



- Straight tip rates are higher due to the nature of the test

Gram Tip Weights* (10 mm)	Grams
Glidewire 0.035" Reg. 3 cm, angled	6.8
Glidewire 0.035" Reg. 3 cm, straight	44.8
Glidewire 0.035" Reg. 1 cm, straight	49.9
Glidewire 0.035" Stiff 3 cm, angled	7.0
Glidewire 0.035" Stiff 3 cm, straight	74.9
Glidewire Gold 0.018" 3 cm 45 degree angle	4.1
Glidewire Advantage 0.035" 5 cm angled	11.1
Glidewire Advantage 0.018" 1 cm angled	2.6
Glidewire Advantage 0.014" 1 cm angled	2.1

\*Standard test method used