

Pinnacle®

Introducer Sheaths

PINNACLE Sheaths ID / OD and Tip Measurements

	Size (Fr)	Sheath OD (mm / in)	Tip ID (mm / in)	Sheath ID (mm / in)
PINNACLE R/O II Radiopaque Marker Introducer Sheath	4	2.11 / 0.083	1.39 / 0.055	1.55 / 0.061
	5	2.44 / 0.096	1.73 / 0.068	1.88 / 0.074
	6	2.77 / 0.109	2.06 / 0.081	2.21 / 0.087
	7	3.10 / 0.122	2.41 / 0.095	2.54 / 0.100
	8	3.48 / 0.137	2.72 / 0.107	2.92 / 0.115
	9	3.89 / 0.153	3.09 / 0.122	3.33 / 0.131

Mini stainless steel spring coil guidewire optional for 6cm and 10cm length sheaths.
Maximum guidewire compatible size: 0.035" for 4 Fr; 0.038" for all others.

	Size (Fr)	Sheath OD (mm / in)	Tip ID (mm / in)	Sheath ID (mm / in)
PINNACLE R/O II HiFlo Introducer Sheath	6	2.77 / 0.109	2.06 / 0.081	2.21 / 0.087
	7	3.10 / 0.122	2.41 / 0.095	2.54 / 0.100
	8	3.48 / 0.137	2.72 / 0.107	2.92 / 0.115

Optional mini stainless steel 0.035" guidewire.

	Size (Fr)	Sheath OD (mm / in)	Tip ID (mm / in)	Sheath ID (mm / in)
PINNACLE Introducer Sheath (Standard)	4	1.96 / 0.077	1.40 / 0.055	1.55 / 0.061
	5	2.29 / 0.090	1.73 / 0.068	1.88 / 0.074
	6	2.62 / 0.103	2.08 / 0.082	2.21 / 0.087
	7	2.95 / 0.116	2.41 / 0.095	2.54 / 0.100
	8	3.33 / 0.131	2.74 / 0.108	2.92 / 0.115
	9	3.73 / 0.147	3.07 / 0.121	3.33 / 0.131
	10	4.39 / 0.173	3.45 / 0.136	3.68 / 0.145
	11	4.83 / 0.190	3.78 / 0.149	4.14 / 0.163

0.035" or 0.038" mini stainless steel spring coil guidewire optional for 10 cm length.
Not available for 25cm length.

	Size (Fr)	Sheath OD (mm / in)	Tip ID (mm / in)	Sheath ID (mm / in)
PINNACLE TIF TIP Introducer Sheath	5	2.29 / 0.090	1.77 / 0.070	1.89 / 0.074
	6	2.62 / 0.103	2.10 / 0.083	2.22 / 0.087
	7	2.95 / 0.116	2.43 / 0.096	2.55 / 0.100
	8	3.32 / 0.131	2.79 / 0.110	2.92 / 0.115

Complete with 0.035" or 0.038" mini stainless steel spring coil guidewire.



Phone: 800.862.4143



www.terumo.com



Fax: 800.411.5870