

Glidewire®

Hydrophilic Coated Guidewire

MASTER YOUR CHOICE FOR PROCEDURE SUCCESS

The TERUMO GLIDEWIRE® product line that you know and trust offers an expanded range of tip shapes and wire configurations, making it possible for you to find the right tools for your challenging procedures.

1. Stiff GLIDEWIRE®

- Extra support and guidewire pushability designed to help straighten tortuous vasculature for better device delivery

2. 1 cm Taper GLIDEWIRE®

- Shorter taper designed to help cross stenotic lesions

3. Long Taper GLIDEWIRE®

- Long taper designed to both help maintain catheter tip shape and gain deeper vessel purchase to support guidewire advancement

4. J-Tip GLIDEWIRE®

- J-curve designed to provide atraumatic tip and avoid selection of vessel side branches



PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients**.



Consistent, reliable TERUMO GLIDEWIRE® performance designed to get you there

With TERUMO Glide Technology™ hydrophilic coating proven to have superior lubricity retention^{1-3*}, and a super-elastic nitinol alloy core for optimal performance and kink resistance⁴, the TERUMO GLIDEWIRE® offers reliable performance you know and trust.

	PRODUCT CODE	WIRE TYPE	DIAMETER	FLEXIBLE TIP LENGTH (TAPER)	TIP SHAPE	TOTAL LENGTH		
Stiff Shaft GLIDEWIRE®	GS3505	Stiff	0.035"	3 cm	Angle	80 cm		
	GS3506					150 cm		
	GS3508					180 cm		
	GS3509					260 cm		
	GS3501					150 cm		
	GS3503		180 cm					
	GS3504		260 cm					
	GS3804		0.038"		Angle	80 cm		
	GS3805					150 cm		
	GS3807					180 cm		
GS3801	150 cm							
GS3803	180 cm							
Long Taper GLIDEWIRE®	GS3507	Stiff	0.035"	5 cm	Angle	150 cm		
	GS3510					260 cm		
	GS3502				0.038"	Angle	150 cm	
	GS3806	150 cm						
	GS3802	150 cm						
	1 cm Taper GLIDEWIRE®	GR3507	Standard		0.035"	8 cm	Angle	150 cm
GR3502		150 cm						
GR3515		Standard	0.035"	1 cm			Angle	150 cm
GR3516								180 cm
GR3517								260 cm
GR3518							150 cm	
GR3519	Straight				180 cm			
GR3520		260 cm						
J-Tip GLIDEWIRE®	GR3510	Standard	0.035"	3 cm	3.0 mm J-tip	150 cm		
	GR3810		0.038"			180 cm		
	GR3525		0.035"		1.5 mm J-tip	180 cm		
	GR3526					260 cm		

FIND OUT MORE ☎ Phone: 800.862.4143 🌐 terumo.com 📠 Fax: 800.411.5870

*When compared to ZIPwire™, HiWire®, AQUATRACK®, and Laureate®.

References:

1. Shah A, Lau C, Stavropoulos SW, et al. Comparison of physician-rated performance characteristics of hydrophilic-coated guidewires. *J Vasc Interv Radiol.* 2008 Mar;19(3):400-405.
2. Niazi K, Farooqui F, Devireddy C, et al. Comparison of hydrophilic guidewires used in endovascular procedures. *J Invasive Cardiol.* 2009 Aug;21(8):397-400.
3. Niazi K. Guidewire Satisfaction Study. Data on File 2012, Terumo Medical Corporation.
4. Data on file.

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2021 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. PM-04701

