


**Navicross**<sup>®</sup>

Support Catheters

**MASTER YOUR  
PROCEDURES**  
WITH PRECISE  
PERFORMANCE



**NAVICROSS<sup>®</sup> 0.018"**  
Support Catheter

**TERUMO**  
INTERVENTIONAL  
SYSTEMS

# PUSHING BOUNDARIES


Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**



## NAVICROSS® 0.018"

PRODUCT CODE	WIRE COMPATIBILITY	CATHETER LENGTH	TIP SHAPE
NC18650	0.018"	65 cm	Straight
NC18651	0.018"	65 cm	Angle
NC18900	0.018"	90 cm	Straight
NC18901	0.018"	90 cm	Angle
NC18130	0.018"	135 cm	Straight
NC18131	0.018"	135 cm	Angle
NC18150	0.018"	150 cm	Straight
NC18151	0.018"	150 cm	Angle



### Double-braided stainless steel<sup>1</sup>

Engineered for optimal pushability and torque control for lesion crossing

### Telescoping Capability<sup>1</sup>



NAVICROSS® 0.018" is designed to telescope through NAVICROSS® 0.035"

### Straight and Angled Tips<sup>1</sup>

Tips allow access to vascular branches, including BTK collaterals

### Three Radiopaque Markers<sup>1</sup>

- Markers facilitate accurate assessment of position
- Unique spacing provides easy measurement of common stent and balloon sizes
- Initial marker is 1 mm from distal tip; 40 mm and 60 mm spacing

**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](https://www.terumo.com)

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

**Reference:**

1. Data on file.

©2024 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. PM-07798

 **TERUMO**  
INTERVENTIONAL  
SYSTEMS