

***Metacross***<sup>®</sup>


PTA Balloon Dilatation Catheter

***Crosstella***<sup>®</sup>

PTA Balloon Dilatation Catheter

***Croserio***<sup>®</sup>

PTA Balloon Dilatation Catheter



MAXIMIZE  
YOUR OPTIONS  
FOR **RX AND OTW**  
**PERFORMANCE**

 **TERUMO**  
INTERVENTIONAL  
SYSTEMS



# PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**

**TERUMO**  
INTERVENTIONAL  
SYSTEMS

# CHOOSE FROM A FULL LINE OF RX AND OTW OPTIONS

## Designed for optimal performance

**METACROSS®** PTA Balloon Dilatation Catheter (0.035")

**CROSSTELLA®** PTA Balloon Dilatation Catheter (0.018")

**CROSPERIO®** PTA Balloon Dilatation Catheter (0.014")

Providing end-to-end performance with less effort and more options, our rapid exchange (RX) and over-the-wire (OTW) PTA balloon dilatation catheters let you choose the optimal solution for you and your patients.

- 3 guidewire sizes available (0.014", 0.018", and 0.035")<sup>1-6</sup>
- Smaller balloon profiles help ease crossability<sup>7,a</sup>
  - Supports re-crossability within the same procedure<sup>1-6</sup>
- Well-designed shaft and stiff tapered core wire engineered for maximizing pushability<sup>b</sup>
  - Effective transition of force from proximal to distal end

# RAPID EXCHANGE: BETTER HANDLING BY DESIGN

## Achieve optimal results with Rapid Exchange

- Enables faster, more convenient exchange of devices
- Supports the use of shorter guidewires
- Encourages procedural efficiency, reducing fluoroscopy time and radiation dose
- Optimizes control, allowing a single operator to mitigate lost wire position



Balloon Catheter	Guidewire Compatibility (in)	Usable Shaft Length (cm)	Shaft Diameter		Proximal Tapered	Proximal Marker (mm)	Center (mm)	Distal Marker (mm)	Distal Tapered (mm)	Distal (mm)
<b>Metacross® RX</b>	0.035	135	Proximal 3.9 Fr (1.30 mm)	Distal 5.7 Fr (1.90 mm)	1.7	1.7	1.6	1.7	1.4	1.02
<b>Metacross® OTW</b>	0.035	65/90/135	5 Fr compatible: 5.4 Fr (1.80 mm) 6/7 Fr compatible: 5.7 Fr (1.90 mm)		1.7	1.7	1.6	1.7	1.4	1.07
<b>Crosstella® RX</b>	0.018	150	Proximal 3.4 Fr (1.12 mm) or 4.1 Fr (1.35 mm)	Distal 4 Fr (1.32 mm)	1.28	1.32	1.28	1.36	1.29	0.62
<b>Crosstella® OTW</b>	0.018	150	4 Fr compatible: 3.9 Fr (1.30 mm) 5 Fr compatible: 4.3 Fr (1.43 mm)		1.19	1.30	1.26	1.30	1.01	0.58
<b>Crosperio® RX</b>	0.014	150	Proximal 3.5 Fr (1.16 mm)	Distal 3.5 Fr (1.16 mm)	0.85	0.85	0.82	0.86	0.83	0.47
<b>Crosperio® OTW</b>	0.014	65/150	3.8 Fr (1.28 mm)		0.90	0.87	0.81	0.83	0.76	0.46

\* Diameter range:  
METACROSS® RX/OTW Balloon Catheter (3-12 mm);  
CROSSTELLA® RX/OTW Balloon Catheter (2-6 mm);  
CROSPERIO® RX/OTW Balloon Catheter (1.5-4 mm).

<sup>b</sup> RX design feature only.

<sup>c</sup> METACROSS® RX/OTW Balloon Catheter illustrations shown.



# Metacross® RX


PTA Balloon Dilatation Catheter

## ENGINEERED FOR PERFORMANCE

### Excellent crossability with innovative catheter design<sup>7</sup>

- Designed to provide best-in-class pushability<sup>7</sup>
- Innovative shaft provides robust strength constructed to cross the lesion<sup>7</sup>
- Smallest balloon profile among tested semi-compliant balloons<sup>7</sup>

### METACROSS RX Balloon Catheter has excellent crossability compared to select competitors



Sample size: METACROSS RX Balloon Catheter n=10, competitors n=5; size used for comparison: 6 mm x 40 mm.

METACROSS® RX

CROSSTELLA® RX

CROSPERIO® RX

See inserts for ordering information.

# Metacross®

PTA Balloon Dilatation Catheter

# Crosstella®

PTA Balloon Dilatation Catheter

# Crosperio®

PTA Balloon Dilatation Catheter



**FIND OUT MORE** 📞 Phone: 800.888.3786 🌐 [terumo.com](http://terumo.com)

#### Indications

CROSPERIO® RX PTA Balloon Dilatation Catheter, CROSPERIO® OTW PTA Balloon Dilatation Catheter, CROSSTELLA® RX PTA Balloon Dilatation Catheter, CROSSTELLA® OTW PTA Balloon Dilatation Catheter, METACROSS® RX PTA Balloon Dilatation Catheter, and METACROSS® OTW PTA Balloon Dilatation Catheter are intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal, and renal arteries, and for the treatment of obstructive lesions of native and synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

**RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.**

#### References:

1. METACROSS® RX PTA Balloon Dilatation Catheter IFU, September 2021.
2. METACROSS® OTW PTA Balloon Dilatation Catheter IFU, September 2021.
3. CROSSTELLA® RX PTA Balloon Dilatation Catheter IFU, April 2021.
4. CROSSTELLA® OTW PTA Balloon Dilatation Catheter IFU, September 2021.
5. CROSPERIO® RX PTA Balloon Dilatation Catheter IFU, September 2021.
6. CROSPERIO® OTW PTA Balloon Dilatation Catheter IFU, April 2021.
7. Data on file.

©2023 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-06596

 **TERUMO**  
INTERVENTIONAL  
SYSTEMS

# Metacross® RX

PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

### METACROSS® RX 0.035" PTA Balloon Dilatation Catheter - 135 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
3	BD-L30020CR	BD-L30040CR	BD-L30060CR	BD-L30080CR	BD-L30100CR	BD-L30120CR	BD-L30150CR	BD-L30200CR	6Fr
4	BD-L40020CR	BD-L40040CR	BD-L40060CR	BD-L40080CR	BD-L40100CR	BD-L40120CR	BD-L40150CR	BD-L40200CR	
5	BD-L50020CR	BD-L50040CR	BD-L50060CR	BD-L50080CR	BD-L50100CR	BD-L50120CR	BD-L50150CR	BD-L50200CR	
6	BD-L60020CR	BD-L60040CR	BD-L60060CR	BD-L60080CR	BD-L60100CR	BD-L60120CR	BD-L60150CR	BD-L60200CR	
7	BD-L70020CR	BD-L70040CR	BD-L70060CR	BD-L70080CR	BD-L70100CR	BD-L70120CR	BD-L70150CR	BD-L70200CR	
8	BD-L80020CR	BD-L80040CR	BD-L80060CR	BD-L80080CR	—	—	—	—	
9	BD-L90020CR	BD-L90040CR	BD-L90060CR	BD-L90080CR	—	—	—	—	7Fr
10	BD-LX0020CR	BD-LX0040CR	BD-LX0060CR	BD-LX0080CR	—	—	—	—	
12	BD-LX2020CR	BD-LX2040CR	—	—	—	—	—	—	



# Metacross® OTW

PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 65 cm									
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
3	—	BD-L30040SW	BD-L30060SW	—	—	—	—	—	5Fr
4	BD-L40020SW	BD-L40040SW	—	—	—	—	—	—	
5	—	BD-L50040SW	BD-L50060SW	—	—	—	BD-L50150SW	—	6Fr
6	BD-L60020SW	BD-L60040SW	BD-L60060SW	—	—	—	BD-L60150SW	—	
7	—	BD-L70040SW	BD-L70060SW	BD-L70080SW	—	—	—	—	7Fr
8	BD-L80020SW	BD-L80040SW	BD-L80060SW	BD-L80080SW	—	—	—	—	
9	BD-L90020SW	BD-L90040SW	BD-L90060SW	—	—	—	—	—	7Fr
10	—	BD-LX0040SW	BD-LX0060SW	BD-LX0080SW	—	—	—	—	
12	—	BD-LX2040SW	—	—	—	—	—	—	

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 90 cm									
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
3	BD-L30020MW	BD-L30040MW	BD-L30060MW	BD-L30080MW	BD-L30100MW	BD-L30120MW	BD-L30150MW	BD-L30200MW	5Fr
4	BD-L40020MW	BD-L40040MW	BD-L40060MW	BD-L40080MW	BD-L40100MW	BD-L40120MW	BD-L40150MW	BD-L40200MW	
5	BD-L50020MW	BD-L50040MW	BD-L50060MW	BD-L50080MW	BD-L50100MW	BD-L50120MW	BD-L50150MW	BD-L50200MW	6Fr
6	BD-L60020MW	BD-L60040MW	BD-L60060MW	BD-L60080MW	BD-L60100MW	BD-L60120MW	BD-L60150MW	BD-L60200MW	
7	BD-L70020MW	BD-L70040MW	BD-L70060MW	BD-L70080MW	BD-L70100MW	BD-L70120MW	BD-L70150MW	BD-L70200MW	7Fr
8	BD-L80020MW	BD-L80040MW	BD-L80060MW	BD-L80080MW	—	—	—	—	
9	BD-L90020MW	BD-L90040MW	BD-L90060MW	BD-L90080MW	—	—	—	—	7Fr
10	BD-LX0020MW	BD-LX0040MW	BD-LX0060MW	BD-LX0080MW	—	—	—	—	
12	BD-LX2020MW	BD-LX2040MW	—	—	—	—	—	—	

METACROSS® OTW 0.035" PTA Balloon Dilatation Catheter – 135 cm									
BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
3	BD-L30020CW	BD-L30040CW	BD-L30060CW	BD-L30080CW	BD-L30100CW	BD-L30120CW	BD-L30150CW	BD-L30200CW	5Fr
4	BD-L40020CW	BD-L40040CW	BD-L40060CW	BD-L40080CW	BD-L40100CW	BD-L40120CW	BD-L40150CW	BD-L40200CW	
5	BD-L50020CW	BD-L50040CW	BD-L50060CW	BD-L50080CW	BD-L50100CW	BD-L50120CW	BD-L50150CW	BD-L50200CW	6Fr
6	BD-L60020CW	BD-L60040CW	BD-L60060CW	BD-L60080CW	BD-L60100CW	BD-L60120CW	BD-L60150CW	BD-L60200CW	
7	BD-L70020CW	BD-L70040CW	BD-L70060CW	BD-L70080CW	BD-L70100CW	BD-L70120CW	BD-L70150CW	BD-L70200CW	7Fr
8	BD-L80020CW	BD-L80040CW	BD-L80060CW	BD-L80080CW	—	—	—	—	
9	BD-L90020CW	BD-L90040CW	BD-L90060CW	BD-L90080CW	—	—	—	—	7Fr
10	BD-LX0020CW	BD-LX0040CW	BD-LX0060CW	BD-LX0080CW	—	—	—	—	
12	BD-LX2020CW	BD-LX2040CW	—	—	—	—	—	—	

**FIND OUT MORE** ☎ Phone: 800.888.3786 🌐 [terumo.com](http://terumo.com)

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2023 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-06596




# Crosstella® RX

PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

### CROSSTELLA® RX 0.018" PTA Balloon Dilatation Catheter - 150 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
2	BD-F20020LR	BD-F20040LR	BD-F20060LR	BD-F20080LR	BD-F20100LR	BD-F20120LR	BD-F20150LR	BD-F20200LR	4Fr
2.5	BD-F25020LR	BD-F25040LR	BD-F25060LR	BD-F25080LR	BD-F25100LR	BD-F25120LR	BD-F25150LR	BD-F25200LR	
3	BD-F30020LR	BD-F30040LR	BD-F30060LR	BD-F30080LR	BD-F30100LR	BD-F30120LR	BD-F30150LR	BD-F30200LR	
4	BD-F40020LR	BD-F40040LR	BD-F40060LR	BD-F40080LR	BD-F40100LR	BD-F40120LR	BD-F40150LR	BD-F40200LR	
5	BD-F50020LR	BD-F50040LR	BD-F50060LR	BD-F50080LR	BD-F50100LR	BD-F50120LR	BD-F50150LR	BD-F50200LR	5Fr
6	BD-F60020LR	BD-F60040LR	BD-F60060LR	BD-F60080LR	BD-F60100LR	BD-F60120LR	BD-F60150LR	BD-F60200LR	





# Crosstella® OTW



PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

### CROSSTELLA® OTW 0.018" PTA Balloon Dilatation Catheter – 150 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)								
	20	40	60	80	100	120	150	200	
2	BD-F20020LW	BD-F20040LW	BD-F20060LW	BD-F20080LW	BD-F20100LW	BD-F20120LW	BD-F20150LW	BD-F20200LW	4Fr
2.5	BD-F25020LW	BD-F25040LW	BD-F25060LW	BD-F25080LW	BD-F25100LW	BD-F25120LW	BD-F25150LW	BD-F25200LW	
3	BD-F30020LW	BD-F30040LW	BD-F30060LW	BD-F30080LW	BD-F30100LW	BD-F30120LW	BD-F30150LW	BD-F30200LW	
4	BD-F40020LW	BD-F40040LW	BD-F40060LW	BD-F40080LW	BD-F40100LW	BD-F40120LW	BD-F40150LW	BD-F40200LW	
5	BD-F50020LW	BD-F50040LW	BD-F50060LW	BD-F50080LW	BD-F50100LW	BD-F50120LW	BD-F50150LW	BD-F50200LW	5Fr
6	BD-F60020LW	BD-F60040LW	BD-F60060LW	BD-F60080LW	BD-F60100LW	BD-F60120LW	BD-F60150LW	BD-F60200LW	



**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](http://terumo.com)

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2023 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-06596



# Crosperio® RX


PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

### CROSPERIO® RX 0.014" PTA Balloon Dilatation Catheter – 150 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	BD-B15020LR	BD-B15040LR	BD-B15080LR	BD-B15120LR	—	—
2	—	BD-B20040LR	BD-B20080LR	BD-B20120LR	BD-B20150LR	BD-B20200LR
2.5	—	BD-B25040LR	BD-B25080LR	BD-B25120LR	BD-B25150LR	BD-B25200LR
3	—	BD-B30040LR	BD-B30080LR	BD-B30120LR	BD-B30150LR	BD-B30200LR
3.5	—	BD-B35040LR	BD-B35080LR	BD-B35120LR	BD-B35150LR	BD-B35200LR
4	—	BD-B40040LR	BD-B40080LR	BD-B40120LR	BD-B40150LR	BD-B40200LR

4Fr



# Crosperio® OTW

PTA Balloon Dilatation Catheter

## ORDERING INFORMATION

### CROSPERIO® OTW 0.014" PTA Balloon Dilatation Catheter – 65 cm



BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	—	—	—	—	—	—
2	—	—	BD-B20080SW	—	—	BD-B20200SW
2.5	—	—	BD-B25080SW	—	—	BD-B25200SW
3	—	—	BD-B30080SW	—	—	BD-B30200SW
3.5	—	—	BD-B35080SW	—	—	BD-B35200SW
4	—	—	BD-B40080SW	—	—	BD-B40200SW

4Fr

### CROSPERIO® OTW 0.014" PTA Balloon Dilatation Catheter – 150 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)					
	20	40	80	120	150	200
1.5	BD-B15020LW	BD-B15040LW	BD-B15080LW	BD-B15120LW	—	—
2	—	BD-B20040LW	BD-B20080LW	BD-B20120LW	BD-B20150LW	BD-B20200LW
2.5	—	BD-B25040LW	BD-B25080LW	BD-B25120LW	BD-B25150LW	BD-B25200LW
3	—	BD-B30040LW	BD-B30080LW	BD-B30120LW	BD-B30150LW	BD-B30200LW
3.5	—	BD-B35040LW	BD-B35080LW	BD-B35120LW	BD-B35150LW	BD-B35200LW
4	—	BD-B40040LW	BD-B40080LW	BD-B40120LW	BD-B40150LW	BD-B40200LW

4Fr

**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](http://terumo.com)

**RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.**

©2023 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-06596

