


R2P Value Analysis Committee Kit



RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

Table of Contents:



| <u>Procedural and Economic Advantages</u> | <u>Product Portfolio Overview</u> | <u>Materials Management</u> |
|--|---|---|
| <ul style="list-style-type: none">• <u>What are Femoral Access Challenges?</u>• <u>When Radial Access is Considered</u>• <u>Paradigm Shift in Treatment of Peripheral Arterial Disease</u>• <u>R2P Procedural Advantages</u>• <u>Reimbursement Codes</u> | <ul style="list-style-type: none">• <u>GLIDESHEATH SLENDER Hydrophilic Coated Introducer Sheath</u>• <u>R2P DESTINATION SLENDER Guiding Sheath</u>• <u>R2P SLENGUIDE Guiding Catheter</u>• <u>GLIDEWIRE Hydrophilic Coated Guidewire</u>• <u>GLIDECATH Hydrophilic Coated Catheter</u>• <u>NAVICROSS Support Catheter</u>• <u>R2P METACROSS RX PTA Balloon Dilatation Catheter</u>• <u>R2P CROSSTELLA RX PTA Balloon Dilatation Catheter</u>• <u>R2P MISAGO RX Self-expanding Peripheral Stent</u>• <u>TR BAND Radial Compression Device</u> | <ul style="list-style-type: none">• <u>Ordering information and Brochure</u>• <u>510(k) clearance letter</u>• <u>Instructions for use</u><ul style="list-style-type: none">• <u>See IFU for Safety and storage</u>• <u>Service Support</u> |

Procedural and Economic Advantages

What are Femoral Access Challenges?

When Radial Access is Considered

Paradigm Shift in Treatment of Peripheral Arterial Disease

R2P Procedural Advantages

Reimbursement

What are the Femoral Access Challenges?

- ◆ Femoral approach compromised for peripheral procedures
 - ◆ Fem-Fem Bypass
 - ◆ Flush/proximal occlusion
 - ◆ Tortuous/Calcified Iliac
 - ◆ Previously deployed Grafts/Stents

- ◆ How has your physician tackled these patients previously?



Previous femoral-femoral bypass



Totally occluded left iliac,
limited access to right femoral



Significant iliac tortuosity
and calcification



Previous abdominal
aortic (AAA) graft



A **different approach** can support
your **patient challenges**

When Radial Access is Considered


- ◆ To overcome complex or challenging anatomy
- ◆ To mitigate or reduce access site complications
- ◆ To increase procedural and operational efficiencies
- ◆ Robust case planning & appropriate patient selection is essential to increase procedural success

Consider the clinical advantages of choosing **radial to peripheral** intervention¹⁻³

- Lower risk access site for antegrade interventions
- Can address both lower extremities from a single access point*
- Proven preference for transradial access versus alternative non-femoral access sites


Achieve procedural success with robust case planning support

PTA Balloon across lesion in distal SFA from Radial approach




Radial access can offer an *alternative solution* for treating Peripheral disease

Paradigm Shift in Treatment of Peripheral Arterial Disease



R2P Procedural Advantages

- ◆ Optimize entry site management with Slender Technology
- ◆ Enables the potential for **early ambulation with radial access**, improved patient comfort, and patient satisfaction
- ◆ Early ambulation provides potential cost reduction and **economic benefit** through various avenues.⁷
 - ◆ These avenues include expediting room turnover and increased throughput both through the procedural area and same day recovery unit.⁷



Reimbursement Codes (HCPCS Codes)

- ◆ GLIDESHEATH SLENDER
 - ◆ C1887: Guiding Catheter
- ◆ R2P DESTINATION SLENDER
 - ◆ C1887: Guiding Catheter
- ◆ R2P SLENGUIDE
 - ◆ C1887: Guiding Catheter
- ◆ R2P GLIDEWIRE
 - ◆ C1769: Guidewire
- ◆ NAVICROSS
 - ◆ C1887: Guiding Catheter
- ◆ GLIDECATH
 - ◆ C1887: Guiding Catheter
- ◆ R2P METACROSS RX
 - ◆ C1725: Catheter, Transluminal Angioplasty, Non-Laser
- ◆ R2P CROSSTELLA RX
 - ◆ C1725: Catheter, Transluminal Angioplasty, Non-Laser
- ◆ R2P MISAGO
 - ◆ C1876: Stent, Non-coated/Non-covered with Delivery System
- ◆ TR BAND
 - ◆ None; Does not fall into the “Pass Through” Categories of Devices

C-Code Disclaimer - Health economic and reimbursement information provided by Terumo Medical is gathered from third party sources and is subject to change without notice as a result of complex and frequently changing laws, regulations, rules and policies. This information is presented for full illustrative purposes only and does not constitute reimbursement or legal advice. Terumo Medical encourages providers to submit accurate and appropriate claims for services. It is always the provider's responsibility to determine medical necessity, the proper site for delivery of any services and to submit appropriate codes, charges, and modifiers for services that are rendered. Terumo Medical recommends that you consult with your payers, reimbursement specialists and/or legal counsel regarding coding, coverage and reimbursement matters. Terumo Medical does not promote the use of its products outside their FDA cleared or approved label.

Product Overview

Product Description
Features and Benefits

Product Description:

Glidesheath Slender[®]

Hydrophilic Coated Introducer Sheath



Indication: GLIDESHEATH SLENDER is indicated to facilitate placing a catheter through the skin into the radial artery.

- 7/6 Fr size
- 10 cm and 16 cm length
- Proprietary thin-wall technology and hydrophilic coating facilitate ease of insertion and removal during peripheral procedures

Features and Benefits:

- 1 *Terumo Glide Technology hydrophilic coating enables easy insertion and removal*
- 2 *Perform diagnostic and interventional procedures without the need to switch to a larger sheath*
- 3 *Ultra thin-wall technology*
- 4 *Low profile sheath to reduce arteriotomy size¹*

Product Description:

R2P *Destination Slender*[™]
Guiding Sheath



Indication: R2P DESTINATION Slender is intended to be used for the introduction of interventional and diagnostic devices in the lower extremities of the peripheral vasculature through an access site, including but not limited to the radial artery.

- 6 Fr size
- 119 cm and 149 cm length
- Full hydrophilic coating allows for smooth transition within the radial artery and Slender Technology provides optimal performance during peripheral procedures

Features and Benefits:

- 1 *Unique thin wall design*
- 2 *Provides a conduit for lower extremity interventions from a radial approach*
- 3 *Enables smaller arterial access while maintaining 6Fr compatibility*
- 4 *Designed to reduce potential for arterial spasm during introduction and removal¹*

Product Description:

R2P SlenGuide™
Guiding Catheter



Indication: The R2P (radial to peripheral) SlenGuide is designed for the introduction of interventional and diagnostic devices into the peripheral vasculature of the lower extremities.

- 7 Fr size
- 120 cm and 150 cm lengths
- Guiding catheter with distal hydrophilic coating and Slender Technology for optimal performance during peripheral procedures

Features and Benefits:

- 1 *Unique thin wall design*
- 2 *Enables smaller arterial access while maintaining 6Fr compatibility*
- 3 *Provides a conduit for lower extremity interventions from a radial approach*
- 4 *Designed to reduce potential for arterial spasm during introduction and removal¹*

Product Description:

Glidewire[®]

Hydrophilic Coated Guidewire



Indication: GLIDEWIRE Hydrophilic Coated Guidewires are designed to direct a catheter to the designed anatomical location during diagnostic or interventional procedures

- Straight, Angled, 1.5 mm and 3 mm JTip
- 350 cm, 400 cm and 450 cm lengths
- Includes multiple procedural options in shaft stiffness, tip shapes and wire configurations for easier vessel navigation

Features and Benefits:

- 1 *Terumo Glide Technology hydrophilic coating*
- 2 *Excellent crossability and trackability for challenging lesions as well as enhanced trackability through tortuous anatomy.*
- 3 *High durability with kink resistance and shape preservation as well as flexibility for easy artery navigation.*
- 4 *Excellent pushability for lesion access and notable torquability for procedure control.*

Product Description:

Glidecath[®]

Hydrophilic Coated Catheter



Indication: GLIDECATH is intended for use in angiographic procedures. It delivers radiopaque media and therapeutic agents to selected sites in the vascular system. It is also used to lead a guide wire or a catheter into the target site.

- 4 Fr size
- 150 cm length
- Offers excellent pushability and torque control with multiple tip shapes for selectivity and access to the peripheral vasculature

Features and Benefits:

- 1 *Provides excellent trackability and torque control⁸*
- 2 *Designed for optimal vessel selectivity and access to peripheral vasculature*
- 3 *Multiple tip shapes offered; Straight, Angle, PV Multicurve, Pigtail, Cobra 2, JR 4.0 & IMA*
- 4 *Double braided stainless steel for increased pushability*

Product Description:

Navicross[®]

Support Catheters



Indication: NAVICROSS is intended to guide and support a guide wire during access of the vasculature, allow for wire exchanges and provide a conduit for the delivery of saline or diagnostic contrast agents.

- Straight, Angled
- 135cm and 150cm lengths
- Double tapered tip provides the smallest crossing profile, giving a near seamless catheter-to-guidewire transition to aid in crossing simple or complex lesions in the peripheral vasculature

Features and Benefits:

- 1 *Best-in-class pushability*⁷
- 2 *Double tapered tip provides the smallest crossing profile*

- 3 *Double- braided stainless steel and three uniquely spaced radiopaque markers*
- 4 *Excellent wire-to-catheter transition enhancing crossing profile*

Product Description:

R2P™ *Metacross*® RX PTA Balloon Dilatation Catheter



Indication: The R2P METACROSS RX PTA balloon dilatation catheter is intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.

- 3-8 mm diameter x 20-200 mm length
- 200 cm shaft length
- The longest radial to peripheral capability and Rapid Exchange (RX) technology are designed to mitigate excessive device management and use of contrast media during a peripheral procedure

Features and Benefits:

- 1 *Tapered core wire from distal tube tip to proximal end of catheter*
- 2 *Unique technology stands apart from competition with the longest radial access to peripheral capability.²*
- 3 *Designed to mitigate excessive device management during procedures.²*
- 4 *Designed to mitigate amount of contrast media that may be used during the procedure.²*

Product Description:

R2P™ Crosstella® RX



Indication: The R2P CROSSTELLA RX PTA balloon dilatation catheter is intended to dilate stenosis in the iliac, femoral, ilio-femoral, popliteal, infra-popliteal and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for stent post-dilatation in the peripheral vasculature.


- 2-6 mm diameter x 40-200 mm length
- 200 cm shaft length
- The longest 0.018" Rapid Exchange (RX) PTA Balloon designed for use in radial to peripheral procedures

Features and Benefits:

- 1 *RX system is designed for ease of use by increasing procedural efficiency and reducing fluoroscopy time.²*
- 2 *CROSSTELLA RX PTA Balloon Dilatation Catheter provides the lowest deflation time results when tested against a competitor.²*
- 3 *Engineered to provide best-in-class pushability to help in crossing the lesion.²*
- 4 *CROSSTELLA RX PTA Balloon Dilatation Catheter provides excellent crossability when tested against competitor.²*

Product Description:

R2P™ Misago®
RX Self-expanding Peripheral Stent



Indication: The MISAGO RX self-expanding Peripheral Stent is indicated to improve luminal diameter in symptomatic patients with de novo or restenotic native lesions or occlusions of the superficial femoral artery and/or proximal popliteal artery with reference vessel diameters ranging from 4mm to 7mm and lesion length up to 150mm

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

- 6Fr sheath compatibility, 6-8 mm diameter x 40-150 mm length
- 200 cm shaft length with 6 gold reference markers
- Rapid Exchange design offers you optimized control and handling.

Features and Benefits:

- 1 Nitinol construction with a direct link design.
- 2 Rapid Exchange (RX) or "monorail" delivery system is designed for ease of use by increasing procedural efficiency and reducing fluoroscopy time.⁶
- 3 Minimizes jumping, recoil or foreshortening with a simplified thumbwheel system that allows for single operator deployment.⁶
- 4 The continuous spine-free stent is designed to promote optimal blood flow and eliminate high-strain stress zones that can lead to fracture.⁶

Product Description:

TR Band Radial Compression Device



Indication: TR Band is a compression device to assist hemostasis of the radial artery after a transradial procedure.

- Regular 24 cm and Long 29 cm
- The #1 preferred radial hemostasis device⁵ on the market provides a more precise way of applying pressure to the radial artery

Features and Benefits:

- 1 *Precise compression to achieve radial hemostasis*
- 2 *Air titration provides a more precise way of applying pressure to the radial artery⁴*
- 3 *Dual balloon technology provides precise compression of the radial artery without compromising local nerve structure³*
- 4 *Hemostasis achieved at low pressures, minimizing the chances of applying occlusive force*

[BACK](#)

Materials Management


Ordering Information and Brochure

510(k) Clearance Letter and PMA Approval Letter

Instructions for Use


Service and Support

Ordering Information Form and Brochure


Radial to Peripheral ORDERING FORM


| AMT | PRODUCT CODE | DESCRIPTION | AMT | PRODUCT CODE | DESCRIPTION |
|--|--------------|--|-----|---|--------------------------------------|
| SHEATHS | | | | | |
| GLIDESHEATH SLENDER® Hydrophilic Coated Introducer Sheath Nitinol Kit (5 per box) | | | | | |
| | 50-1060 | 6 Fr, 10 cm, 0.021" GW | | GR3534 | Standard, 0.035", 350 cm, 3 mm J Tip |
| | 50-1070 | 7 Fr, 10 cm, 0.021" GW | | GR3538 | Standard, 0.035", 400 cm, 3 mm J Tip |
| GLIDESHEATH SLENDER® Hydrophilic Coated Introducer Sheath Access Kit (5 per box) | | | | | |
| | 60-1060 | 6 Fr, 10 cm, 0.021" GW | | GR3545 | Standard, 0.035", 450 cm, 3 mm J Tip |
| | 60-1660 | 6 Fr, 16 cm, 0.021" GW | | GS3536 | Stiff, 0.035", 350 cm, Straight |
| | 60-1070 | 7 Fr, 10 cm, 0.021" GW | | GS3541 | Stiff, 0.035", 400 cm, Straight |
| | 60-1670 | 7 Fr, 16 cm, 0.021" GW | | GS3546 | Stiff, 0.035", 450 cm, Straight |
| | 60-1065 | 6 Fr, 10 cm, 0.025" GW | | GS3535 | Stiff, 0.035", 350 cm, 45° Angle |
| | 60-1075 | 7 Fr, 10 cm, 0.025" GW | | GS3540 | Stiff, 0.035", 400 cm, 45° Angle |
| GLIDESHEATH SLENDER® Hydrophilic Coated Introducer Sheath Stainless Steel Kit (5 per box) | | | | | |
| | 80-1060 | 6 Fr, 10 cm, 0.021" GW | | GS3543 | Stiff, 0.035", 450 cm, 45° Angle |
| | 80-1660 | 6 Fr, 16 cm, 0.021" GW | | GS3533 | Stiff, 0.035", 350 cm, 1.5 mm J Tip |
| | 80-1070 | 7 Fr, 10 cm, 0.021" GW | | GS3539 | Stiff, 0.035", 400 cm, 1.5 mm J Tip |
| | 80-1670 | 7 Fr, 16 cm, 0.021" GW | | GS3544 | Stiff, 0.035", 450 cm, 1.5 mm J Tip |
| | 80-1065 | 6 Fr, 10 cm, 0.025" GW | | CATHETERS | |
| | 80-1075 | 7 Fr, 10 cm, 0.025" GW | | GLIDECATH® Hydrophilic Coated Catheter (5 per box) | |
| R2P™ DESTINATION SLENDER™ Guiding Sheath (1 per box) | | | | | |
| | GS-R65T1C12W | 119 cm, 6 Fr, Straight | | CG430 | 150 cm, 4 Fr, PV Multicurve |
| | GS-R65T1C15W | 149 cm, 6 Fr, Straight | | CG431 | 150 cm, 4 Fr, Straight |
| R2P™ SLENGUIDE™ Guiding Catheter (1 per box) | | | | | |
| | WG-575T1L23A | 120 cm, 7 Fr, Straight | | CG432 | 150 cm, 4 Fr, Angle |
| | WG-575T1L53A | 150 cm, 7 Fr, Straight | | CG433 | 150 cm, 4 Fr, Pigtail |
| GUIDEWIRES | | | | | |
| GLIDEWIRE® Hydrophilic Coated Guidewire (5 per box) | | | | | |
| | GR3536 | Standard, 0.035", 350 cm, Straight | | CG434 | 150 cm, 4 Fr, Cobra 2 |
| | GR3541 | Standard, 0.035", 400 cm, Straight | | CG435 | 150 cm, 4 Fr, JR 4.0 |
| | GR3546 | Standard, 0.035", 450 cm, Straight | | CG436 | 150 cm, 4 Fr, IMA |
| | GR3535 | Standard, 0.035", 350 cm, 45° Angle | | NAVICROSS® Support Catheter (1 per box) | |
| | GR3540 | Standard, 0.035", 400 cm, 45° Angle | | NC35130 | Straight Tip, 0.035" x 135 cm |
| | GR3543 | Standard, 0.035", 450 cm, 45° Angle | | NC35131 | Angle Tip, 0.035" x 135 cm |
| | GR3533 | Standard, 0.035", 350 cm, 1.5 mm J Tip | | NC35150 | Straight Tip, 0.035" x 150 cm |
| | GR3539 | Standard, 0.035", 400 cm, 1.5 mm J Tip | | NC35151 | Angle Tip, 0.035" x 150 cm |
| | GR3544 | Standard, 0.035", 450 cm, 1.5 mm J Tip | | PTA BALLOON | |
| R2P™ METACROSS® RX PTA Balloon Dilatation Catheter (1 per box) | | | | | |
| | BD-P30020ER | 200 cm, 6 Fr, 3 mm x 20 mm | | BD-P30040ER | 200 cm, 6 Fr, 3 mm x 40 mm |
| | BD-P30060ER | 200 cm, 6 Fr, 3 mm x 60 mm | | BD-P30080ER | 200 cm, 6 Fr, 3 mm x 80 mm |
| | BD-P30100ER | 200 cm, 6 Fr, 3 mm x 100 mm | | BD-P30120ER | 200 cm, 6 Fr, 3 mm x 120 mm |
| | BD-P30150ER | 200 cm, 6 Fr, 3 mm x 150 mm | | BD-P30200ER | 200 cm, 6 Fr, 3 mm x 200 mm |

TERUMO
INTERVENTIONAL SYSTEMS




510(k) Clearance Letters & R2P MISAGO PMA Letter

Glidesheath Slender®
Hydrophilic Coated Introducer Sheath




K173831

R2P Destination Slender®
Guiding Sheath




K193125

R2P SlenGuide®
Guiding Catheter




K161546

Glidewire®
Hydrophilic Coated Guidewire




K152740

Glidecath®
Hydrophilic Coated Catheter




K090040

Navicross®
Support Catheters




K110540

R2P Metacross® RX
PTA Balloon Dilatation Catheter




K163479

R2P Crosstella® RX
PTA Balloon Dilatation Catheter




K152873

R2P Misago®
RX Self-expanding Peripheral Stent



P140002/S002

TR Band®
Radial Compression Device



K152525

Instructions for use





Glidesheath Slender®
Hydrophilic Coated Introducer Sheath

R2P Destination Slender®
Guiding Sheath

R2P SlenGuide®
Guiding Catheter

Glidewire®
Hydrophilic Coated Guidewire

Glidecath®
Hydrophilic Coated Catheter








Navicross®
Support Catheters

R2P Metacross® RX
PTA Balloon Dilatation Catheter

R2P Crosstella® RX
PTA Balloon Dilatation Catheter

R2P Misago®
RX Self-expanding Peripheral Stent

TR Band®
Radial Compression Device



Service and Support

- ◆ Terumo Medical Corporation- Corporate Headquarters
265 Davidson Ave
Suite 320
Somerset, NJ 08773
- ◆ United States Customer Care
Phone:800-888-3786
- ◆ Inside Sales
Phone:800-862-4143
- ◆ Accounts Receivable
Phone:800-283-7866