

Navicross®

Support Catheters

**MASTER YOUR
PROCEDURES**
WITH PRECISE
PERFORMANCE



NAVICROSS® 0.018"
Coming Soon

 **TERUMO**
INTERVENTIONAL
SYSTEMS

PUSHING BOUNDARIES


Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**



NAVICROSS® 0.018" - Coming Soon

| PRODUCT CODE | WIRE COMPATIBILITY | CATHETER LENGTH | TIP SHAPE |
|--------------|--------------------|-----------------|-----------|
| NC18650 | 0.018" | 65 cm | Straight |
| NC18651 | 0.018" | 65 cm | Angle |
| NC18900 | 0.018" | 90 cm | Straight |
| NC18901 | 0.018" | 90 cm | Angle |
| NC18130 | 0.018" | 135 cm | Straight |
| NC18131 | 0.018" | 135 cm | Angle |
| NC18150 | 0.018" | 150 cm | Straight |
| NC18151 | 0.018" | 150 cm | Angle |



Double-braided stainless steel

Engineered for optimal pushability and torque control for lesion crossing

Telescoping Capability

NAVICROSS® 0.018" is designed to telescope through NAVICROSS® 0.035"

Straight and Angled Tips

Tips allow access to vascular branches, including BTK collaterals

Three Radiopaque Markers

- Markers facilitate accurate assessment of position
- Unique spacing provides easy measurement of common stent and balloon sizes
- Initial marker is 1 mm from distal tip; 40 mm and 60 mm spacing

FIND OUT MORE |  Phone: 800.862.4143  terumo.com  Fax: 800.411.5870

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2018 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. PM-01071

 **TERUMO**
INTERVENTIONAL
SYSTEMS