Pinnacle Destination Guiding Sheath

GUIDE YOUR
PROCEDURES
THROUGH
CONSISTENT
RELIABILITY





Terumo Interventional Systems is committed to your success with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve better outcomes for more patients.



Pinnacle Destination

Guiding Sheath

ALLOWS SMOOTH TRANSITION AND EASY PENETRATION

Atraumatic tip¹

• Minimizes potential for vessel damage

Smooth transitions¹

• Guidewire-to-dilator and dilator-to-sheath

Provides lesion access¹

• Exceptional trackability, even in the most challenging anatomy

Stainless Steel Coil

PTFE Inner Layer

Atraumatic tip

Gold Coil Marker

SHAPES AND SIZES TO MEET YOUR DIAGNOSTIC CHALLENGES

Maintains shape and positioning¹

• Outstanding durability, even in the longest procedures

Enhances visualization for precise positioning¹

 Triple radiopacity (sheath, dilator, and gold coil marker)

Resists kinking¹

 Stainless steel coil for consistent reliability

Provides smooth movement¹

• Coil reinforced tubing and a PTFE inner layer that minimizes friction



Tuohy-Borst valve



Accommodates anatomical variations in vessel access

- Straight
- Hockey Stick
- Multi-Purpose
- RDC (Renal Double Curve)
- LIMA (Left Internal Mammary Artery)





Designed to work with Navicross and Glidewire Advantage

Nylon Outer Layer

Hydrophilic Coating

Pinnacle Destination

Guiding Sheath

5Fr 0.076" (1.92 mm)







PINNACLE® DESTINATION® Guiding Sheaths

Packed 1 per box and includes dilator. All dilators are 0.038" wire compatible.

French Size	Product Code	Length (cm)	Hydrophilic Coating Length (cm)	Valve Type	Curve Style	Curve Shape
5 Fr	54-54501	45	35	Cross-Cut	Straight	
	54-54502	45	5	Cross-Cut	Hockey Stick	
	54-54503	45	5	Cross-Cut	Multi-Purpose	
	54-54504	45	5	Cross-Cut	RDC	
6 Fr	54-64501	45	35	Cross-Cut	Straight	
	54-64506	45	35	Tuohy-Borst	Straight	
	RSR01	45	5	Cross-Cut	Straight	
	RSR02	45	5	Cross-Cut	Hockey Stick	
	RSR03	45	5	Cross-Cut	Multi-Purpose	
	RSR07	45	5	Tuohy-Borst	Straight	
	RSR08	45	5	Tuohy-Borst	Hockey Stick	
	RSR09	45	5	Tuohy-Borst	Multi-Purpose	
	RSR13	45	5	Cross-Cut	RDC	
	RSR14	45	5	Cross-Cut	LIMA	
	RSR17	45	5	Tuohy-Borst	RDC	(
	RSR18	45	5	Tuohy-Borst	LIMA	
	RSP01	65	35	Cross-Cut	Straight	
	RSP03	65	35	Tuohy-Borst	Straight	
	RSC01	90	15	Tuohy-Borst	Straight	
	RSC03	90	15	Tuohy-Borst	Multi-Purpose	\
	RSC05	90	15	Cross-Cut	Straight	
	RSC07	90	15	Cross-Cut	Multi-Purpose	
7Fr	54-74501	45	35	Cross-Cut	Straight	
	54-74506	45	35	Tuohy-Borst	Straight	
	RSR04	45	5	Cross-Cut	Straight	
	RSR05	45	5	Cross-Cut	Hockey Stick	
	RSR06	45	5	Cross-Cut	Multi-Purpose	~
	RSR10	45	5	Tuohy-Borst	Straight	
	RSR11	45	5	Tuohy-Borst	Hockey Stick	
	RSR12	45	5	Tuohy-Borst	Multi-Purpose	
			5	Cross-Cut	RDC	
	RSR15 RSR16	45 45	5	Cross-Cut		,
			5		LIMA RDC	
	RSR19 RSR20	45 45	5	Tuohy-Borst	LIMA	,
				Tuohy-Borst		
	RSP02	65	35	Cross-Cut	Straight	
	RSP04	65	35	Tuohy-Borst	Straight	
	RSC02	90	15	Tuohy-Borst	Straight	
	RSC04	90	15	Tuohy-Borst	Multi-Purpose	
	RSC06	90	15	Cross-Cut	Straight	
	RSC08	90	15	Cross-Cut	Multi-Purpose	
8Fr	54-84501	45	35	Cross-Cut	Straight	
	54-84506	45	35	Tuohy-Borst	Straight	
	54-86501	65	35	Cross-Cut	Straight	
	54-86506	65	35	Tuohy-Borst	Straight	
	54-89001	90	60	Cross-Cut	Straight	
	54-89006	90	60	Tuohy-Borst	Straight	

FIND OUT MORE | Phone: 800.888.3786 terumois.com





RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

Reference:

1. Data on File.

