




Glidewire GT-R™

Guidewire

RESHAPING PRECISION

Designed with enhanced shapeability¹⁻²

- Soft, flexible tip to track through visceral anatomy with ease
- Superior tip shape retention
- Kink-resistant, nickel-titanium core
- Available in up to 220cm for radial access

*0.016" and 0.018" diameters compatible with
the PROGREAT® family of microcatheters*

 **TERUMO**
INTERVENTIONAL
SYSTEMS

PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**

GLIDEWIRE GT-R™ Guidewire

PRODUCT CODE	DIAMETER	LENGTH	TIP SHAPE	TAPER LENGTH	RADIOPAQUE LENGTH
GT*ES1618A	0.016"	180 cm	Straight/ Shapeable	25 cm	3.5 cm
GT*ES1818A	0.018"	180 cm		30 cm	
GT*ES1622A	0.016"	220 cm		25 cm	
GT*ES1822A	0.018"	220 cm		30 cm	

Superior hydrophillic coating



3.5cm gold coil
at distal tip
for optimal visibility

FIND OUT MORE



US: 800.888.3786
CA: 833.883.7866



US: terumo.com
CA: terumocanada.ca

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

References:

1. Data on file.
2. Glidewire GT-R 510(k).

©2025 Terumo Medical Corporation and Terumo Medical Canada Inc. All rights reserved.
All brand names are trademarks or registered trademarks of their respective owners. PM-09548

 **TERUMO**
INTERVENTIONAL
SYSTEMS